



**LET'S GET STARTED**  
**Guardian Alert 911 Plus**





Congratulations!

You're a few steps away from staying connected to emergency services

## Contents [ [hide](#) ]





- 1 Device Features
- 2 INCLUDED ACCESSORIES
- 3 DEVICE ACTIVATION
- 4 CONFIRMING BATTERY STATUS
- 5 Tips and Reminders
- 6 Customer Support
- 7 Documents / Resources
  - 7.1 References

## Device Features

	Built-in fall detection technology
	Two-way voice communication
	Call emergency services at the push of a button
	Powered by 4G LTE Technology



## STATUS LIGHT CHART

-  **SOLID RED:** Fully Charged
-  **SOLID BLUE:** Placing a Call
-  **FLASHING RED:** Charging
-  **FLASHING BLUE:** Fall

Detection is ON

## INCLUDED ACCESSORIES

Charging Base

Lanyard

Belt Clip

Setting Up a New Device

## DEVICE ACTIVATION

If your device was purchased anywhere other than LogicMark. com, please call our dedicated customer support team at 1 [800-519-2419](tel:800-519-2419). They will assist you in activating your device and purchasing any necessary subscriptions.

### SETTING UP A DEVICE

**STEP 1:** Plug the charging base into a power outlet.

**STEP 2:** Place your device into the charging base with the charging contacts touching inside the base. Guardian Alert 911 Plus will say “Your Device is Powering On”. It will then say “Your Device is Ready”.

**STEP 3:** Charge Guardian Alert 911 Plus fully.

**STEP 4:** Perform a test call by pressing and holding the emergency button until it vibrates. After a short pause you will speak with an emergency services representative.

## CONFIRMING BATTERY STATUS

**STEP 1:** Place your device in the charging base until it announces “Charging”.

**STEP 2:** Remove it from the base and press the battery test button.

**STEP 3:** Your device will announce “Battery Low” or “Battery Okay”.

## Tips and Reminders

Charge your Guardian Alert 911 Plus device nightly to ensure it is always charged in an emergency.

It is essential to test your Guardian Alert 911 Plus every month to ensure that it is functioning properly.

Your Guardian Alert 911 Plus has fall detection included as a feature. No fall detection device can detect 100% of falls. You should always push the emergency button to call

for help whenever possible.

## Customer Support

### CONTACT US VIA PHONE, EMAIL OR ONLINE

To change device settings Customer Support must be contacted.

Call: 1 [800-519-2419](tel:800-519-2419)

Email: [info@logicmark.com](mailto:info@logicmark.com)

Go to: <https://www.logicmark.com/customer-support/>

### WARRANTY INFORMATION

To register your FREE 1-year limited warranty, go to:

<https://www.logicmark.com/warranty-manuals/>

### DISCLAIMERS

Privacy Policy: <https://www.logicmark.com/privacy-policy/>

Terms and Conditions: <https://www.logicmark.com/terms/>




Guardian Alert 911 Plus

Connected. Safe. Independent.

Copyright © 2025 LogicMark, Inc.

## Documents / Resources

	<a href="#">Logic Mark 911 Plus Medical Alert Devices [pdf]</a> Instruction Manual 911 Plus Medical Alert Devices, 911 Plus, Medical Alert Devices, Alert De vices
---	--

## References

- [User Manual](#)



Leave a comment

Your email address will not be published. Required fields are marked \*

Comment \*

Name

Email

Website

☐ Save my name, email, and website in this browser for the next time I comment.

**Post Comment**

**Search:**

e.g. whirlpool wrf535swhz

**Search**

[Manuals+](#) | [Upload](#) | [Deep Search](#) | [Privacy Policy](#) | [@manuals.plus](#) | [YouTube](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.