



---

[Home](#) » [Linux](#) » **Linux LX26 2 Head Deep Tissue Massage Instruction Manual** 

## Contents [ [hide](#) ]

- [1 Linux LX26 2-Head Deep Tissue Massage](#)
- [2 Introduction](#)
- [3 What's In The Box?](#)
- [4 Terms of Use](#)
- [5 WARRANTY](#)
- [6 Benefits](#)
- [7 Precautions](#)
- [8 BATTERY SAFETY GUIDELINES](#)
- [9 Product Overview](#)
- [10 Operations](#)
- [11 Control Pad](#)
- [12 Care, Maintenance & Cleaning Instructions](#)
- [13 Troubleshooting](#)
- [14 Linux Wellness Products Collection](#)
- [15 Documents / Resources](#)
  - [15.1 References](#)



## Linux LX26 2-Head Deep Tissue Massage



Extend your guarantee to 2 years—just text your order number at [307-461-9622](tel:307-461-9622), and you're all set! It's quick, easy, and hassle-free.

### **US Based Customer Support:**

[307-461-9622](tel:307-461-9622) [SUPPORT@LUNIXINC.COM](mailto:SUPPORT@LUNIXINC.COM)

Lunix Replacement Parts available at: [www.lunixinc.com](http://www.lunixinc.com)

### **Congratulations on your purchase of the Lunix LX26 2-Head Deep Tissue Massager!**

This advanced massager is designed to squeeze away the stress and tension of everyday life, offering you a powerful and unparalleled massage experience.

Your satisfaction is our top priority, and we are here to assist you with any questions or concerns.

We are a small American family business in California rooted in the core value of enhancing our customers' lives. Our dedicated team strives daily to create thoughtful, high-quality products, giving you confidence in your purchase. After all, our products are crafted with care for those who matter most: OUR FAMILY & FRIENDS!

### **Introduction**

Thank you for purchasing the LX26 2-Head Deep Tissue Massager.

This manual contains instructions on operating and maintaining your model so that you can enjoy it for years. However, please be advised that mishaps could happen in time. Rest assured that LUNIX is a team of professionals committed to providing the highest factory support possible.

We know you're excited about using your LX26 Massager, but you must take some time to read through the User Manual. This manual contains all the necessary setup and operating procedures that will allow you to use your new device. Also, follow the precautions and warnings in this manual and any labels or tags attached to your model. They are there to educate you on how to operate your unit safely and also get maximum life and performance from it.

Thank you again for choosing LUNIX. We work hard every day to ensure the highest level of customer satisfaction possible. We genuinely want you to enjoy your new LX26 2-Head Deep Tissue Massager!

You can find more information on our website at

[www.lunixinc.com](http://www.lunixinc.com)

If you have any questions about your model or its operation, contact LUNIX Customer Support by email :

[support@lunixinc.com](mailto:support@lunixinc.com)

## **What's In The Box?**



PRODUCT NAME: 2-Head Deep Tissue Massager

MODEL: LX26

POWER: 10W

WEIGHT: 5 lbs / 2.3kg

VOLTAGE: 5V 2A

BATTERY: 11.1V 2500mAh

PRODUCT DIMENSIONS: 53.54×5.9×5.19 in / 136x15x13.2 cm

## Terms of Use

This product is supplied by LUNIX to the purchaser with the understanding that the purchaser accepts the responsibility that using the LX26 2-Head Massager with its accessories in a careless, improper, or unsafe manner can result in serious injury or death. Also, the purchaser assumes all liability resulting from misuse, unsafe handling, failure to follow instructions, or any action violating applicable laws or regulations. LUNIX shall not be liable for personal injury, property loss, or life resulting from using this product under any circumstances, including intentional, reckless, negligent, or accidental behavior. LUNIX shall also not be liable for any special, indirect, incidental, or consequential damages arising from the assembly, installation, or use of this product or any accessory or chemical required to use this product. By the act of operating & using

the product, the user accepts all resulting liability and releases LUNIX of any and all liability associated with its use.

As a user, you have the power to decide whether to accept the liability for using the LX26 2-Head Massager. If you do not, we kindly request that you refrain from using this product. Please do not open any of the enclosed materials. Instead, return the model to LUNIX. However, please be aware that we cannot accept an item for return or exchange after it has been used or is no longer in as-new condition.

All information in this User Manual is subject to change without notice. For the most updated version, please visit our website at <https://lunixinc.com/pages/user-manuals>.

LUNIX is committed to product improvement and reserves the right to make changes without incurring any obligation to incorporate such improvements into products previously sold. If you have any questions about your model or its operation, please feel free to contact our dedicated Customer Support team by E-mail: [support@lunixinc.com](mailto:support@lunixinc.com)

## **FCC Compliance**

This device contains a module that complies with the limits for a Class B digital device as described in part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

The limits for a Class B digital device are designed to provide reasonable protection against harmful interference in residential settings. This product generates, uses and can radiate radio frequency energy, and, if not operated in accordance with the instructions, may cause harmful interference to radio communications. The user is cautioned that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## **WARRANTY**

To guarantee you have the best experience, we offer a one-year warranty covering

defects in materials, workmanship, and assembly on all our products\*. Please read and understand this manual to ensure you uphold the warranty terms.

\*If a defect in materials, workmanship, or assembly is found BEFORE the product has been used, we will exchange the entire product if it is returned complete and in perfect condition in its original packaging. \*\*

IT IS THE PARENTS' OR GUARDIANS' RESPONSIBILITY TO ENSURE MINORS RECEIVE APPROPRIATE GUIDANCE AND SUPERVISION WHEN USING THIS PRODUCT.

LUNIX DOES NOT operate a new-for-old warranty once the product has been used. Warranty is only valid if the items were purchased directly from LUNIX (purchases from third-party sellers will not be covered under the warranty).

## **COVERED BY WARRANTY**

– 1-year warranty for the massager (excludes damages caused by the user, including but not limited to incorrect maintenance.)

THE LX26 IS A COMPLEX HIGH-PERFORMANCE MASSAGER AND MUST BE TREATED WITH CARE AND ATTENTION.

LUNIX has made every effort to ensure your product is as solid and durable as possible. However, it is still possible that some mishaps could happen. If the program starts to malfunction within the warranty period, please contact us at [support@lunixinc.com](mailto:support@lunixinc.com), and we will provide a replacement.

## **NOT COVERED BY WARRANTY**

- All other parts not mentioned above
- Lack of maintenance, improper use or abuse
- Water damage

Extend your guarantee to 2 years—just text your order number at [307-461-9622](tel:307-461-9622), and you're all set! It's quick, easy, and hassle-free.

## **Benefits**

## **Comprehensive Muscle Relief and Tailored Relaxation**

The LX26 2-Head Deep Tissue Massager features three interchangeable heads for precise muscle targeting, acupoint relief, and muscle shaping. With adjustable modes and intensity levels, it offers a customizable experience for gentle relaxation and deep tissue relief, promoting overall muscle health.

## **Convenient, Portable, and Versatile**

The massager's durable design and wireless USB charging provide convenient use at home, the office, or on the go. Its versatility ensures professional-quality massage and consistent relief wherever you are.

## **Promotes Well-Being and Empowers Your Wellness**

The LX26's deep tissue massage aids muscle recovery, stress relief, and enhances your wellness routine. Its intuitive controls and versatile functions make self-care easy and effective.

## **A Thoughtful and Practical Gift**

The LX26's user-friendly, multifunctional design makes it more than a massager—it's a thoughtful gift for relaxation and well-being, perfect for enhancing any self-care routine.

## **Massage Guidelines**

Limit your session to 15 minutes per area, then move to another spot to avoid overworking. After the massage, drink warm water to aid metabolism. Start with a low intensity and gradually increase as your body adjusts. If you feel unwell, consult a doctor immediately.

## **Important Considerations**

The massager's effects can vary from person to person; some may experience significant relief, while others might find the results less noticeable. If you have specific health conditions, consult a healthcare professional before using the device. It's crucial to follow the guidelines and contraindications provided, as overuse or incorrect use can lead to potential issues. Always use the massager in moderation and seek professional advice to ensure a safe and beneficial experience.

## **Precautions**

The LX26 2-Head Massager is a personal care product not intended to diagnose or treat health conditions. The information provided here is not a substitute for professional medical advice, diagnosis, or treatment. Pay attention to your body, stay hydrated, and consult a healthcare professional, especially if you have specific health conditions or concerns, are pregnant, or are concerned about its suitability for your specific situation. Always follow the instructions and precautions outlined in the manual to ensure safe operation. You are solely responsible for following these instructions and precautions.

1. **Initial Inspection:** After unpacking, carefully inspect the massager and its accessories for any damages. Do not use the device if you notice any issues. Read Before Use: Thoroughly read the user manual for specific instructions and safety warnings before using the massager.
2. **Health and Medical Considerations:**
  1. If you are pregnant, have diabetes, a pacemaker, heart disease, osteoporosis, or any condition affecting skin sensitivity or pressure tolerance, consult a healthcare provider before use.
  2. Discontinue use immediately if you experience any pain or discomfort, and consult a healthcare provider.
3. **Usage Precautions:**
  1. Start with the lowest setting and gradually adjust as your body adapts to the pressure and heat.
  2. Limit each session to 15 minutes per body area, and avoid overuse.
  3. This massager is designed for seated use only; standing on it may cause damage or injury.
4. **Charging Instructions:**
  1. Use only the provided USB cable for charging, ensuring it fits securely into the USB port.
  2. Never charge the device with a USB plug or cable exceeding 5V 2A, as this could damage the unit and cause fire.
  3. Avoid charging overnight or for extended periods during the day.
  4. Do not charge the device on wooden surfaces or near liquids to prevent spills and damage.
  5. Do not leave the massager unattended while plugged in.
  6. Never use the massager while charging, as this can affect the performance and



damage the unit.

#### **5. Hygiene and Maintenance:**

1. Regularly inspect the device for any wear or damage, particularly the electrical cords.
2. Clean the massager after each use, especially if used by multiple people.

#### **6. General Safety Tips:**

1. Avoid falling asleep while using the massager to prevent prolonged pressure on your body.
2. Keep liquids away from the device to prevent spills and damage.
3. Do not wear jewelry while using the massager.

#### **7. Operating Conditions:**

1. Do not use the massager on wet skin, humid environments, or near heated surfaces, as this could cause electric shock or malfunction.
2. Avoid using the device if you've taken painkillers, sedatives, or alcohol.
3. Never attempt to disassemble or repair the unit or its charging cord, as this will void the warranty.
4. Do not use the massager in wet or humid rooms, at high temperatures, or near corrosive gases.

8. **Use with Caution:** Use the massager only for its intended purpose, and always apply common sense to ensure safe and effective use.

### **SAFETY WARNING:**

#### **Check Health Guidelines:**

- Do not use this massager as a substitute for medical treatment. If you are experiencing unexplained pain, receiving treatment for circulation or vascular issues, or using at-home medical equipment, refer to your healthcare provider before use.
- Seek medical advice before using the massager if you have any physical health concerns or if you have any of the following conditions:
  - Pregnancy
  - Pacemakers or other medical implants
  - Heart or vascular disease
  - Irregular blood pressure or are receiving treatment to manage blood pressure.

- Osteoporosis
- Cancer
- Recent surgery

## **Usage Restrictions:**

- Do not use this product on swollen or inflamed skin, open wounds, or if you have a physical condition that limits mobility.
- Never use the massager on individuals with sensory deficiencies that affect their bodies.
- Do not leave the appliance unattended, especially around children.
- Limit usage to 15-minute sessions, with no more than 60 minutes of continuous use.
- For those with sensitive skin, start with the lowest intensity setting.
- General Safety Guidelines:
- Carefully read and follow all the manufacturer's instructions, warnings, and precautions before using the massager.
- Children under 16 should not use the massager without the supervision of a responsible adult. Remember, this product is not a toy.

## **BATTERY SAFETY GUIDELINES**

This massager is a powerful device that demands responsible and careful use to ensure your safety. Any form of improper handling or careless operation can lead to serious injury or even death, underlining the gravity of the situation.

- The massager contains a lithium battery, a potentially hazardous component if not handled correctly. It is of utmost importance to always adhere to the provided instructions and warnings regarding battery care and use.
- Prior to charging, it is crucial to ensure that the charger and plug settings are compatible with the battery type and specification (5V, 2A). Failure to do so can result in overcharging and potential fire hazard.
- Always charge the battery in a well-ventilated area, away from flammable materials, and never leave it unattended during charging. Keep the device out of reach of small children while charging.
- Do not disassemble, crush, short-circuit, or expose the battery to fire or other ignition

sources.

- If the device becomes hot to the touch during charging (above 110°F / 43°C), disconnect it immediately and stop charging.
- Keep a Class D fire extinguisher nearby in case of an emergency.
- Avoid exposing the device to water or moisture to prevent electrical hazards.

### **WARNING: RISK OF INJURY!**

We strongly recommend that people who have very sensitive skin use the lowest intensity.

### **! CAUTION! DANGER! FIRE HAZARD!**

CHARGING OR DISCHARGING BATTERIES POSES A NUMBER OF SIGNIFICANT HAZARDS WHICH CAN INCREASE THE RISK OF FIRE, EXPLOSION, SERIOUS INJURY, AND PROPERTY DAMAGE IF NOT PERFORMED AS PER THE INSTRUCTIONS. BEFORE USE, READ AND FOLLOW ALL MANUFACTURER'S INSTRUCTIONS, WARNINGS AND PRECAUTIONS.

NEVER ALLOW CHILDREN UNDER THE AGE OF 16 TO USE THE MASSAGER WITHOUT THE SUPERVISION OF A RESPONSIBLE, KNOWLEDGEABLE ADULT. THIS PRODUCT IS NOT A TOY.

### **Product Overview**



## Operations

The LuniX LX26 2-Head Massager is designed as a personal care product with your comfort and well-being in mind. It is not intended to diagnose or treat any health condition. Operate the LuniX Massager sensibly and with care, and it will be exciting, safe, and relaxing for you. Failure to operate your model safely and responsibly may result in property damage and serious injury. We strongly recommend that people with very sensitive skin use the lowest intensity.

If you have physical health concerns, please consult your physician before using the LuniX LX26 2-Head Massager.

### **To use the 2-Head Massager, follow these steps:**

1. Fully charge the massager before first use. Only use the charger provided; do not use the device while charging, as it may affect performance or damage the unit. The power button flashes while charging and will have a solid green light when fully

charged.

2. Choose the preferred Massage Head from the 3 heads included in the package. Align the required head shaft with the socket and firmly press until the shaft is fully inserted. To remove it, pull the attachment and withdraw it from the socket.
3. Press the Power button for 3 seconds to turn on the massager. The default mode is Massage Mode 1 and Speed Intensity Level 1. Press the Massage Mode and Intensity Level buttons to adjust the settings.
4. The Massager is pre-programmed to run for 15 minutes and will automatically shut off after the session. If an additional session is desired, restart the massager.

### **Caution:**

Limit massage to 15 minutes on any body part and slowly transition to other areas. After use, drink warm water to help metabolism.

Strictly Not for Children: This device is not suitable for children under the age of 16. It must always be used under the supervision of a responsible adult. Consult your doctor immediately if you feel unwell after use.

### **How to Determine the Battery Capacity of the Massager:**

When the unit is off, short press the power button and check the light indicator:

- Flashing Red Light: 0-5% (unit won't turn on)
- Steady Red Light for 5 seconds: 6-20% (a prompt to charge the unit)
- Steady Yellow Light for 5 seconds: 21-50%
- Steady Green Light for 5 seconds: More than 50%

When the unit is on, look at the power button; the color display corresponds to the battery capacity.

Red Light: Under 20% Yellow Light: 20-50% Green Light: Over 50%

### **Control Pad**



**Note:** The massager has 3 interchangeable massage heads that can be used in different parts of the body:

- A. Flat Massage head for shaping various muscle parts.
- B. Double Ball Head for massaging the acupoints on the body.
- C. W-shaped PU Self-Formed Leather Massage Head for massaging arms, waist, back buttocks, and thighs.

## Care, Maintenance & Cleaning Instructions

- After use, wipe the massager with a soft, damp cloth and non-abrasive cleaners, then store it in a cool and dry place.
- Do not place the unit into any liquid.
- Keep away from all solvents and harsh detergents.



## Troubleshooting

When operational problems occur, please refer to the chart below.

PROBLEM	POSSIBLE CAUSE	SOLUTION
My unit suddenly stopped working	Power outage OR massage more than 15min OR long working time, the unit has been shut off due to overheating protection	The unit automatically shuts off after 15min of use. Use the unit after the motor has been naturally cooled
Unit stopped working and will not turn on after pressing the power button.	The unit is out of power or the battery shuts off due to overcurrent protection.	If the unit is out of power, charge it again before using. If it still has battery and will not turn on, charge it for a few seconds to activate.
Unit stopped charging	Unit not receiving the correct amount of voltage	Make sure you use the charger included in the package. If issue persists, contact us at support@lunixinc.com
For more troubleshooting guides, please scan the QR code.		

If you have any questions, concerns or if you would like to order any spare parts or replacements. Please contact us at [support@lunixinc.com](mailto:support@lunixinc.com) and we will get back to you within 24 hours

## Lunix Wellness Products Collection







LX8 2 Layers Orthopedic Bed Wedge Pillow Set



LX10 Cordless Full Leg Compression Massager



LX11 5Pcs Orthopedic Bed Wedge Pillow Set



LX12 3Pcs Orthopedic Wedge Pillow Set for Kids



LX13 6Pcs Orthopedic Bed Wedge Pillow Set



LX14 1Pc Orthopedic Bed Wedge



LX15 14pcs Modular Couch for Kids



LX16 Red Light Therapy Belt



LX17 Infrared Sauna Blanket



LX18 Vibration Plate



LX19 Red Light Therapy for Knee, Elbow or Shoulder




LX20 Shiatsu Foot Massager





## Documents / Resources

	<p><a href="#">Lunix LX26 2 Head Deep Tissue Massager [pdf]</a> Instruction Manual</p> <p>LX26, LX26 2 Head Deep Tissue Massager, LX26, 2 Head Deep Tissue Massager, Deep Tissue Massager, Tissue Massager, Massager</p>
---	--

## References

- [User Manual](#)

Linux

2 Head Deep Tissue Massager, Deep Tissue Massager, Lunix, LX26, LX26 2 Head Deep Tissue Massager, massager, Tissue Massager

## Leave a comment

Your email address will not be published. Required fields are marked \*

Comment \*

Name

Email

Website

☐ Save my name, email, and website in this browser for the next time I comment.

**Post Comment**

**Search:**

**Search**

[Manuals+](#) | [Upload](#) | [Deep Search](#) | [Privacy Policy](#) | [@manuals.plus](#) | [YouTube](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.