WARNING

READ ALL INSTRUCTIONS FOR THE SCALE BEFORE USING THIS QUICK REFERENCE GUIDE.

If you no longer have the instruction manual in your possession, contact Beurer customer service. This quick guide is not a substitute for reading and understanding the instruction manual.

1 Device description

1. Display
2. Conductors
3. "Down" key ▼
4. "SET" key
5. "UP" key ▲
6. Body fat %
7. Muscle percentage
8. Body Mass Index
9. Body water %
10. Bone mass
11. Active metabolic rate
12. Weight
13. Memory position

2 Battery installation

1. If present, remove the battery insulating strip or remove the battery's protective film and insert the batteries according to the polarity.
2. Replace battery cover.
3. Replace batteries when "Lo" is displayed.

3 Setting measurement unit

You can switch to "pound", "stone" and "kilogram" (lb, st, lb, kg) by switching on the scale and immediately pressing the ▼ or ▲ button.

4 Measuring weight only

1. Stand on the scale without moving and distribute your weight evenly on both legs. The scale immediately begins to measure your weight. After flashing, the weight is established.
2. The scale switches off a few seconds after you step down.

5 Setting user data

Your personal parameters must be stored in order to measure body fat and other physical data.

1. Switch on your scale (tap the weighing area) and wait until "0.0" is displayed.
2. Then press "SET". Now the first memory preset flashes in the display.

3. Select the desired preset by pressing the buttons "up" ▲ or "down" ▼ and confirm the desired preset with the "SET" button.
 4. The following steps will allow you to select your desired settings with the "up" ▲ or "down" ▼ buttons and confirm with "SET":
 - Memory position
 - Height in inches or centimeter
 - Age in years
 - Gender setting for "Woman" or "Man"
 - Activity mode
- The scale is ready for measuring.

6 Measuring weight, body fat, body water, muscle percentage, bone mass, BMI and active metabolic rate

Your personal parameters must be stored (see no. 5 Setting user data).

1. Switch on your scale and click "SET".
2. Select the memory preset where your basic personal data is stored by pressing the "up" ▲ or "down" ▼ button repeatedly. The data is displayed consecutively until "0.0" appears.
3. Step onto the scale barefoot and make sure you're standing on both conductors.
4. First your body weight is calculated and displayed.
5. Remain on the scale and avoid moving; your body fat is being analyzed. This may take a few seconds.

The following data will be displayed:

3. "Down" key ▼
4. "SET" key ("SET")
5. "UP" key ▲
6. Body fat %
7. Muscle %
8. Body Mass Index
9. Body water %
10. Bone mass
11. Active metabolic rate

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