
CONGRATULATIONS ON YOUR NEW POWER PERCUSSION DEEP TISSUE MASSAGER

Please read all instructions and warnings prior to use. Keep this as it contains important information for future reference.

WHAT YOU GET

- a. Power Percussion Deep Tissue Massager
- c. Carrying case
- d. Adapter
- b. Interchangeable massage nodes e. Attachment case










CARE INSTRUCTIONS

- a damp cloth.
- · Air-dry only. Never place in a dryer.
- · Spot clean only. Unplug and wipe with · Allow to dry fully before plugging back in.
 - · Do not submerge in water.
 - Do not use detergent or soap when cleaning.

01 To charge the massager, plug the adapter into the massager and then plug into an outlet. Do not use the massager while it is

charging.




• The image shown here is for illustrative purposes only and may be subject to change



the desired massager node.

See reverse side for interchangeable node guide.



03

Press and hold the power button to power on the massager. Press the power button again to cycle through the 6 intensity levels. Press and hold the power button to turn off the massager.

You can change the nodes to customize your massage. Only swap nodes when the massager is off.



PRINTED IN MEXICO Elbrd: J.E. SC: 11-05-2024 SHNMDL: 2437599A_B 1P OBPN: 1014747 8 1 4 1 0 0 3 0 2

04

Use for 60 seconds per region. Do not use on neck, head or any hard or bony areas of the body.

Enjoy your massage for up to 10 minutes.




BALL

Impact large areas for overall comfort.









ROUND

Constructed from durable foam but designed for more precision.









FLAT

The hammer-like impact hits deeper into muscles.




GOOD FOR Pecs Quads



MUSCLE IMPACT Medium

FORK

Gently rock from side to side to work out knots and tense spots.



GOOD FOR Traps, Quads Calves, Feet




High

MUSCLE IMPACT

BAR

The curved edge provides comfort with surprising strength.





GOOD FOR

MUSCLE IMPACT High





SHARPER IMAGE®

INTERCHANGEABLE NODES GUIDE

SAFETY WARNINGS: PLEASE READ BEFORE USE

- · Not a children's item. Keep out of the reach of children.
- Always wear a shirt/clothing when using. Do not use on bare skin to avoid possible irritation.
- As with any massage, bruising may occur. Test your comfort level and check targeted areas often. Stop massage immediately if you feel any increase in physical discomfort or pain.
- · Do not use attachments of any kind other than those that may be included.
- · Do not fall asleep while using.
- · Do not leave this item unattended when it is plugged in.
- · Always allow item to cool if it becomes too hot.
- · Do not expose to extreme temperatures or moisture.
- · Keep away from sources of water and heated surfaces.
- · Never touch the power plug with wet hands.
- · Do not carry this item by its cord.
- · For indoor use only.
- Thoroughly inspect products before every use. If the product appears worn, frayed/splintered, cracked, or broken in any manner, discontinue use and discard immediately.
- · This product is not intended to treat, cure, or prevent medical conditions.
- Do not use this item if you are pregnant or have an injury or illness. Consult your
 physician if you have a medical condition that may prevent safe use of this item.

RECHARGEABLE BATTERY WARNING

- · This product contains batteries that are non-replaceable.
- There is a risk of fire and personal injury if battery is punctured, damaged, or misused.
- · Do not dismantle, crush, incinerate or short-circuit battery.
- · Never expose battery to extreme temperatures or moisture.
- · Keep away from flammable materials and direct sunlight.
- · Only charge using the provided charger, cable, and/or adapters.
- · Rechargeable battery is only to be charged by an adult.
- · Adult supervision is required while charge the battery.
- · Allow battery to cool to room temperature before charging.
- · Charge on a non-flammable surface and keep away from flammable products.
- Discontinue charging if battery becomes hot, swells, smokes, or gives off a strong odor.

- Battery chargers used with the product are to be regularly examined for damage to the cord, plug, enclosure, and other parts. In the event of such damage, chargers must not be used until the damage has been repaired.
- · Recycle or dispose of battery according to federal, state, and local laws.

FCC STATEMENT

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1. This device may not cause harmful interference, and
- this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

WARNING: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.



To reduce the risk of electric shock, always unplug this appliance from the electrical outlet immediately after using and before cleaning.

PRINTED IN MEXICO Elbrd: J.E. SC: 11-05-2024 1P

OBPN: 1014747 SHNMDL: 2437599A B

8 1 4 0 0 0 3 0 2

