




QUICK START GUIDE


IPL Hair Removal Device

90-Day Hassle-Free Return
2-Year Warranty
eusupport@ulike.com



The hair removal system **U**like


Product Introduction



The Power of IPL

IPL stands for Intense Pulsed Light. It is a type of FDA-cleared light treatment used in the hair removal industry for years.

IPL uses intense pulses of a broad spectrum of polychromatic light. It targets the melanin in hair follicles and sets the hair root dormant, bringing permanent hair removal results.




■ What's Difference Between IPL And Laser?

Both IPL and laser are light treatments. The critical difference between laser and IPL is the type of light used. IPL uses a broadband pulsed light source, whereas laser is a monochromatic coherent light source.

Both methods target the melanin in the hair follicle and can achieve permanent results.

The laser can be more precise than IPL and suitable for more skin tones. But since it's high power output, most medical experts recommend doing laser treatments at professional facilities.




Ulike Patented Sapphire Ice-cooling Technology

Since in IPL hair removal treatment, the melanin in the hair follicle will absorb the IPL light and heat up, the heat will set the hair follicle dormant and stop hair from regrowing. So sometimes, during the treatments, people will feel slight pain like a rubber band snapping. Some at-home IPL devices' light window can even reach temperatures of 200 °F during use.

To ease the burning sensation, many experts suggest cooling to soothe the skin and prevent/help the burning.

Ulike is one step ahead in the game. We have developed our patented Sapphire ice-touch technology, which allows the ice-cooling system works along with the IPL system. So the treated area gets cooling applications and effective IPL light simultaneously. Effectively avoiding burning and pain during use.

In Air 3, Ulike has upgraded its exclusive Sapphire ice-touch technology to a higher level. It can effectively lower the epidermal temperature of the effective area to 68 °F and below while emitting high-energy light for 20 minutes straight.




How to Use

Step 1

Prepare


- Remove all visible hair by shaving
- Cleanse the treatment area and pat dry
- DON'T use any gels or creams with this device!



Step 2

Connect the power


- Find a comfortable place
- Connect adapter to the handset
- Plug in the power



Step 3

Goggle up


- Wear the professional IPL Safety Glasses comes in the box



Step 4

Turn on the device


- Short press to turn on the handset
- Once the handset is on, short press to change the intensity levels



Step 5

Full contact


- Place the handset in a vertical position against skin.
- Make sure the skin sensor indicator is light up



Step 6

Start Ice-cooling painfree treatment


- Press the treatment button. Light will pulse from the device. Enjoy the treatment.



Tips on How to Use the Hair Removal Device

■ Finds the Right Level for U

- We strongly suggest you start with Soft Mode and do a test patch.
- Note your skin's reaction and proceed accordingly.
- Increase the level of intensity as needed.



- The higher the level, the faster you will achieve results.

■ Choose the Treatment Method


Auto-Glide Method

- Press and hold for 2 seconds to activate
- Glide the light window on skin, the light will automatically flashes
- Great for large areas such as arms, legs



Spot Method

- Short press the treatment button, one flash at a time
- Suitable for small areas such as fingers and upper lips



■ Follow the Agenda

To get the best results ASAP:

Week 1&2 Treat every other day	Week 3&4 Treat twice per week
--	---

10+ minutes for complete body treatment.

Follow instructions and consistent use for 3-4 weeks


You will definitely notice the difference.

■ Why do I Need to Shave before Using Ulike?


Shaving off the hair will let the treatment have better performance. So the IPL light will work only toward the hair root but not distracted by the surface hair.



Please only use shaving to prep the skin because other methods will remove the hair root, affecting the treatment performance.

If you smell some burning during the treatment, it is probably because the IPL light burns the surface hair stubble/residue.



What to Expect in 4 Weeks

	Shaving	With Ulike	
Week 1	Shave Smooth for a day or two Shave again	Shave and prepare Start treatment every other day Not much change	
Week 2	Shave Smooth for a day or two Shave again	Shave, and prepare. Start treatment every other day. The hair is thinner and softer.	

	Shaving	With Ulike	
Week 3	Shave Smooth for a day or two Shave again	Shave, and prepare. Start treatment twice per week. The hair is growing way much slower, and they are even more thin and soft.	
Week 4	Shave Smooth for a day or two Shave again	Shave and prepare. Start treatment twice per week. The hair is almost not growing back, and there is not much left.	

The mode suitable for skin and hair color

Will IPL Hair Removal Work for Me?

For IPL hair removal to work, there must be a contrast between the skin tone and the hair colour. Users with fair skin and dark hair will achieve the best results. Because the IPL light is absorbed by the melanin and heated up, the hair follicle can be killed, resulting in long-term hair removal.

Red, white, grey, and light blonde hair, which lacks melanin, may not yield effective hair removal results. Because of the lack of melanin, the follicle will not absorb enough light to become dormant.

At the same time, deep dark skin with a lot of melanin may be at risk of burning because the skin absorbs too much light. This is a lamentable limitation of IPL technology.

Take a look at the chart below to see if IPL hair removal will work for you.

Skin color Hair color						
	White	Beige	Light brown	Medium brown	Dark brown	Brown black and darker
White/Grey	×	×	×	×	×	×
Red	×	×	×	×	×	×
Light blond	×	×	×	×	×	×
Dark blond/light brown	regular mode power mode	regular mode power mode	regular mode power mode	soft mode regular mode	×	×
Brown	regular mode power mode	regular mode power mode	regular mode power mode	soft mode regular mode	×	×
Dark brown	regular mode power mode	regular mode power mode	regular mode power mode	soft mode regular mode	×	×
Black	regular mode power mode	regular mode power mode	regular mode power mode	soft mode regular mode	×	×

(Note: X means not applicable)

Warnings

■ Do not Use the Device

- If you are under 18.
- If your skin tone is darker than the recommended tones.
- If your hair color is lighter than the recommended tones.
- If you have skin diseases such as allergic dermatitis, eczema, etc.
- If you have a known skin disorder.
- If you have photosensitive epilepsy or photosensitivity.
- If you have been burned from sun overexposure.
- If you have been naturally or artificially deep-tanned to deep dark skin.
- If you are on your period.
- If you are pregnant/breastfeeding. Although there is no evidence showing that IPL has adverse effects on pregnant or breastfeeding women, we still don't recommend using our handset to be safe.
- On the scalp, eyebrows, or near or around the eyes. Use in these areas may permanently damage your eyes.
- Over blisters, scars, moles, tattoos, and dark brown or black spots.
- On damaged skin, scabs, healing wounds, cancer, or hemangioma.
- On any artificial aesthetic, prosthetic, or plastic area.

If you are still unsure whether this device is safe for you to use, please consult your physician or dermatologist.

If you take any medication regularly and over a long time, please consult your physician for any potential impact on skin sensitivity to light.

■ Eye Safety

Do not look directly into the light window or attempt to activate the device toward the eyes.

For eye safety, the device will only pulse light when the light window is fully covered by skin.

If the device activates when the light window is NOT in contact with the skin, STOP use immediately and contact your retailer.

Possible eye injury (potentially leading to loss of vision or even blindness) or skin injury if instructions are not followed.

Protect the eye from exposure. Wear the included IPL safety glasses before treatment.

Do not use the appliance over the eyelids or close to the eye.