



**MILLION DOLLAR GROWTH FACTOR
MICRO-INFUSION SYSTEM**

YOUR SKIN WILL LOOK LIKE A MILLION DOLLARS!

**ANTI-AGING REJUVENATION SYSTEM
USER GUIDE & MANUAL**

**THANK YOU FOR PURCHASING OUR
MILLION DOLLAR GROWTH FACTOR
MICRO-INFUSION SYSTEM!**

We Want to Ensure Your Best Results!

IMPORTANT:

**Please Read All of the Care & Usage Instructions
Completely Before Beginning Your Treatment.**

Watch Our Treatment Video.

Micro-Infusion Care Instructions

INDEX

What's Included?	3
Micro-Infusion Device	3
Beautify® Million Dollar Growth Factor Serum	4 - 5
Treatment Instructions & Treatment Video	6
Treatment Step 1: Preparing The Serum Vial	7
Treatment Step 2: Priming The Needles	7
Treatment Step 3: Applying The Serum To Your Skin	7
Post-Treatment Care Instructions	8
Important Safety & Usage Guidelines Disclaimer	9
Frequently Asked Questions – FAQ's	10-12
Important Safety & Usage Guidelines	13
Micro-Needle Tip Use Precautions – PLEASE READ	14

IMPORTANT: Takes Pics!

Be Sure to Take Before, During & After Pics:

- **Day 1**
- **Day 15**
- **Day 30**

This allows you to better see your actual progress & results!

Here are some helpful tips when taking pics of your face:

- Take pics at the same distance away from your face each time
- Same lighting each time
- Same camera angle each time
- Use the same expression each time
- Don't wear makeup

You can take Selfies or have someone take your pic.



MILLION DOLLAR GROWTH FACTOR MICRO-INFUSION SYSTEM

Each Micro-Infusion Kit Contains:

- Two 5mL sealed infusion vials of our Beautify® Million Dollar Growth Factor Serum
- Two sterilized, medically-sealed disposable 64 micro-needle Roller tips with caps
- *Optional:* One 1 oz. bottle of our Beautify® Million Dollar Growth Factor Serum for everyday use & optimal results.

24K Gold-Plated Micro-Needles

- 0.5mm needle length for gentle and effective micro-infusion. Our nerves, blood vessels, oil glands are down 0.6mm to 4mm is the dermis layer, *so our 0.5mm micro-needles cannot reach them.* This allows for a virtually pain-free treatment.
- Each device features a Roller Tip with 64 surgical-grade 24k gold-plated stainless steel micro-needles for precision and safety. **Note: Roller tip is only good for 1 treatment.**
- These micro-needles are thinner than a human hair and create tiny, invisible micro-channels on your skin's surface, increasing absorption by 900% stimulating natural Collagen & Elastin production. This reduces lines & wrinkles and firms and tightens skin.

Infusion Vial:

- Contains 4mL of serum to deliver active ingredients during treatment.
- The serum infusion vials & Roller Tips are not re-usable. *Bacteria could get on them and cause treatment issues if re-used and cause serious skin issues.* Please dispose of them properly – Do Not Flush.

**** KEEP OUT OF THE REACH OF CHILDREN ****

Beautify® Million Dollar Growth Factor Serum Ingredients

Water (Aqua), Aloe Vera Gel, Vegetable Glycerin, Butylene Glycol, Carbomer, Polysorbate 20, Hydrolyzed Collagen, Super Peptide Palmitoyl Oligopeptide, Super Peptide Palmitoyl Tetrapeptide-7, Lecithin, Acetyl Glutamine, Epidermal Growth Factors: (sH-Oligopeptide-1, sH-Oligopeptide-2, sH-Polypeptide-1, sH-Polypeptide-9, sH-Polypeptide-11) Hyaluronic Acid, Sodium Hyaluronate, Cyclopentasiloxane, Dimethicone, Phenyl Methicone, Glucosamine Phosphate, Epsom Salt Magnesium Sulfate (Skin Mineral), Salvia Miltiorrhiza Leaf (Red Sage) Extract, Dimethyl Sulfonyl Methane, Folic Acid (Skin Vitamin B9), Magnesium Aspartate (Skin Mineral), Niacinamide (Skin Vitamin B3), Copper Gluconate (Skin Mineral), Dexa Panthenol (Skin Vitamin B5), Hydrolyzed Oat Protein, Hydrolyzed Elastin, Manuka Honey, Citric Acid and Lactobacillus Ferment/Radish Root Ferment (Natural Preservative), Lemon Basil (Natural Scent). (*Beautify® Micro-Infusion treatment serum has no fragrance*)

Note: The Beautify® Micro-Infusion treatment serum is a bit more liquidy (**than our everyday serum**) to allow it to flow better for treatment and has no scent, **because** it could irritate the **open** skin during treatment.

Our 1 oz. bottle of our Beautify® Million Dollar Growth Factor Serum is for everyday use only *and is not to be used for treatment.*

Main Active Ingredients

- ◆ **5 Growth Factors:** This Growth Factors is so Good; **It Actually Won a Nobel Prize!**
- ◆ **2 Super-Peptides:** Clinically-Proven to Grow Collagen Reducing Lines & Wrinkles
- ◆ **Hydrolyzed Collagen:** 75% of Skin is Comprised of Collagen and We Lose 1% a Year, Starting at Age 25. Tightens & Firms Skin
- ◆ **Hydrolyzed Elastin:** 10% of Skin is Comprised of Elastin – Improves Skin Elasticity
- ◆ **2 Hyaluronic Acids:** Plumps & Super-Hydrates the Skin By Attracting & Binding Water 1,000 Times Its Weight in Water. Two HA's are better than One!
- ◆ **3 Special Skin Care Vitamins:** Nourishes the Skin with Anti-Aging Nutrients
- ◆ **3 Special Skin Care Minerals:** Improves Skin Health Giving Skin A Healthy Glow

Five Growth Factors:

Growth Factor sh-Oligopeptide-1: **It's so exceptional; it was awarded a Noble prize.** Clinical evaluations over a 3-Month period, showed a statistically significant improvement in the appearance of fine lines and wrinkles, skin texture, pore size. Known as EGF (Epidermal Growth Factor), it's comprised of 53 Amino acids, which promotes skin repair processes, increasing the rate of self-regeneration, plus, it enhances blood vessel generation around damaged parts of the skin.

Growth Factor sh-Oligopeptide-2: It enhances the biosynthesis of Hyaluronic Acid and other components of the Extracellular Matrix (ECM) developing a smoother, youthful, and healthy appearance of the skin, and it stimulates cell regeneration and skin healing. Known as IGF-1, this special ingredient consists of 70 amino acids, which minimizes the appearance of fine lines and deep wrinkles and reduces blotchiness. Thanks to its anabolic fat-burning effect, sh-Oligopeptide-2 effectively refines texture glides, and slims the face. Provides conditioning and helps maintain skin in good condition. Protects skin from external influences. Helps build up a stronger skin barrier.

Growth Factor SH-Polypeptide-1: Stimulates the growth of Fibroblasts, which are dermal cells that produce Collagen & Elastin. This results in improving the elasticity of the skin and smoothes out lines & wrinkles and it creates firmer, bouncier skin. Contains 288 amino acids which aids in the restoration of the aging tissue. Interacting with EGF and IGF-1 sh-Polypeptide-1 boosts the production of Hyaluronic Acid, exposing a smoother, youthful, and healthy appearance of the skin.

Growth Factor sH-Polypeptide-9: It's everything you would want for Anti-Aging and Rejuvenation, including anti-wrinkle, skin rejuvenation, elasticity enhancement, dermal-moisturizing and cell vitalization properties. Known as VEGF, it stimulates blood vessel formation and increases nutrient flow to fibroblasts, which grow Collagen, while restraining cell death.

Growth Factor sh-Polypeptide-11: A Bioactive Growth Factor that significantly boosts Collagen production and skin regeneration. This special ingredient enhances skin elasticity, reduces fine lines & wrinkles, and improves overall skin texture, and fights Free-Radicals contributing to a more youthful appearance. It also strengthens the skin's structural framework, promoting firmness and resilience. It boosts the proliferation of blood vessel endothelial cells, improving vascularization and blood circulation. **Lastly, it reduces Free Radicals, which are the #1 cause of skin issues, plus, it inhibits the aging of skin cells blocking cell programmed death and function failures. It reinforces the skin's natural defenses.**

Growth Factors Are Clinically Proven to Grow Collagen.

They Noticeably Improve the Appearance of Fine Lines & Wrinkles and Improve the Volume & Plumpness of Skin, Making You Look More Youthful!

Our Million Dollar Serum Also Contains Two Super-Peptides at a Very High Potency Needed to Deliver Real Results!

Super Peptides signal the skin to produce Collagen, reducing lines & wrinkles and work to rebuild and repair damaged cells slowing down the skin aging process. They also have Anti-Oxidant properties to neutralize the silent damage that occurs from free radical damage, UV rays, pollutants, and other toxins.

Super-Peptide Matrixyl Synthe' 6: The #1 Peptide for Growing Collagen!

Super Peptide Palmitoyl Oligopeptide: AKA *Palmitoyl Tripeptide-1*. This latest version stimulates 6 major skin components **showing a 31% decrease* in wrinkles** over 2 months.

**Source: Matrixyl Clinical Trial Study.*

Super-Peptide Argireline: The #1 Peptide for Wrinkle & Fine Line Reduction!

Referred to as a Natural Botox®. This Peptide relaxes the muscles responsible for the formation of expression lines, **reducing wrinkles & lines by 17%* in 15 days.** **Source: Argireline Clinical Trial Study*

Argireline (Palmitoyl Tetrapeptide-7) is known as the "youth hormone molecule" and it provides four major skin benefits:

- **Increases Skin Firmness**
- **Improves Skin Elasticity**
- **Produces Collagen Production**
- **Helps Reduce Inflammation**

This Super-Peptide aids the reduction of fine lines & wrinkles, and promotes a smoother, more youthful complexion. It protects the skin from degradation, preserving skin's natural youthfulness and boosting its replenishment.

**This is What's Needed To Deliver Positive Skin Improvement!
Creating a More Beautiful, More Youthful-Looking You!**

Beautify® Micro-Infusion Treatment Instructions



Scan QR Code for Treatment Video Instructions:



PRE-CARE INSTRUCTIONS: **IMPORTANT**

PLEASE READ ALL INSTRUCTIONS FOR SAFETY & BEST RESULTS

- No Botox or Injectables 2 weeks before and after a Micro-Infusion Treatment.
- Do not lay out in the sun three (3) days before, or seven (7) days (*minimum*) after your Micro-Infusion treatment. If your sunburnt, do not do a treatment.
- If you use any Glycolic/Lactic Acid or Prescription Retin-A Products, discontinue use three (3) days before, *and for three (3) days after treatment*.
- If you use Non-Prescription Retin-A Products, discontinue use at least one (1) day, preferably three (3) days, before treatment and discontinue use for seven (7) days after treatment.
- Cold Sores or Herpes Simplex: If you are highly prone to cold sores, a micro-infusion treatment can sometimes cause a break-out, especially if we are doing a treatment on your lips. We recommend that you take an anti-viral drug such as; Zovirax, Famvir, Valtrex or L-lysine as prescribed by your Physician. *This is at your discretion*. Take the medication as prescribed for 1 week prior and one week after the lip or full face Micro-Infusion Treatment.

IMPORTANT: Takes Before, During & After Pics!

Be Sure to Take Before, During & After Pics on Day 1 - Day 15 & Day 30.
This allows you to better see your actual progress & results!

PRE-TREATMENT - IMPORTANT

1. Tie your hair back and keep it out of your face.
2. **Three Cleanses:** First; double cleanse your Face & Neck (*Two Cleanses*) to ensure all makeup and products are fully removed on skin. *Use makeup remover if necessary*. Then use 70% isopropyl alcohol on Face & Neck with a makeup cotton round. **It's very important to remove every bit of oil, dirt, and bad bacteria to prevent irritation or adverse reactions during treatment as you will have open wounds in your skin.**
3. **Make sure your hands and the area you're using are clean.** If you touch anything besides your skin or the device, like a bottle or brush, it could cross-contaminate. Stick to just touching your face and the tool during the treatment.

STEP 1: PREPARING THE SERUM VIAL FOR TREATMENT

1. Your sealed glass vial is already filled with our Million Dollar Growth Factor Serum. Hold the glass vial firmly: Keep a steady grip on the body of the vial with one hand. Look for the two arrows on the cap.
2. Break the seal: Simply lift up on the clear plastic piece that covers the Gold cap from where the arrow points. Once the seal is broken, the plastic piece and Gold cap should easily pop off. Remove the remaining Gold seal from the top of the vial and dispose of it.
3. Remove the rubber stopper: To access the serum. *Keep Rubber Stopper – Do Not Throw Away.*
4. Open and remove the Roller tip from its sterile package - (*Leave the Roller cap on that covers the micro-needles*). Insert the Roller tip by firmly pressing down into the bottle and snap it into the top of the glass vial...it should be snug and not at an angle.

STEP 2: PRIMING THE MICRO-NEEDLE ROLLER TIP

1. Turn the device upside down, remove the clear cap and *gently* press the Roller needle head onto a clean paper towel laid on a table or Counter and roll it gently on the paper towel. (*There is a spring mechanism that moves up and down that allows the serum flow out*). **Do Not press the roller down too hard as it may bend the needles.** *Keep the Clear Cap handy.*
2. Look for bubbles: You should see bubbles forming in the vial, which shows that the serum is beginning to flow. Keep pressing until you see the serum flowing normally on to the paper towel. While rolling; gently push the vial up & down to pump the serum; this allows it to flow out more easily.
Note: Serum does not flow out excessively.

STEP 3: TREATMENT – MICRO-NEEDLING THE SERUM INTO YOUR SKIN

1. **Using a mirror; It's best to start on the forehead:** Total actual treatment time is about 3 to 5 minutes.
 - Do vertical passes up & down. Keep the vial upside down (*tip down*) when treating your skin.
 - Push down firmly into the skin, *but not too hard and not too fast - do a slow to medium speed.*
 - **Note: Serum does not flow out excessively.**
 - **Roll up & down doing vertical passes. While rolling; push the vial 'up & down' into the skin to help pump out the serum, this allows the serum to flow out more easily.**
 - One pass is all you need...if you have a problem area, you can do another pass over it.
 - Do the entire right side of your face down to the jawline, then do the left side of the face.
 - **Be very careful around your eyes;** a slower, more controlled pace is recommended around the eye area. You can do the upper brow area. **DO NOT TREAT EYELIDS!**
 - You can treat your lips. Be sure to treat your Neck and Jawline as well.
 - You can treat your upper Décolleté as well, if so desired.
 - Once your treatment is completed; **you will have some serum left in the vial:**
 - Put the clear cap back on the Roller tip and carefully pull (remove) the Roller tip out of the bottle by wiggling it back and forth and dispose of it properly - **Do Not Flush.**
 - **Pour some of the remaining serum into your hand and gently apply the serum evenly to your face & neck. Do not press too hard. No need to use all the Serum.**
 - If you still have some serum is still left, put the gray rubber cap back into the bottle. Apply any remaining serum **12 hours later** on a clean face & neck. And again 12 hours later, if there is any left.
 - Properly dispose of the glass vial & microneedle in the trash when done – **Do Not Flush.** The vial and Roller tip are designed for one-use only. **DO NOT REUSE TIP OR INFUSION VIAL.** *Bacteria growth would cause treatment issues if re-used and could cause serious skin issues.*
BE SURE CHILDREN OR PETS CAN NOT GET ACCESS TO MICRONEEDLE TIP OR INFUSION VIAL; THEY CAN BE A CHOKING HAZARD.

2. **Only perform one treatment every two weeks. Do not treat every week.** Your Body will begin Collagen production immediately and your skin will improve over time. **Perform your second treatment two weeks later: Collagen growth peaks at Day 14, so performing another treatment on Day 14 will help maintain that Collagen peak, producing maximum results! For best results; use the Beautify Million Dollar Growth Factor Serum twice a day on face & neck – Save \$111! MillionDollarSerum.com (QR Code below)**

POST-CARE INSTRUCTIONS:

- **Do not touch the treated area as it may be more prone to infection; this includes children & kissing pets.** It's important to keep it as clean as possible. If your skin itches; lightly tap the itch with a new Q-tip. **Change your pillow case after treatment - Sleep on a clean pillow case.**
- **Do Not Wash Your Face after your Micro-Infusion treatment** before going to sleep. The following day after your treatment, you may gently cleanse the treated area. Be sure to clean (sanitize) your hands prior. **DO NOT SCRUB YOUR FACE.**
- **Do not apply ANY make-up cosmetics, or skin care products for a minimum of 24 hours** following your Micro-Infusion treatment. We recommend that you do not use any make-up cosmetics on the treated area for 48 hours after your Micro-Infusion treatment as make-up may contain bacteria. Clean your make-up brushes with Alcohol. It is best if you use Mineral make-up.
- **DO NOT lay in direct sunlight for at least one (1) week.** If you are going to be exposed to the sun, wear a hat & apply an SPF 30 or higher Natural Sun Screen starting on the 2nd day, after your Micro-Infusion treatment. Additionally, it is important to avoid Tanning Bed exposure for at least one week post treatment. *For lasting results, reduce the amount of direct sun exposure and tanning beds as part of your healthy skin care regimen.*
- **Do not use any Glycolic/Lactic Acid or Retin-A Products for a minimum of three (3) days** after Micro-Infusion Treatment.
- **Do not sweat heavily, workout, swim, use a steam room or sauna, or do hot yoga for one day** after Micro-Infusion.
- **Do not get a Massage, a Chiropractic adjustment, do Physical Therapy or go to the Dentist for a minimum of one day after your Treatment.** It's best if you do not travel for 24 hours after treatment. Wait 3 days before doing an LED Skin Rejuvenation Treatment.
- **Do not or have any other facial therapy performed for a minimum of (7) seven days. Wait two (2) weeks to get Botox or Injectables after your Micro-Infusion Treatment.**



You Can Purchase Beautify® World's Best Natural Sunscreen Here:



Our Sunscreen Anti-Aging Formula Contains 7 Anti-Oxidants to Fight Free-Radicals:

It even contains a Copper Peptide, which helps grow Collagen! Contains Super Anti-Oxidant Grapeseed and Vitamins B2, B5, C, E & Anti-Oxidant Green Tea. Also contains Beautifying & Skin-Soothing Aloe-Vera, Shea Butter, Jojoba, Seaweed, Irish Moss & Chamomile. Even Contains Hyaluronic Acid! Uses 19% Zinc Oxide.



**You Can Purchase Beautify® Million Dollar Growth Factor
1 oz. Daily Serum Here: Save \$111!**



IMPORTANT SAFETY & USAGE GUIDELINES

This guide outlines essential information for the safe and effective use of the Beautify® Micro-Infusion Device. Users are responsible for thoroughly reviewing all content in this guide before using the Beautify® Micro-Infusion System.

Always follow the Pre-care & Post Care guidelines provided to ensure a positive treatment experience. It is up to the user or buyer to be fully informed about any local regulations and restrictions regarding micro-infusion treatments – *if any*.

Be sure to read this guide carefully and completely. Keep this guide in a place where it's easy to reference whenever needed. Watch out treatment video for best results.

DISCLAIMER

For a safe and successful treatment; please take time to read and understand the disclaimers, usage & care instructions prior to using the Beautify® Micro-Infusion System. Beautify® and its Affiliates are not responsible or liable for the treatment outcomes, as results may vary from person to person.

By using the Beautify® Micro-Infusion System, you are consenting to our pre-care and post-care Instructions. *This treatment is cosmetic only and not designed to diagnose or treat medical conditions* and should be strictly used according to the instructions provided.

20 Micro-Needles vs. 64 Micro-Needles

Are you sick & tired of buying skin care products & devices that don't work as promised? To truly affect positive skin care improvement, you need more than 20 micro-needles, plus you need a serum with the finest ingredients all at the highest potency! Our Beautify® Roller micro-needle system has over 3 times the micro-needles with 64 micro-needles. This creates a much greater number of micro-channels allowing for more product to be infused into your skin at up to 900% better absorption! Plus, by rolling it's way faster and much easier! Coupled with one of the finest serums in the world; you see better results! Rolling is much better than out dated 'stamping' with only 20 needles, which is considered outdated & archaic. *Think of a pogo stick vs. a bicycle!* Our system has more than 3 times the micro-needles delivering 900% absorption, yielding significantly better results!

Stamping vs. Micro-Roller

Our Beautify® Micro-Infusion Roller System is much better than stamping, which only has 20 needles. Our Beautify® Micro-Roller has 64 micro-needles that is much easier and way faster to better create micro-channels to infuse our Clinically-Formulated Million Dollar Growth Factor Serum into the skin's surface. Think of a pogo stick vs. a bicycle! Our system nourishes your skin and amplifies absorption by over 900%, resulting in noticeably visible improvements. The results; more youthful and beautiful looking skin! Stamping with only 20 needles is harder to perform, much slower, and considered to be a dinosaur; archaic.

Virtually Pain-Free

Our micro-needles are finer than human hair, and provide a virtually pain-free treatment! *When we think of needles, we automatically think; 'PAIN!'* But that is simply not the case...our micro-needles are only 0.5mm long, which is 1/64" long. **The nerves, blood vessels, oil glands are down 0.6mm to 4mm is the dermis layer, so our 0.5mm micro-needles cannot reach them.** This allows for a virtually pain-free treatment, yet it creates the micro-channels to infuse our Million Dollar Growth Factor serum directly into your skin at a 900% better absorption rate for maximum results! Treatment can be a bit tingly or spicy. It will not produce any droplets of blood. *If you have severe allergies, sometimes it can be a bit painful.*

FAQ's - Frequently Asked Questions

When I will see the results?

You should notice some results fairly quickly—After your first treatment, you may notice your skin is a bit tighter and firmer often within a day or two. It may also look brighter, smoother, and more refreshed.

Your own Body will begin Collagen production will begin immediately and your skin will improve over time. Perform your second treatment two weeks later! **Collagen growth peaks at Day 14, so performing another treatment on Day 14 will help maintain that Collagen peak, producing maximum results!**

How long will my results last?

The more treatments you do, the better the results will be. Remember, you're fighting the numbers of years you are in age...so more is more. However, only do a treatment every two weeks, *not every week*; let your Body regenerate. Typically, results can last around up to 6 months, depending on how many treatments you do. Then do a treatment every 6 to 8 weeks for maintenance.

To improve your results and make the results last even longer; use our Beautify® Million Dollar Growth Factor Serum every day, twice a Day – Morning & Night. It only takes about 6 seconds to apply! **Save 62% Off!** **(Order link on Page 8)**

Does this treatment hurt?

Our micro-needles are finer than human hair, and provide a virtually pain-free treatment! *When we think of needles, we automatically think; 'PAIN!' But that is simply not the case...*our micro-needles are only 0.5mm long, which is 1/64" long. **The nerves, blood vessels, oil glands are down 0.6mm to 4mm is the dermis layer, so our 0.5mm micro-needles cannot reach them.** This allows for a virtually pain-free treatment, yet it creates the micro-channels to infuse our Million Dollar Growth Factor serum directly into your skin at a 900% better absorption rate for maximum results!

Is this treatment safe?

Yes, it is quite safe, provided you are responsible and follow our guidelines and instructions. Our system only creates micro-channels using 0.5mm micro needles, and even though they are virtually invisible; however the skin is open and do close up for up to 24 hours. It is your responsibility to follow the pre-care and post-care instructions. *Beautify cannot be responsible for issues.* This is a simple, easy at-home treatment; however, it is of paramount importance that you follow the pre-care & post-care instructions.

Is using a Roller bad for your skin bad?

No, not with our Beautify® Micro-Infusion Roller. **Derma-Rollers, which are much larger, are considered bad for your skin because they create a 'V' in the skin.** This occurs with a larger type Roller because it pushes the needles into the skin at an angle, then the Roller keeps turning and comes out at the opposite angle creating a 'V', *which is too invasive and damaging for the skin.* Our Beautify® Micro-Infusion Roller is very small, so our micro-needles don't get pushed into the skin at a severe angle like large Derma-Rollers.

How often can I do a treatment?

Only do a treatment every two weeks, Do NOT do a treatment every week; Allow your Body regenerate for optimal results. **Collagen growth peaks at Day 14, so performing another treatment on Day 14 will help maintain that Collagen peak, producing maximum results!**

Keep doing treatments every two weeks until desired results are obtained. Then do one treatment every 6 to 8 weeks for maintenance to keep your results. For best results; use the Beautify Million Dollar Growth Factor Serum twice a day on face & neck.

Can I reuse the Mini-Roller Tip?

No, DO NOT REUSE TIP. Roller tip is designed for one-use only. *Bacteria growth could cause serious skin issues if re-used.* **Keep Roller Tip away from Children & Pets. Small parts may present a choking risk. Dispose of properly - Do Not Flush.**

When is the best time to do micro-infusion?

We recommend doing the treatment in the late afternoon or evening. Try to avoid the sun and skip makeup for 24 hours afterward to let your skin recover fully. **BE SURE TO SLEEP USING A CLEAN PILLOW CASE.**

Can you use it around the eyes, lips and the neck?

Eye Area: Yes. Just be very gentle and use light pressure when rolling around the under-eye area as the skin is thinner and more sensitive. Only do 1 pass to start, and watch your skin closely since everyone's skin can react differently. You can do the upper Brow area. **DO NOT ROLL YOUR EYELIDS. DO NOT ROLL TOO CLOSE TO YOUR EYES.**

Neck Area: Use our Beautify® Micro-Infusion Growth Factor system on your neck area to help reduce fine lines and rejuvenate the skin. It's great for tackling neck creases, horizontal lines, and signs of aging, making the skin look smoother, tighter & firmer, and more youthful. *The Neck area is often neglected and may show your age, therefore, it may need more treatments.*

Lips: Yes, however, do not press too hard and only do one pass. It makes lips appear plumper, fuller and healthier, hydrates dry lips, reduces fine lines and wrinkles, Promote natural, long-lasting fullness. *Note: Lip results do not last as long due to constant movement of the lips.*

Can I use any serum every day after treatment?

We would *not* recommend it. So many skin care products have harmful, dangerous ingredients including proven carcinogens and endocrine disruptors, which can cause cancer and disrupt your hormones. Further, no serum has the power to improve skin like Beautify Million Dollar Growth Factor Serum. It contains 5 Clinically-Proven Growth Factors + 2 Super Peptides Formulated at the Highest Potency Possible! Plus Collagen and 2 Hyaluronic Acids, and Special Skin Care Vitamins & Minerals to Nourish Skin! NO Dangerous or Harmful Ingredients! Takes Just 6 Seconds – Twice A Day To Realize Your Best Skin! **Save up to 62% Off! (Order link on Page 8)**

What are the benefits of Beautify® Micro-Infusion?

Our 7 Super Ingredients Have Clinically-Proven Effectiveness And Work Synergistically Together and are infused into the skin with 900% absorbability...Producing Beautiful, More Youthful Looking Skin!

- **Reduces Lines & Wrinkles:** Clinically-Proven to Stimulate Collagen & Hyaluronic Acid
- **Anti-Oxidant Protection:** Fights Free Radicals; The #1 Cause of Skin Issues!
- **Firms Sagging Skin:** Seven Clinically-Proven Ingredients to Tighten & Firm Skin
- **Increases Cell Turnover:** Makes Skin Look More Radiant, Brighter & Smoother
- **Improved Skin Texture:** Visibly Refined Appearance, Produces Glowing Skin
- **Reduced Pore Size:** Creates a Smoother, More Even Looking Complexion
- **Reverses Sun Damage & Pigmentation:** Helps Fade Brown Age Spots
- **Reduces Inflammation:** Revealing Tighter, Firmer, More Sculpted Looking Skin
- **Reduces Skin Roughness:** Makes Skin Softer & More Supple
- **Improved Skin Hydration:** Plumps Skin & Gives Better Elasticity
- **Enhances Blood Vessel Generation:** Produces Healthier, Radiant, Fresh-Looking Skin
- **Strengthens Skin Barrier:** Protects Skin from Harmful Toxins & Substances
- **Total Skin Rejuvenation:** Produces a Beautiful Complexion!

How many treatments do I get out of one Micro-Infusion Kit?

You get two treatments with our 2-bottle kit, and it contains two Micro-Roller tips in sterile packs. (One for each treatment). It's considered a 30-Day kit, because you should do your first treatment on Day One, and your second treatment on Day 15, and then you will see the best results by Day 30.

Do I need the Beautify® Daily Anti-Aging Serum?

Yes, if you want the best optimal results! This is one of the finest serums on the planet and contains the same exact strength & formula used in your Beautify® Micro-Infusion kit! **Save up to 62% Off! (Order link on Page 8)**

Our 'First-of-its-Kind' Anti-Aging Rejuvenation System Contains 5 Clinically-Proven Growth Factors + 2 Super Peptides Formulated at the Highest Potency Possible! Plus Special Skin Care Vitamins & Minerals! NO Dangerous or Harmful Ingredients! Just 6 Seconds, Twice A Day, To Realize Your Best Skin!

I have some leftover serum after my treatment, what should I do with it?

Once you've completed your treatment, take any leftover serum from the vial and pat it onto your face and neck right after you finish your treatment. Leave it on overnight to let your skin infuse into your open microchannels and optimal results! If you have any serum still left, then apply it again all over your clean face and neck 12 hours after your treatment. One bottle is meant for a single treatment session. Dispose of it properly & Keep Out of the Reach of Children or Pets. *Do not flush down the toilet.*

What is size of the needles?

Our micro-needles are finer than human hair, and provide a virtually pain-free treatment! *When we think of needles, we automatically think; 'PAIN!'* But that is simply not the case...**our micro-needles are only 0.5mm long and our nerves, blood vessels, oil glands are down 0.6mm to 4mm in the dermis layer, so our 0.5mm micro-needles cannot reach them.** This allows for a virtually pain-free treatment.

What's type of micro-needling is better? Stamping or Rolling?

Our Beautify® Micro-Infusion Roller System is much better than stamping, because stamping only has 20 needles. Our Beautify® Micro-Roller has 64 micro-needles that is much easier and way faster to better create micro-channels to infuse our Clinically-Formulated Million Dollar Growth Factor Serum into the skin's surface. Think of a pogo stick vs. a bicycle! Our system nourishes your skin and amplifies absorption by over 900%, resulting in noticeably visible improvements. The results; more youthful and beautiful looking skin! Stamping with only 20 needles is harder to perform, much slower, and considered to be a dinosaur; archaic.

How does the 60-day risk-free money back guarantee work?

You Try it Risk-Free: 60 Day Money Back Guarantee on All Orders! We give you a full 60 days because we know it works, and we want to give you enough time to see it work. Your best results will be realized a full 30 days *after your second treatment because of Collagen growth*, so this gives a chance for it to see some skin improvement. *Remember, you're fighting the number of years you are in age, so you may need more than 2 treatments to realize better results*, but you should see some results after 2 treatments.

If you didn't see results and want a refund, you must submit a Before pic from day one and an After pic from day 45 to 60 and provided you have performed 2 treatments within a 20 day time period (15 days is optimal – meaning Treatment on Day 1 and a 2nd treatment on day 15). This will allow us to see what you're seeing, and if approved, you will be contacted and asked to mail back your used bottles and we will give you a 100% full refund. 60 days starts from the date of delivery.

Why isn't the serum flowing out of the device when I'm Rolling?

Take a paper towel and place it on the counter. Turn the serum bottle upside down and gently press the micro-roller tip down on the paper towel. *(Do not press down too hard as it may bend the micro-needles).* The Micro-needle tip has a spring in it should allow some air into the serum bottle, *(You should see a bubble arise thru the serum showing it is allowing air into the bottle)*, which should help prime it and start the flow of the serum when you gently press the roller tip on your skin and roll. **During treatment; push the vial 'up & down' to pump the serum; this allows it to flow out more easily.**

Important Safety & Usage Guidelines

PRECAUTIONS – PLEASE READ

The Beautify® Micro-Infusion System is designed for personal use only.

Keep the device and its components away from Children & Pets

Small parts may present a choking risk. Dispose of properly – Do Not Flush.

DO NOT REUSE MICRO-NEEDLE TIP. Roller tip is designed for one-use only.

Bacteria growth could cause serious skin issues if re-used.

If you experience any discomfort or adverse effects, discontinue use immediately and consult your healthcare provider. Only use skincare products that are fresh and uncontaminated after treatment. A patch test is recommended prior to starting a full treatment session.

- **Consult a Professional:** If you're unsure whether the Beautify® Micro-Infusion System is right for you, seek advice from a healthcare professional or dermatologist before starting treatment.
- **Clean Skin First:** Make sure your skin is thoroughly cleansed and free from any makeup, oils, or impurities before beginning your micro-infusion session.
- **Pigmentation Caution:** Use extra care if you have a history of post-inflammatory hyperpigmentation (PIH), which may cause darkening of the skin. Those with PIH are more prone to pigmentation changes following micro-infusion treatments. Stay out of the sun for at least one week. Only wear natural sun block - no chemical sun block!
- **Product Compatibility:** Avoid applying products to your skin during or after treatment that are not specified for use with the Beautify® Micro-Infusion System. Only use Beautify® Million Dollar Growth Factor Serum.
- **Protect Treated Skin:** Safeguard newly treated skin from harsh environments, including sun exposure, UV rays, and any physical contact that could cause irritation. Wear only natural sun block...no chemical sun block.
- **Targeted Use Only:** Only use the Beautify® Micro-Infusion device on the facial areas specified in this guide. Avoid using it on sensitive zones or areas with active skin conditions.
- **Treatment Frequency:** Limit micro-infusion sessions to a maximum of one time per every two weeks to prevent skin sensitivity or over-treatment.
- **Needle Head Hygiene:** Always ensure that needle heads do not touch dirty surfaces or come into contact with unclean objects, as this can introduce bacteria and cause skin issues.
- **Proper Disposal:** Dispose of used the Micro-needle Tip and Infusion Vial responsibly after use. Do not reuse or attempt to sterilize them for future use. Do Not Flush.
- **Keep of the Reach of Children and Pets:** Do not let children play with the Micro-needle Tip and Infusion Vial. Do not leave them laying around for children to get access to or you pets as well – these items can cause a choking hazard.

Micro-Needle Tip Use Precautions – PLEASE READ

- **Needle Head Safety:** The Beautify® Micro-Infusion System utilizes specialized needle heads that, if improperly handled, may pose a risk of injury. Using contaminated or damaged needle heads can lead to serious health risks, including infection or transmission of diseases.
- **Check for Damage:** Exercise caution when managing needle heads, as any damage to the needles can cause injury. If a needle head appears bent, broken, or incomplete, stop use immediately and replace it with a new one to avoid potential harm.
- **Built-in Safety Mechanism:** The Beautify® Micro-Infusion System is designed with a built-in safety mechanism that prevents needles from penetrating deeper than 0.5 mm, ensuring a safe and controlled application.
- **Do Not Reuse:** Never attempt to reuse a needle head. Reusing needle cartridges can compromise safety and reduce the effectiveness of the treatment. Roller tip is designed for one-use only. ***Bacteria growth could cause serious skin issues if re-used.***
- **Proper Attachment & Removal:** Only attach or detach the needle head when the protective cover is securely in place. Ignoring this precaution can result in accidental injury.
- **Avoid Direct Contact:** Avoid direct contact with the needles at all times. Touching the needle tips can introduce contaminants or cause injury.
- **Use Only Approved Parts:** Only use Beautify® certified needle heads with the Micro-Infusion System. *Using unapproved needles may damage the device or result in improper treatment.*
- **Keep Away From Children & Pets:** Keep this device and its components away from children as small parts may present a choking risk.
- **Proper Disposal:** Dispose of used or contaminated needle heads immediately after each use. Proper disposal in household waste is recommended to prevent accidental injury. Do Not Flush.

DISCLAIMER

For a safe and successful treatment; please take time to read and understand the disclaimers, usage & care instructions prior to using the Beautify® Micro-Infusion System. Beautify® and its Affiliates are not responsible or liable for the outcomes of treatments as results may vary from person to person.

By using the Beautify® Micro-Infusion System, you are consenting to our pre-care and post-care Instructions. This treatment is cosmetic only and is not designed to diagnose or treat medical conditions and should be strictly used according to the instructions provided.

Store this guide in a place where it's easy to reference when needed.

Your Best Skin Now.

Give yourself a 'Med Spa quality' treatment
in Only 5 Minutes Without the 'High Med Spa prices'!
Utilize the same technology used
by Aestheticians & Dermatologists
to improve your overall skin appearance, quickly & easily.

No Appointments - No Waiting - No High Cost - No Hassles.

MillionDollarSerum.com

Beautify® Skin Care 2025
Beautify I AM LLC