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NAIPO

USER MANUAL

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Contents

English-----01-06

Deutsch-----07-12

Français-----13-18

Español-----19-24

Italiano-----25-30

Aukley International Limited

Web: www.naipocare.com

Email: wholesale@naipocare.com

Add: No.102, Building P09, China South City, Pinghu Town, Longgang Distrct, Shenzhen, China

CE

RoHS

Made in China

Neck & Back Massager with Heat and Vibration

Model: MGBK-C2030

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Introduction

Thank you for purchasing the Naipo Neck & Back Massager with Heat and Vibration. Please read this user manual carefully and keep it in a safe place for future reference. If you need any assistance, please contact our support team with your product model number and Amazon order number.

Included

1 x Neck & Back massager

1 x User manual

1 x Adapter

1 x Car adapter

Product Specifications

Model	MGBK-C2030
Name	Naipo Neck & Back Massager with Heat and Vibration
Input	DC 12V=4A
Power consumption	48 W
Default timer	15 minutes

Important Safeguards

Children, pregnant women, and the old should not use the product.

People with osteoporosis, spinal fracture, sprain, pulled muscles or other symptoms should not use the product.

While the product is not in use, please disconnect it from the power supply.

Use the product only with the matched adapter or car adapter.

Please do not use the product in a humid environment (like a bathroom) to avoid electric shock or other safety hazards.

If abnormal or severe pain is felt during the use, please stop using the product and contact our official customer service promptly.

Please do not use the product if damage is found on its outer materials or the power cable.

- 01 -

- 02 -

①b.

②a.

②c.

①a.

①b.

③c.

③a.

③b.

④c.

④a.

④b.

⑤

①a.

②a.

②c.

①a.

①b.

③c.

③a.

③b.

④c.

④a.

④b.

⑤

①a.

②a.

②c.

①a.

①b.

③c.

③a.

③b.

④c.

④a.

④b.

⑤

① Neck

a. The "Back Kneading" button: Press this button, and the massage nodes massage in one direction (indicator light: flickering). Press this button again to activate medium-level vibration function (indicator light: flickering rapidly). Press this button for the third time, and the massage would stop.

b. The "Up" button: Long press this button to move the neck massage nodes up.

c. The "Down" button: Long press this button to move the neck massage nodes down.

② Back

a. The "Full" button: Press this button, and the back massage nodes would massage up and down on the whole back (indicator light: on). Press this button again to stop the massage (indicator light: off).

b. The "Upper" button: Press this button, and the back massage nodes would massage up and down on the upper back (indicator light: on). Press this button again to stop the massage (indicator light: off).

c. The "Lower" button: Press this button, and the back massage nodes would massage up and down on the lower back (indicator light: on). Press this button again to stop the massage (indicator light: off).

③ Massage on Selected Area

a. The "▲" button: Long press this button, and the back massage nodes would move up (indicator light: flickering). When reaching the top, the back massage nodes will stop moving. After loosening the press (indicator light: on), the product will begin massage on selected areas: the back massage nodes would massage on the area between where you press this button and where you loosen the press.

b. The "▼" button: Long press this button, and the back massage nodes would move down (indicator light: flickering). When reaching the bottom, the back massage nodes will stop moving. After loosening the press (indicator light: on), the product will begin massage on selected areas: the back massage nodes would massage on the area between where you press this button and where you loosen the press.

c. The "Reset" button: Press this button to activate set massage mode (indicator light: on), and the massage nodes would massage up and down on the certain part in a small range. Press this button again to stop set massage mode (indicator light: off).

- 03 -

① Hip Massage

The "Back Kneading" button: Press this button to activate the vibration function with low-level by default (indicator light: flickering slowly). Press this button again to activate medium-level vibration function (indicator light: flickering rapidly). Press this button for the third time to activate high-level vibration function (indicator light: on). Press this button for the fourth time to turn off the function (indicator light: off).

② Hip Massage

The "Back Kneading" button: Press this button to activate the vibration function with low-level by default (indicator light: flickering slowly). Press this button again to activate medium-level vibration function (indicator light: flickering rapidly). Press this button for the third time to activate high-level vibration function (indicator light: on). Press this button for the fourth time to turn off the function (indicator light: off).

- 04 -

Questions & Answers

Problem	Solution
The product could not be turned on.	Please check whether the adapter has been connected to the product properly. Please check whether the adapter has been connected to the power supply properly. The product would shut off automatically after 15 minutes of continuous use.
The product stops working during the use.	The product is overloaded. Please use after disconnecting the back. The built-in overheat protection has been activated. Please wait for at least 30 minutes before using the product again.
The height is not proper.	Please adjust the height of neck massage nodes following the product Operation Guide.
The product does not work satisfactorily on the driver seat.	It's recommended to place the product on the back seat for use.
Pain is felt during the use.	Please avoid direct contact of your skin with the product.
The heating function does not work.	After activating the heating function, it needs 2-3 minutes to warm up.

- 05 -

Product Maintenance

1) Please disconnect the product from the power supply before cleaning.  
2) Wipe the dirt on the controller with a dry towel.  
3) Please wipe the surface of the PU leather with a mild damp sponge or dry cloth.  
4) After using, please disconnect the product from the power supply, and place it in a dry and airy place for storage.

Warranty & Customer Support

If you need help regarding product issues, after-sales service, or other relevant problems, please contact us via the email on the cover.

\*Note: Products purchased directly from Naipo benefit from our after-sales service. If the purchases are made from other selling teams, we recommend you contact them for warranty issues.

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