

材质：70g书纸
印刷四色
尺寸：85*140mm
工艺：风琴折

出血线
折叠线
裁切线

<div><div>Welcome to use our Breo Laser Hair Growth Comb</div><div><div><div><div></div><div></div></div><div><div></div><div></div></div></div><div><p>This product is a laser instrument. Users are required to read the manual carefully before operation to ensure efficient and safe use of all the functions and maintenance of the product.</p><p>Warnings:</p><ol style="list-style-type: none">Do not make direct eye contact with the laser light diodes or its reflection in a mirror because it can temporarily irritate your eyes! you stare at the light for too long, it could harm your eyes. Always wear eye protection when using the product because it can cause temporary irritation to your eyes.Keep out of reach of children at all times.If the component or structure is damaged of the device, do not use the device because you could get an electrical shock. If the component or structure gets damaged, contact our company.This is an ESD sensitive device, do not touch charging contacts.If you are taking any medications or products that make you sensitive to light, you should test the device to make sure you do not have a reaction. Some of the things that can cause increased sensitivity to the device could be medicines for colds, allergies, pain medicines and drugs to treat infections. To test the device sensitivity, please refer to "Section Testing the Laser Comb for Sensitivity".Handle the device with care, keep away from heat and liquids because this could cause the device not to work properly and may cause an electrical shock.Ensure to keep the light emitting ports clean. Otherwise, the laser light intensity output will be affected, adversely affecting the therapeutic effects of the comb.Do not charge the comb while it is being used.Turn off the comb while cleaning it or after use.Only use the included charger to charge the comb.If the comb is not used for a long period of time, charge it at least every 4-6 months to preserve its battery life.Only use the comb as directed. Do not try to change the way the comb operates. Changes may result in harmful laser exposure.</div><div><div>Version V1.0 Date 2024-6-30</div></div></div></div>	<div><div>Table of Contents</div><div><div>I. Working Principle</div><div>II. Performance Parameter</div><div>III. Package Content</div><div>IV. Structure and Components</div><div>V. Indications and Contraindications</div><div>VI. Operation</div><div>VII. Testing the Breo Laser Hair Growth Comb for Sensitivity</div><div>VIII. Product Maintenance</div><div>IX. Environmental Statement</div><div>X. Electromagnetic Compatibility</div><div>XI. Product Service</div><div>XII. Symbol Explanation</div><div>XIII. Contact Us</div></div></div>	<div><div>I. Working Principle</div><p>Breo Laser Hair Growth Comb works by providing laser energy to stimulate hair follicles. For optimal results, the laser must not be blocked by the hair and must have an unobstructed path to the scalp. Breo Laser Hair Growth Comb has teeth that part the hair and allow the maximum amount of laser energy to reach the scalp.</p></div>	<div><div>II. Performance Parameter</div><div><div>Technical Information</div><div>Trade name/Common name: Breo Laser Hair Growth Comb</div><div>Laser medium: GaAlAs Semiconductor</div><div>Laser wavelength: 650nm ± 10nm</div><div>Power: 15 mW</div><div>Laser diodes quantity: 7PCS</div><div>Input voltage: 5Vdc</div><div>Input power: less than 5V1A</div><div>Working environment:</div><div>● Temperature: 5℃~30℃</div><div>● Relative humidity: 15%~80%RH</div><div>● Atmospheric pressure: 860hpa~1060hpa</div><div>Battery capacity: 2000mAh 3.7V</div></div></div>	<div><div>III. Package Content</div><div><div><div><div></div><div>Breo Laser Hair Growth Comb</div></div><div><div></div><div>Storage Bag</div></div><div><div></div><div>User Manual</div></div><div><div></div><div>Charging Cable</div></div></div></div></div>	<div><div>IV. Structure and Components</div><div><div><div>Main Structure</div><div>The treatment instrument consists of host, comb cover, power charger jack. Massage head is applied part</div></div><div><div><div>Massage head</div><div>Laser light</div><div>Power/ Mode button</div></div><div><div></div><div>charging Cable</div></div></div></div></div>	<div><div>V. Indications and Contraindications</div><div><div>1. Indications for use</div><div>Breo Laser Hair Growth Comb is indicated to treat: Androgenetic Alopecia, and promote hair growth in females who have Ludwig's(Savin) I-4, II-2, or frontal patterns of hair loss and in males who have Norwood Hamilton Classifications of 1 to 5 V and who both have Fitzpatrick Skin Types I to IV.</div></div><div><div>2. Contraindications</div><div>Pregnancy, cancer and hemorrhagic diseases are prohibited to use this product. Forbidden for use by individuals with pacemakers since the Breo Laser Hair Growth Comb may interfere with the operation of the pacemaker. If you are taking medications or if you have any medical condition which makes your skin sensitive to sun or light or if you suffer from scalp injury, the product should not be used.</div></div><div><div>It should not be used by or on children under 18 years of age.</div><div>Before operation, please note the following:</div><div>> The System should never be connected directly to the wall charger and only be powered by the adapter.</div><div>> How do you know if the Breo Laser Hair Growth Comb is regulator you?</div><div>Breo Laser Hair Growth Comb can be used by men and women with thinning hair (including color treated hair) or pattern baldness caused by a hereditary condition. Doctors use a system known as the Norwood Hamilton Classification (men) and the Ludwig-Savin Classification (women) to describe the degree of hair loss. Below are the pictures of the scales. The shaded areas show the type of hair loss that can be treated with the Breo Laser Hair Growth Comb.</div></div></div>	<div><div>V. Indications and Contraindications</div><div><div>1. Indications for use</div><div>Breo Laser Hair Growth Comb is indicated to treat: Androgenetic Alopecia, and promote hair growth in females who have Ludwig's(Savin) I-4, II-2, or frontal patterns of hair loss and in males who have Norwood Hamilton Classifications of 1 to 5 V and who both have Fitzpatrick Skin Types I to IV.</div></div><div><div>2. Contraindications</div><div>Pregnancy, cancer and hemorrhagic diseases are prohibited to use this product. Forbidden for use by individuals with pacemakers since the Breo Laser Hair Growth Comb may interfere with the operation of the pacemaker. If you are taking medications or if you have any medical condition which makes your skin sensitive to sun or light or if you suffer from scalp injury, the product should not be used.</div></div><div><div>It should not be used by or on children under 18 years of age.</div><div>Before operation, please note the following:</div><div>> The System should never be connected directly to the wall charger and only be powered by the adapter.</div><div>> How do you know if the Breo Laser Hair Growth Comb is regulator you?</div><div>Breo Laser Hair Growth Comb can be used by men and women with thinning hair (including color treated hair) or pattern baldness caused by a hereditary condition. Doctors use a system known as the Norwood Hamilton Classification (men) and the Ludwig-Savin Classification (women) to describe the degree of hair loss. Below are the pictures of the scales. The shaded areas show the type of hair loss that can be treated with the Breo Laser Hair Growth Comb.</div></div></div>	<div><div>VI. Operation</div><div><div>This product cannot be used while charging. After charging, please unplug the power adapter before use.</div><div><div>Power on/ Mode Selection</div><div>Press and hold the power/ mode button for about 2 seconds, and the prompt sound is long "dt", which means the product is turned on. After turning on, the default mode is mode one. After turning on the machine, short press the power button to switch the mode, and the prompt sound short "dt" means switching succeeded.</div><div>Mode 1: Slow massage</div><div>Mode 2: Laser light</div></div><div><div>Skin contact start/stop:</div><div>For safety, the laser switch limit is added. When the comb is within a safe distance from the scalp, the laser starts. When the safe distance is exceeded, the laser turns off.</div><div><div>Treatment frequency:</div><div>3 times per week, every time 15 mins, spaced out every other day, as little as 16 weeks.</div><div>When the ambient temperature is 20 °C, the intended use of the ME device can be achieved from the lowest storage temperature after use to the readiness of the ME device.</div><div>When the ambient temperature is 20 °C, the intended use of the ME device can be achieved from the maximum storage temperature to the ambient temperature of the ME device 30°.</div><div>The patient is expected to be the operator.</div></div></div></div></div>	<div><div>VII. Testing the Breo Laser Hair Growth Comb for Sensitivity</div><div><div>Testing the Breo Laser Hair Growth Comb for sensitivity before beginning treatment, be sure to test the Breo Laser Hair Growth Comb for sensitivity. Turn the device on, hold the Breo Laser Hair Growth Comb on one spot, shining the light on your forehead. After two minutes, turn the device off and set aside. Look at the spot where you held the laser light over the next five minutes. If you see any reaction such as redness or rash, do not use the Breo Laser Hair Growth Comb and contact Agent. If you don't experience any reaction, proceed with the treatment instructions.</div></div></div>	<div><div>VIII. Product Maintenance</div><div><div>Maintenance</div><div>(1) Wipe the dust off of the device's surface with a soft dry cloth.</div><div>(2) If stained with dirt or oil, use a soft brush to gently remove the stain; then wipe down with a dry soft cloth.</div><div>(3) To disinfect the device, moderately wet a cotton ball/square with disinfecting ethyl alcohol until moist. Use the cotton ball/square to wipe down the surface of the device and allow to air dry until the ethyl alcohol smell is gone.</div><div>(4) Keeps stored in a cool and dry area away from direct sunlight.</div><div>(5) Storage and transport environment:</div><div>a) Ambient temperature: 0℃~40℃</div><div>b) Relative humidity: 10%-80% RH</div><div>c) Atmospheric pressure: 860hpa~1060hpa</div><div>2) Cleaning your Breo Laser Hair Growth Comb</div><div>(1) Keep your Breo Laser Hair Growth Comb clean.</div><div>(2) Turn off the Breo Laser Hair Growth Comb by pressing and holding the power button (lower button) until the light turns off.</div><div>(3) Unplug the power cord from the device and from the electrical wall outlet.</div><div>(4) Clean the laser window with a soft, slightly moistened cloth.</div><div>(5) To clean the comb teeth, gently remove any hairs or dust with a slightly moistened towel.</div><div>(6) Don't use any acetone, soap or detergents, only warm water or ethyl alcohol.</div><div>(7) Visual inspection for residual soil and an additional repeat cleaning instruction in the event that residual soil is visible on the device.</div></div><div><div>Warning: Do not modify this equipment without authorization of the manufacturer.</div></div></div>	<div><div>VIII. Product Maintenance</div><div><div>Maintenance</div><div>(1) Wipe the dust off of the device's surface with a soft dry cloth.</div><div>(2) If stained with dirt or oil, use a soft brush to gently remove the stain; then wipe down with a dry soft cloth.</div><div>(3) To disinfect the device, moderately wet a cotton ball/square with disinfecting ethyl alcohol until moist. Use the cotton ball/square to wipe down the surface of the device and allow to air dry until the ethyl alcohol smell is gone.</div><div>(4) Keeps stored in a cool and dry area away from direct sunlight.</div><div>(5) Storage and transport environment:</div><div>a) Ambient temperature: 0℃~40℃</div><div>b) Relative humidity: 10%-80% RH</div><div>c) Atmospheric pressure: 860hpa~1060hpa</div><div>2) Cleaning your Breo Laser Hair Growth Comb</div><div>(1) Keep your Breo Laser Hair Growth Comb clean.</div><div>(2) Turn off the Breo Laser Hair Growth Comb by pressing and holding the power button (lower button) until the light turns off.</div><div>(3) Unplug the power cord from the device and from the electrical wall outlet.</div><div>(4) Clean the laser window with a soft, slightly moistened cloth.</div><div>(5) To clean the comb teeth, gently remove any hairs or dust with a slightly moistened towel.</div><div>(6) Don't use any acetone, soap or detergents, only warm water or ethyl alcohol.</div><div>(7) Visual inspection for residual soil and an additional repeat cleaning instruction in the event that residual soil is visible on the device.</div></div><div><div>Warning: Do not modify this equipment without authorization of the manufacturer.</div></div></div>
---	---	---	---	---	--	---	---	--	--	---	---

[illegible]