RENPHO

Model: RP-SNM061

Kneading Shoulder Massager

CONTENTS

ENGLISH	1-5
DEUTSCH	6-10
FRANÇAIS	11-15
ITALIANO	16-20
ESPAÑOL	21-25

Thank you for your purchase of RENPHO Shoulder & Neck Massager. The Massager offers a powerful shiatsu massage. While enjoying massage relief, you should use it with care and safety in mind. To ensure safety and avoid any property damage, please read this manual carefully and use it as indicated. Non-compliance may cause serious injuries or death and/or severe property damage.

IMPORTANT SAFETY INSTRUCTIONS

When using an electrical appliance, basic precautions should always be followed, including the following:
Read all instructions before using this appliance.

- For household use only, not for medical or other commercial purposes.
- Any individual who may be pregnant or has recently given birth, has a heart condition or a pacemaker, suffers from diabetes or other sensory impairment, phlebitis and/or thrombosis, joint dysfunction, inflammation or injuries, weak bones, osteoporosis, a high fever, an abnormal or cured sore, is at an increased risk of developing blood clots, has phisicocardiographic joints or other medical devices implanted in his/her body, is being treated for a medical condition or physical injury or is under a doctor's care should consult with a medical professional before using a massaging device designed for domestic use.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.

DANGER - To reduce the risk of electric shock:

- Always unplug the massager from the electrical outlet immediately after using and before cleaning.
- Do not use while bathing or in a shower.
- Do not reach for the massager that has fallen into water.
- Do not place or store appliance where it can fall or be pulled into a tub or sink.
- Do not place in or drop into water or other liquid.

WARNING - To reduce the risk of burns, fire, electric shock, or injury to persons:

- The massager should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
- Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- Close supervision is necessary when this massager is used by, on, or near children, invalids, or disabled persons.
- Use this massager only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this massager if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the massager to a service center for examination and repair.
- Do not carry this massager by supply cord or use cord as a handle.
- Keep the cord away from heated surfaces.
- Never use pins or other metallic fasteners with this appliance.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.

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Manufacturer: Shenzhen Ruiyi Business Technology Co., Ltd.
~~EU REP~~ C&E Connection E-Commerce (DE) GmbH
~~UK REP~~ Richers & Vision Ltd.