Myofascial Physiotherapy Device

Safety Instructions

User Manual

Danger: In order to reduce the risk of electric shock, fire and personal injury, this product must be used in accordance with the following instructions.

For adults only

It can only be used on dry, clean skin surfaces of the body, and cannot be used through clothes. Press and move it gently on the skin.

The duration on each part is about 60 seconds.

Use this device only on soft tissues of the body without causing pain or discomfort.

Do not use it on head or any hard or bony part of the body.

Only use the vibrating head that can produce the best results.

Bruise may be caused whether under moderate pressure or normal pressure.

Please always check the area that receives physiotherapy.

Stop using it immediately once there is any sign of pain or discomfort.

Keep fingers, hair or other parts of the body away from the shaft and back of the vibrating head, because they might be pinched.

Do not place any object in the vent of the device.

Do not immerse it in water or allow water to enter the vent of the device.

Do not drop the device or use it improperly.

Charge the device charger with only 8.4 VDC.

Carefully check the device and battery before each use.

Do not change or modify the device in any way.

Never operate or charge the device in an unattended way.

Instructions

1. Power on / off button (long press to start / long press to turn off)
2. Press the key gear to change
3. Detachable and replaceable vibrator, When disassembling the vibrator, just pull it out with force.
4. Please select 5V charging equipment for charging
5. For the first time, please charge for 3 hours. Please connect the data cable to the 5V adapter for charging.
6. Button type: when charging, 6 lights will flash, when full, 6 lights will be on, and when working, the lights will indicate the gear position.
7. Touch type: when charging, the nixie tube displays the electric quantity, while working, it displays the gear.

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