amazonbasics

Safety Manual

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Sicherheitshandbuch

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Veiligheidshandleiding

Handheld Percussion Massage Gun

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NLPagina 26

1. IMPORTANT SAFETY INSTRUCTIONS

Read these instructions carefully and retain them for future use. If this appliance is passed to a third party, then these instructions must be included.

When using electrical appliance, basic safety precautions should always be followed to reduce the risk of fire, electric shock, and/or injury to persons including the following:

Read all instructions before using this appliance.

⚠ DANGER

To reduce the risk of electric shock:

- ▶ Always disconnect this appliance from the electrical outlet immediately and switch off after using and before cleaning.

⚠ DANGER

- ▶ Do not reach for an appliance that has fallen into water. Disconnect immediately.
- ▶ Do not use while bathing or in a shower.
- ▶ Do not place or store appliance where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid.

⚠ WARNING

Risk of burn and damage!

- ▶ Never operate the appliance with the air openings blocked. Keep the air openings free of lint, hair, and the like.

⚠ WARNING

Risk of burn and damage!

- ▶ Never operate on a soft surface such as a bed or couch where the air openings may be blocked.
- ▶ Do not use massager in close proximity to loose clothing or jewelry.
- ▶ Keep long hair away from massager while in use.
- ▶ Never drop or insert any object into any opening.

⚠ WARNING

To reduce the risk of burns, fire, electric shock, or injury to persons:

- ▶ Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- ▶ Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- ▶ Do not use outdoors.
- ▶ Do not operate where aerosol (spray) appliances are being used or where oxygen is being administered.
- ▶ Close supervision is necessary when this appliance is used by, on, or near children, invalids, or disabled persons.
- ▶ Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the appliance to a service center for examination and repair.
- ▶ Keep the cord away from heated surfaces.

- ▶ For household use only.
- ▶ Do not carry this appliance by supply cord or use cord as a handle.
- ▶ To disconnect, turn all controls to the off position, then remove plug from outlet.
- ▶ For use only with Battery Pack: PT18650-3S 10.8V

- ▶ For use only with Charger: SHENZHEN SHI YINGYUAN ELECTRONIC CO.,LTD SAW20-150-1000XX Series.
- ▶ Do not use the appliance continuously for longer than 15 minutes to avoid overheating. Let the appliance cool down before using it again.
- ▶ The appliance is only to be charged with the supplied power adapter.
- ▶ This appliance should never be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls or who has sensory deficiencies in the lower half of their body.
- ▶ Do not use the appliance directly on:

a	Irritable skin conditions	g	Recent burn
b	Varicose veins	h	Recent surgery
c	Unconscious or sleeping person	i	Frostbite
d	Acute conditions requiring first aid or medical attention	j	Severe unstable hypertension
e	Deep vein thrombosis	k	Acute flare-up of inflammatory conditions such as rheumatoid arthritis
f	Sensitive skin, swollen or inflamed areas, in the presence of poor circulation, where skin eruptions are present, or in the presence of unexplained abdominal pain.		

- ▶ Consult your doctor prior to using the appliance if you are pregnant, if you have a pacemaker or if you have any concerns regarding your health.
- ▶ If you feel any discomfort whilst using this appliance, discontinue use and consult your doctor.

Battery Warnings

WARNING

- ▶ Read all safety warnings and instructions. Failure to follow the warnings and instructions may result in electric shock, fire and/or serious injury.
- ▶ Prevent unintentional starting. Ensure the switch is in the off-position before connecting to battery pack, picking up or carrying the appliance. Carrying the appliance with your finger on the switch or energizing appliance that have the switch on invites accidents.
- ▶ Recharge only with the charger specified by the manufacturer. A charger that is suitable for one type of battery pack may create a risk of fire when used with another battery pack.
- ▶ This appliance contains internal, non-replaceable lithium batteries. Use appliances only with specifically designated battery packs. Use of any other battery packs may create a risk of injury and fire.
- ▶ Do not use a battery pack or appliance that is damaged or modified. Damaged or modified batteries may exhibit unpredictable behavior resulting in fire, explosion or risk of injury.
- ▶ Do not expose a battery pack or appliance to fire or excessive temperature. Exposure to fire or temperature above 265°F may cause explosion.
- ▶ Follow all charging instructions and do not charge the battery pack or appliance outside of the temperature range specified in the instructions. Charging improperly or at temperatures outside of the specified range may damage the battery and increase the risk of fire.
- ▶ Have servicing performed by a qualified repair person using only identical replacement parts. This will ensure that the safety of the product is maintained.

- ▶ Do not modify or attempt to repair the appliance or the battery pack except as indicated in the instructions for use and care.
- ▶ This appliance contains batteries that are non-replaceable. Do not throw the appliance into a fire. Before disposal, the battery must be removed from the appliance. Contact qualified professional.
- ▶ Do not leave the battery charging for a prolonged period of time when the appliance is not in use.
- ▶ Do not charge the battery overnight.
- ▶ If the appliance appears to be bulging or shows other undesirable phenomena (e.g. excess noise), stop using the appliance immediately.
- ▶ Do not cover the appliance while charging the battery. Charge the appliance in a well-ventilated environment.
- ▶ When in long term storage, charge the battery to its 70% capacity once every 3 months to avoid shortening its life.
- ▶ Batteries produce heat when charging and discharging. If the appliance becomes too hot to touch, switch the appliance off and unplug immediately.
- ▶ If the battery leaks, avoid contact with skin and eyes. Rinse affected areas immediately with plenty of clean water, then consult a doctor.

SAVE THESE INSTRUCTIONS

Additional Warnings for Europe:

⚠ WARNING

- ▶ For the purposes of recharging the battery, only use the detachable supply unit provided with this appliance.
- ▶ This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- ▶ The supply cord cannot be replaced. If the cord is damaged the power adapter should be scrapped.

- The appliance is only to be used with the power supply unit provided with the appliance.
- The battery must be removed from the appliance before it is scrapped.
- The appliance must be disconnected from the supply mains when removing the battery.
- The battery is to be disposed of safely.

3. FCC – Supplier's Declaration of Conformity

Unique Identifier	B08L57T4P5- AmazonBasics Handheld Percussion Massage Gun - Professional Deep Tissue Muscle Massager 5 Speeds, 4 Interchangeable Attachment Heads
Responsible Party	Amazon.com Services LLC.
U.S. Contact Information	410 Terry Ave N. Seattle, WA 98109 USA
Telephone Number	(206) 266-1000

3.1 FCC Compliance Statement

1. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
 - (1) this device may not cause harmful interference, and
 - (2) this device must accept any interference received, including interference that may cause undesired operation.
2. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

3.2 FCC Interference Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio / TV technician for help.

4. Canada IC Notice

- This Class B digital apparatus complies with Canadian CAN ICES-3(B) / NMB-3(B) standard.

5. Disposal

The Waste Electrical and Electronic Equipment (WEEE) Directive aims to minimise the impact of electrical and electronic goods on the environment, by increasing re-use and recycling and by reducing the amount of WEEE going to landfill. The symbol on this appliance or its packaging signifies that this appliance must be disposed separately from ordinary household wastes at its end of life. Be aware that this is your responsibility to dispose of electronic equipment at recycling centers in order to conserve natural resources. Each country should have its collection centers for electrical and electronic equipment recycling. For information about your recycling drop off area, please contact your related electrical and electronic equipment waste management authority, your local city office, or your household waste disposal service.

6. Battery Disposal

Do not dispose of used batteries with your household waste. Take them to an appropriate disposal/collection site.

To learn more about recycling batteries, visit:
call2recycle.org/what-can-i-recycle

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