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iOS App

Android App

FitGo - Ideas for health
Health & fitness

ONE-CLICK SYNC WITH

Apple Health

Google Fit

MOBILE & DEVICE

iOS 7.0

Android 4.3

Bluetooth 4.0

FOR FIRST USE

1. BATTERY INSTALLATION

OPEN THE BATTERY DOOR, PUT INTO BATTERIES WITH CORRECT POLARITY.

2. APP INSTALLATION

SCAN THE QR CODE ON THE ABOVE, OR DOWNLOAD "FITGO" FROM APP STORE OR ANDROID GOOGLE PLAY.

3. REGISTRATION

TURN ON BLUETOOTH ON MOBILE DEVICE, RUN THE APP "FITGO", FOLLOW THE APP INSTRUCTIONS.

12 PARAMETERS IN 1 MEASUREMENT

WEIGHT

BMI

BODY FAT

WATER

MUSCLE

BONE

BMR

OBESITY

V-FAT

PROTEIN

AGE

LBM

MULTI USERS

8 USERS SHARE ONE USER ACCOUNT IN THE DEVICE.

TARGET WEIGHT

USERS CAN SET THEIR OWN TARGET WEIGHT TO CUSTOMIZE THEIR OWN HEALTH PLAN.

ATTENTION

1. DO NOT USE THE ACHATÉ SCALE ON PEOPLE WHO HAVE BODY IMPLANTS SUCH AS A PACEMAKER, ARTIFICIAL LIMBS, CONTRACEPTIVE DEVICES, OR METAL PLATES OR SCREWS. IT MAY CAUSE THE DEVICE TO MALFUNCTION OR PRODUCE AN INACCURATE RESULT. WHEN IN DOUBT, CONSULT YOUR PHYSICIAN.
2. PLEASE MAKE SURE GOOD CONTACT BETWEEN BARE FEET AND THE METAL ELECTRODES, OTHERWISE IT WILL DISPLAY YOUR WEIGHT ONLY.
3. ENSURE THE ACHATÉ SCALE IS ON A FLAT, EVEN AND FIRM SURFACE. AVOID CARPET OR SOFT SURFACES.
4. THE SURFACE OF THE ACHATÉ SCALE WILL BE SLIPPERY IF WET. PLEASE MAKE SURE BOTH THE GLASS AND YOUR FEET ARE DRY BEFORE USE.
5. THIS SCALE IS NOT WATER PROOF. DO NOT EXPOSE TO MOISTURE.
6. CLEAN THE ACHATÉ SCALE WITH A SLIGHTLY DAMP CLOTH. DO NOT USE ANY CHEMICAL/ABRASIVE CLEANING AGENTS TO CLEAN THE SCALE.
7. PLEASE REMOVE THE BATTERIES, IF THE ACHATÉ SCALE IS NOT TO BEING USED FOR A LONG TIME.
8. THE APPLIANCE IS INTENDED SOLELY FOR DOMESTIC USE WITHIN THE HOME.
9. CHILDREN SHOULD BE SUPERVISED TO ENSURE THAT THEY DO NOT PLAY WITH THE APPLIANCE.

Achaté