

COSLUS

ACTIVATE YOUR **24-MONTH** WARRANTY FOR FREE



* Please contact customer service if you're unable to provide the order number

www.coslus.com/warranty



Win Free Accessories



Covered by 24-Month Warranty



Friendly Responsive Customer Service

Expire time: 7 days after the delivery date

☎ +1 855 594 3950 (8am-5pm Mon. - Fri. (EST))


✉ support@coslus.com




FLIP →

QUICK GUIDE


Scan QR code to watch the operation video. >>>




Charge Instruction



Install nozzle



Fill water




! Slight bleeding is normal for first-time use or sensitive gums. Recommend using "Massage" mode for a week. If bleeding continues, consult a dentist for possible gum disease and follow their advice.



! For first-time or sensitive gums, use the periodontal tip with Massage mode for gum protection.



Start irrigating



! Align at 90° to the gum line, not on the gums, to prevent irritation.