HEATCONFORT 6000 THERMALPAD

Almohadilla eléctrica/ Electric heating pad

Manual de instrucciones
Instruction manual
Manuel d'instructions
Bedienungsanleitung
Manuale di istruzioni
Manual de instruções
Handleiding
Instrukcja obsługi
Návod k použití
Kullanma kılavuzu

Instrucciones de seguridad	4
Safety instructions	8
Instructions de sécurité	12
Sicherheitshinweise	16
Istruzioni di sicurezza	21
Instruções de segurança	25
Veiligheidsinstructies	29
Instrukcja bezpieczeństwa	33
Bezpečnostní pokyny	37
Güvenlik talimatları	41

ÍNDICE

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SAFETY INSTRUCTIONS

Read these instructions thoroughly before using the appliance. Keep this manual for future reference or new users as it contains important instructions.

- This appliance is designed for domestic use only and is not intended for bars, restaurants, farmhouses, hotels, motels, and offices.
- This appliance is not intended for use in hospitals.
- If the power cable is damaged, it must be replaced by the official Cecotec Technical Support Service or similar qualified personnel to avoid risks.
- This appliance can be used by children aged 8 years and above and people with reduced physical, sensory, or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children must not play with the appliance. Cleaning and user maintenance should not be carried out by unsupervised children.
- Children must be supervised to make sure they do not play with the appliance.
- Check the appliance regularly for signs of wear or damage. If there are such signs, if the appliance has been misused or is not working, return it before switching it on again.
- The appliance must not be used for heating animals.

- The appliance must be repaired or replaced if the cover is worn.
- Do not crumple the appliance by placing objects on top of it during storage.
- Do not operate the heating pad if the controller or the power cord is damaged or if you notice that the appliance is not working properly.
- Before connecting the appliance to the mains, make sure that the voltage indicated on the technical specification label corresponds to the voltage in your home.
- Do not bend or crumple the heating pad during use.
- Do not stick needles, pins, safety pins or other sharp objects into the heating pad.
- This appliance must not be used by people insensitive to heat and other very vulnerable people who cannot react to overheating.
- Children under 3 years of age should not use the appliance, as they are not capable of responding to overheating.
- The appliance can be used by children over 3 years of age and under 8 years of age under supervision and with the controller previously set to the lowest temperature level.
- The cable and controller of the appliance may create a risk of entanglement, strangulation, tripping or tripping hazards if not properly arranged. The user must ensure that ties and surplus power cable are securely arranged.
- When using the appliance for a long period of time, it is recommended to set the controller to the lowest temperature level to avoid the risk of overheating.
- Avoid falling asleep while the heating pad is on.
- Prolonged use of the appliance at maximum temperature may cause skin burns.
- Do not cover the appliance with any other object. The

controller should never be covered or placed on or under the heating pad while it is in operation.

- Do not sit on the heating pad. Instead, place it on the part of the body to be treated. It is designed to be placed generally on the back and neck.
- Do not touch the appliance if it has fallen into water. Unplug it immediately from the socket.
- Keep the power cable away from hot surfaces.
- Do not pull on the power cable or the appliance to disconnect the plug from the socket.
- Do not use the heating pad when it is wet or damp. Only use it in dry rooms (avoid using it in the bathroom or similar places).
- Do not expose the controller and power cable to weather or moisture.
- In case of damage, do not try to repair the appliance on your own. Repairs must only be carried out by qualified personnel.
- The heating pad should only be used with the controller provided with the appliance.
- When storing the heating pad, allow it to cool before folding.
- If you have any health problems, consult your doctor before using the heating pad.
- Do not use the heating pad on parts of the body that are inflamed, injured or swollen.
- If you experience constant pain in muscles and joints, inform your doctor. This type of pain can be a symptom of a serious illness.
- If you experience any pain or discomfort while using the appliance, stop using it immediately.

Symbol	Meaning
	Read the instruction manual
	Class II appliance
	No pins or needles
	Do not use folded or bent
	Not suitable for children under 3
	For domestic use only
	Normal wash, max. 30 °C
	Do not use bleach
	Do not tumble dry

	Do not iron
	Do not dry clean

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1. BEFORE USE

- This appliance is packaged in a way as to protect it during transport. Take the appliance out of its box and remove all packaging materials. You can keep the original box and other packaging elements in a safe place to prevent damage to the device if you need to transport it in the future. In case the original packaging is disposed of, make sure all packaging materials are recycled accordingly.
- Make sure all parts and components are included and in good conditions. If there is any piece missing or in bad conditions, contact the official Cecotec Technical Support Service immediately.
- Do not remove the serial number of the appliance in order to keep a correct traceability of it in case of assistance.

Box content

- Heating pad
- Instruction manual

NOTE:

The graphics in this manual are schematic representations and may not exactly match the device.

2. OPERATION

Controller

Figure 1:

1. Level display
2. Increase button
3. Power button
4. Decrease button

Switch on

Power button: press the power button to switch on the heating pad. The appliance will start to heat up and the level display will switch on, showing the current temperature level.

Temperature settings

Increase button (+): press this button to increase the temperature level. Each press of the button will increase the temperature level by one level, up to a maximum of nine levels.

Decrease button (-): press this button to decrease the temperature level. Each press of the button will reduce the temperature level by one level.

Temperature levels: the heating pad has nine temperature levels: level 1 is the lowest and level 9 is the highest. Adjust the temperature to your preference.

Level display

The controller display will show the current temperature level. Be sure to adjust to a comfortable level to avoid overheating the skin.

Auto shutdown

The heating pad is equipped with an auto shutdown function for safety. After 120 minutes of continuous use, the pad will automatically switch off. If you wish to continue using it, you must switch it on again by pressing the power button.

3. CLEANING AND MAINTENANCE

- Before washing the heating pad, unplug it from the mains and allow it to cool for at least ten minutes.
- The appliance has a removable cable. First disconnect the plug and then remove the cable.
- You can wash the pad in the washing machine using a delicate washing programme. If you prefer, you can wash it by hand. To do this, place the pad in a basin with lukewarm water and a little detergent and then carefully wring it out. Rinse the pad several times to remove all detergent.
- Store the appliance in a clean, dry place, in a horizontal position, making sure that it is not covered by any other object.
- Use the heating pad once it is completely dry.

4. TECHNICAL SPECIFICATIONS

Product reference: EU01_100641

Product: HeatConfort 6000 ThermalPad

Power: 100 W

Voltage: 220 – 240 V~

Frequency: 50 / 60 Hz

Technical specifications may change without prior notification to improve product quality.
Made in China | Designed in Spain

5. DISPOSAL OF OLD ELECTRICAL AND ELECTRONIC APPLIANCES

This symbol indicates that, according to the applicable regulations, the product and/or batteries must be disposed of separately from household waste. When this product reaches the end of its shelf life, you should dispose of the cells/batteries/accumulators and take them to a collection point designated by the local authorities.

Consumers must contact their local authorities or retailer for information concerning the correct disposal of old appliances and/or their batteries.

Compliance with the above guidelines will help protecting the environment.

6. TECHNICAL SUPPORT AND WARRANTY

Cecotec shall be liable to the end user or consumer for any lack of conformity that exists at the time of delivery of the product under the terms, conditions, and deadlines established by the applicable regulations.

It is recommended that repairs be carried out by qualified personnel.

If at any moment you detect any problem with your product or have any doubt, do not hesitate to contact the official Cecotec Technical Support Service at +34 96 321 07 28.

7. COPYRIGHT

The intellectual property rights over the texts in this manual belong to CECOTEC INNOVACIONES, S.L. All rights reserved. The contents of this publication may not, in whole or in part, be reproduced, stored in a retrieval system, transmitted, or distributed by any means (electronic, mechanical, photocopying, recording or similar) without the prior authorization of CECOTEC INNOVACIONES, S.L.

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