# Handleiding

## User manual

# PureFeet

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comfore

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## Warnings and Safety Instructions

- Consult a physiotherapist or doctor if you are in any doubt about the product.
- Do not put objects or other items between the rotating massage buttons.
- Do not pull on the cord to remove the plug from the socket.
- Do not immerse or wash the device with water.
- Bruising can result from misuse. Stop immediately the moment it hurts or feels uncomfortable.
- Check and examine the foot massager for any damage before use
- Do not remove any stuck parts yourself.
- Do not use the massager, or consult your doctor first if: A pregnancy has been established, you have a pacemaker, artificial joints or electronic implants, you suffer from one or more of the following ailments or objections: circulatory disorders, varicose veins, open wounds, bruises, torn skin, inflammation of the veins.

## Storaging

- Store the device in a dry, clean place. Do not store the device in extremely high or low temperatures.

## Maintenance

- Clean the surface with a dry cloth. Do not use aggressive abrasives or cleaning products. Please note, the massager must NOT come into contact with water. The fabric part where you rest your feet can be detached. This fabric part can be washed in the washing machine.

## Operating Instructions

- ❗ Do not use the device for more than 20 minutes at a time.
- ❗ The device is only suitable for massaging feet. Do not insert other body parts or objects into the device.
- ❗ Do not stand in the device, the device is not built to support body weight.

## Control Panel

Timer

Massage modes

On/off

Heater button

Massage strength

## **On/off**

Turn the device on or off. After 15 minutes of continuous use, the device will turn itself off automatically (unless you manually set the device for 20 minutes or more).

## **Massage modes**

Press the button to switch between 3 different modes.

P1: The massager under the soles of the feet and the air compression massage are on. The pressure can be adjusted and the heating function can be turned on or off.

P2: The massager under the soles of the feet and the air compression massage are on. The massage buttons on the bottom now move in a different direction.

P3: Only the air compression massage is on. The pressure can be adjusted with the strength knob. Heating function can be turned on and off.

## **Strenght**

Press the button to switch between 5 different massage strenghts. The device is automatically set to L1.

## **Timer**

Press the button to set the timer. You can choose between 10, 15, 20 and 30 minutes. The default timer setting is 15 min.

## **Heater**

The heating function is switched on automatically. The massage buttons light up red to indicate that the heating function is on. Press the button to turn the function off.

## Specifications

|                         |                             |
|-------------------------|-----------------------------|
| Size                    | 44 x 24,4 x 37,8 cm         |
| Output                  | 30W                         |
| Input voltage           | 12V/2.5A                    |
| Material                | ABS                         |
| Automatic shutdown      | Na 10, 15, 20 of 30 minuten |
| Number of massage modes | 5                           |
| Intensity levels        | 5                           |
| Overheating protection  | Ja                          |
| Weight                  | 3,8 kg                      |
| Cord length             | 1,7 m                       |
| Heat Levels             | 1                           |

## Troubleshooting

| Problem                                              | Cause                                                                   | Solution                                               |
|------------------------------------------------------|-------------------------------------------------------------------------|--------------------------------------------------------|
| The device does not turn on.                         | The power adapter is not connected.                                     | Connect the power adapter.                             |
|                                                      | The plug is connected incorrectly.                                      | Insert the plug correctly.                             |
| The device turns off during use.                     | The device turns off automatically after 15 minutes of continuous use.  | Turn the device off and let it cool down for one hour. |
| The device makes a faint sound when it is turned on. | The motor is making noise. The rotating massage nodes are making noise. | This is normal, the product is not defective.          |

## Warranty

We grant 1 year warranty on the massage device. The customer must be able to show proof of purchase and the warranty must be made valid within a period of 1 year from the date of purchase. Please note, the warranty does not apply: for damage caused by incorrect use or repair; for wearing parts; for defects of which the customer was aware at the time of purchase; in case of defects caused by the customer's own fault, in case of damage by third parties.

comfore

© Comfore