



Contraindications

DO NOT use this device to treat any other conditions apart from those listed in Section 1.0 of User Manual.

DO NOT use this device if you are pregnant, planning on becoming pregnant or breastfeeding.

DO NOT use the Maysama PRANA LED Light Therapy Mask if you suffer from light induced headaches.

DO NOT use this device if you suffer from any genetic conditions of the eyes. This device has not been tested under/for these conditions and therefore the risk is unknown.

DO NOT use this device if you suffer from lupus erythematosus, photosensitive eczema, or albinism – severe skin reaction can occur.

DO NOT use this device if you suffer from any photosensitive disorder (sensitization to light) – severe skin reaction may occur.

DO NOT use this device if you are taking any medication that can cause photosensitivity – can cause severe skin reaction.

DO NOT use this device, without consulting your healthcare provider, if you suffer from physical, sensory or mental illness or are undergoing treatment for a dermatology condition.

DO NOT use Maysama PRANA LED Light Therapy Mask on the same day as skincare product containing photosensitive ingredients such as Alpha Hydroxy Acids.

DO NOT use Maysama PRANA LED Light Therapy Mask without first cleaning your face and removing any emollients, face oil or SPF.

DO NOT use Maysama PRANA LED Light Therapy Mask with skincare products other than those recommended by Maysama for use with Maysama PRANA LED Light Therapy Mask. Other products may reduce the efficacy of your mask or cause reaction.

Photosensitivity is a common side effect of various medications. These can include certain antibiotics, chemotherapy drugs and diuretics. If you are unsure about any medication you may be taking, consult your healthcare provider.

Other substances not listed above can also cause photosensitivity. Common examples of these substances are: St John's Wort, coal tar, deodorant, antibacterial soaps, artificial sweeteners, naphthalene (mothballs), petroleum products, brightening agents found in laundry detergent and cadmium sulphide (a chemical injected into the skin during tattooing)