



# Mito Red Light MitoADAPT Red Light Therapy Floor Stand User Guide


[Home](#) » [Mito Red Light](#) » Mito Red Light MitoADAPT Red Light Therapy Floor Stand User Guide 

## Contents

- [1 Mito Red Light MitoADAPT Red Light Therapy Floor Stand](#)
- [2 Product Information](#)
- [3 Key Features](#)
- [4 Product Usage Instructions](#)
- [5 Quick set-up without using the Mite Red Light opp](#)
- [6 Quick set-up using the Mito Red Light opp](#)
- [7 Usage Guidelines](#)
- [8 Troubleshooting](#)
- [9 Documents / Resources](#)
  - [9.1 References](#)



**Mito Red Light MitoADAPT Red Light Therapy Floor Stand**



## Product Information

The MitoADAPT is a user-friendly device designed to provide red light therapy for various health and wellness benefits. It utilizes advanced LED technology to emit specific wavelengths of red light that can penetrate the skin and stimulate cellular activity.

## Key Features

- Compact and portable design
- Adjustable intensity settings
- Multiple treatment modes
- Easy-to-use interface
- Safe and non-invasive

## Product Usage Instructions

Before using the MitoADAPT, please read the user manual and watch the instructional videos available at <https://mitoredlight.com/pages/user-manual>.

1. Connect the power adapter to the device and plug it into a suitable power outlet.
2. Press the power button to turn on the MitoADAPT.
3. Select the desired treatment mode using the mode button. Refer to the user manual for information on different modes and their specific benefits.
4. Adjust the intensity level using the intensity control buttons. Start with a lower intensity and gradually increase it as desired.
5. Position the device at a suitable distance from your body or target area. Follow the recommendations provided in the user manual for optimal treatment distance.
6. Place the device directly facing the skin or area to be treated. Ensure that the red light is evenly distributed over the target area.
7. Set the treatment duration using the timer function if available. Follow the recommended treatment times provided in the user manual.
8. During the treatment, relax and allow the red light to penetrate the skin. You may choose to close your eyes or wear protective goggles if recommended.
9. After the treatment session, press the power button to turn off the device.
10. Store the MitoADAPT in a cool, dry place when not in use.

For any further assistance or inquiries, please visit our website at [www.mitoredlight.com](http://www.mitoredlight.com) or follow us on social media:

- Facebook: [@mitoredlightofficial](#)
- Instagram: [@MitoRedLight](#)
- Twitter: [@mitoredlightofficial](#)

Congratulations on your new MitoADAPT panel! Below is a quick startup guide to get you on your way to your red light therapy pronto! If you need more detailed instructions or need more information, check out our main user manual at [www.mitoredlight.com/pages/user-manual](https://www.mitoredlight.com/pages/user-manual)

### Quick set-up without using the Mito Red Light app



Unbox your panel, plug in the power cord, and turn on the power switch.

This is the Home Screen. Using the '+' and '-' buttons you can adjust the time for the current session. If all the other settings look good, press 'GO' to start and stop your session.



From the Home Screen press 'SET' to navigate to the Settings Screen. To adjust the Default Time, Mode, Bluetooth connectivity, and Brightness settings press 'SET'.

Press + /- to change any specific setting. Press 'SET' to toggle between each setting.



Once you have adjusted these settings, press 'GO' to go back to the Home Screen, and then press 'GO' again to start/stop your session.

Please note, changing the Default Time on the Settings Screen will not affect the current session. The panel must be turned off and back on to display the new Default Time.



## Quick set-up using the Mito Red Light app



Yes, we have an app so you can control your MitoADAPT from your mobile device.



Unbox your panel, plug in the power cord, and turn on the power switch. If you have multiple devices please make sure the panels are connected before starting the session.

Make sure the Bluetooth function is enabled on your primary device.



Download the Mito Red Light app using the QR codes below! Start the app and fill out the brief onboarding questions. Accept permissions and make sure the Bluetooth on your mobile device is turned on.



Select one of the 11 modes, and follow the prompts. If you need to stop your session, press 'Stop' on your app or press 'Go' on the main panel.

Note: Modes 9, 10, and 11 are near infrared light only. NIR light is invisible to the naked eye so it may appear that the session is not running. Please rest assured that significant light energy is being emitted from the panel.



## Usage Guidelines

For first time users, we recommend starting out with 2-3 minute treatments per area, and slowly working up to a 10-minute session over the course of 2-3 weeks as your body acclimates to the therapy.

## Modes

- We recommend starting with Mode 1 as you get acclimated to the light.
- Then we recommend rotating your sessions through Modes 1-6 (i.e. one day do a session with Mode 1, the next session use Mode 2, etc).
- Modes 1 through Mode 6 are various combinations of Red (peak 630nm & 660nm) and NIR (peak 810nm & 850nm) light. See the user manual for a summary of the modes.
- Modes 7-10 have only 50% of the LEDs lit and can be used for gentler sessions. Mode 11 is 100% NIR and is available only through the app.

## Troubleshooting

- If the control panel is not working, make sure you are connected to power and that the power switch on the back of the device is turned on.
- When resetting the Default Time on the Settings Screen, after the time is selected hit the 'GO' button then turn the power switch off then back on.
- If there is no light coming from the panel when the session is started, confirm that the brightness is not set to 0% (also note that NIR light 830nm/ 850nm is invisible to the naked eye).
- If app will not connect, make sure you are within 5 feet of the panel, close the app on your device and restart. Turn the panel off and on.
- Make sure your Bluetooth is enabled on your mobile device and your panel.

<http://www.mitoredlight.com>

User Manual and Instructional Videos This user manual and instructional videos can be found here:

<https://mitoredlight.com/pages/user-manual> Updated June 18th 2023

## Documents / Resources



[Mito Red Light MitoADAPT Red Light Therapy Floor Stand](#) [pdf] User Guide  
MitoADAPT Red Light Therapy Floor Stand, MitoADAPT, Red Light Therapy Floor Stand, Therapy Floor Stand, Floor Stand, Stand

## References

- [Mito Red Light | Red Light Therapy Devices for your Home](#)
- [Mito Red Light | Red Light Therapy Devices for your Home](#)
- [User Manuals and Video Tutorials | Mito Red Light](#)