



Qumax XMTZC05HM Bluetooth Smart Scale User Manual

[Home](#) » [Qumax](#) » Qumax XMTZC05HM Bluetooth Smart Scale User Manual 

Qumax XMTZC05HM Bluetooth Smart Scale User Manual





Contents

- [1 Check required specification for Phone or Tablet](#)
- [2 Install the batteries](#)
- [3 Install App](#)
- [4 Pair the devices, register and complete user information](#)
- [5 Weighing](#)
- [6 View the test result](#)
- [7 Specification](#)
- [8 Instruction for OKOK App](#)
- [9 Documents / Resources](#)
- [10 Related Posts](#)


Check required specification for Phone or Tablet

Before installing OKOK App, kindly check if the phone or tablet meets following requirement:

-  Compatible with iOS 7.0 and above with Bluetooth 4.0 and above
 - Phone 4s, iPhone 5, iPhone 5c, iPhone 5s
 - iPhone 6, iPhone 6 Plus, iPod touch 5, iPad 3
 - iPad 4, iPad Air, iPad Air 2, iPad mini
 - iPad mini 2, iPad mini 3
-  Compatible with Android 4.3 and above with Bluetooth 4.0 and above
 - **SAMSUNG:** galaxy s5, galaxy s4, galaxy s3, galaxy note4, galaxy note3, galaxy note2, galaxy grand2
 - **XIAOMI:** XIAOMI Note, XIAOMI4, XIAOMI3, XIAOMI2/2s, XIAOMI tablet, HONGMI Note 4G, HONGMI 1s
 - **HUAWEI:** Ascend P8, Ascend P7, Honor 6Plus, Honor 6, mate 7, Ascend G7

Install the batteries

1. Open the cover of battery
2. Install the batteries in the right direction
3. Close the cover



Install App

1. Switch on Bluetooth
2. Phone or Tablet connects to Net (WiFi, 3G or 4G)



Pair the devices, register and complete user information

- Open APP, begin to weigh and pair



- Register as a user



- Complete personal information



Weighing

- Enter into “Dynamic Weight” interface



- Stand on the scale and synchronize weight data onto APP



View the test result

- Enter into “Weight Trend” interface



- View test report



Specification

Display:	All data displayed when power on
Capacity:	3kg-180kgs
Division:	1D=100g/0.2lb
Unit:	KG/LB/ST
Voltage:	4.5V
Battery:	3 *AAA battery
Low battery indication:	2.4V±0.2V, display “Lo” and power off
Current:	<25mA
Standby current:	<10 μ A
Overload indication:	Overload more than 2.5kgs, total 182.5kgs display “Err”
Best working temperature:	5~35 °C
Storage temperature:	-10~60 °C
Moisture Environment:	<90%RH

Instruction for OKOK App

- **Instantly weighing**

Open the OKOK APP, it will remind you to switch on Bluetooth. Turn Bluetooth on and the auto-pair process completed automatically. User begin weighing, the weight data will be synchronizing to APP. If the Bluetooth disconnected, the weight data will be saved into memory of the scale, then synced to APP while the Bluetooth connected.

- **Type weight data into APP**

The user can type weight data into APP to keep a track

- **Current weight**

Current weight interface displays weight information of the user and family members, and generates the report

and give health advice based on the measurement. Including current weight interface, weight list interface, date report and notice.

- **Weight trend**

Including weekly report, monthly report and quarterly report, indicating the changes of weight through curves.

Notice: For fat scale, there are 8 data: weight, body fat, visceral fat, muscle, bone, BMR, age, BMI.

- **Management of scale**

Management of scale including two parts: scale management and connection management.

Scale management

Manage the process of pairing, bind and unbind.

Connection management

Including Bluetooth status indication and Bluetooth connection.

- **Account and user**

Including account management and user management.

Account management

Including registration, login in, change password, forget password, links to Facebook, Twitter...etc

User management

Manage user information, including nickname, birthday, height, gender...etc, also managing family member account and information (under the same registration account).

Warnings

1. Please keep the scale away from wet floor, don't use the scale if the glass is wet in case of slip.
2. Please stand on the middle of the scale in case of turnover.
3. The platform of the scale is made of glass, please avoid sharp or harsh things scratching the surface, excessive forces or sudden shocks, unexpected fallen-off will break the glass and injure the user.
4. The Bluetooth scale is not waterproof. Keep the surrounding environment dry, do not submerge in water.
5. Use soft cloth with neutral cleaning liquids to clean the surface, keep away contacting with organic solvent nor chemicals.
6. When the batteries are in low power status, please remove the batteries in case of leaking liquids corrode the scale.

Notices

1. Please stand on the scale in the right position.
2. Please stand on the scale no less than 10 seconds if you measure body fat.
3. Please confirm gender, height, weight data is correct.





[Qumax XMTZC05HM Bluetooth Smart Scale](#) [pdf] User Manual
XMTZC05HM, Bluetooth Smart Scale