

# TALEC WS2601 Bluetooth Smart Scale User Manual

## Contents

- 1 Bluetooth Smart Scale chinese and english
  - 1.1 Step 1 Check required specification for Phone or Tablet
  - 1.2 Step 2 Install the batteries
  - 1.3 Step 3 Install App
  - 1.4 Step 4 Pair the devices, register and complete user information
  - 1.5 Step 5 Weighing
  - 1.6 Step 6 View the test result
  - 1.7 Specification
  - 1.8 Instruction for OKOK App
    - 1.8.1 Scale management
    - 1.8.2 Connection management
    - 1.8.3 Account Management
    - 1.8.4 User management
  - 1.9 Warnings
  - 1.10 Notices
- 2 Documents / Resources
  - 2.1 References
- 3 Related Posts

## Bluetooth Smart Scale chinese and english



Sync your progress wirelessly to any Bluetooth smart ready device.



**Smartphones not required for scale use.**




### Step 1 Check required specification for Phone or Tablet



Before installing OKOK App, kindly check if the phone or tablet meets following requirement:

 <p><b>Compatible with iOS7.0 and above with Bluetooth 4.0 and above</b></p>	<p>iPhone 4s, iPhone 5, iPhone 5c, iPhone 5s          iPhone 6, iPhone 6 Plus, iPhone 6s, iPhone 6s Plus          iPhone 7, iPhone 7 Plus, iPhone 8, iPhone 8 Plus, iPhone X          iPad 1, iPad 2, iPad 3, iPad 4, iPad 5, iPad mini          iPad mini 2, iPad mini 3, iPad mini 4, iPad Air          iPad Air 2, iPad Pro, iPad Pro 2</p>
 <p><b>Compatible with Android 4.3 and above with Bluetooth 4.0 and above</b></p>	<p>SAMSUNG: galaxy s8, galaxy s7, galaxy s6, galaxy s5, galaxy s4          galaxy s3, galaxy note8, galaxy note7, galaxy note6, note5          galaxy note4, galaxy note3, galaxy note2, galaxy grand2XIAOMI: Hongmi Note 4x, Xiaomi 5c, Hongmi 4x, Xiaomi Note 2          Xiaomi 6, Xiaomi Max 2, Hongmi 5x, Hongmi Note 5A, Xiaomi Mix 2          Xiaomi Note 3, Hongmi 4, Hongmi 5, Hongmi 5 Plus, Xiaomi 4</p> <p>HUAWEI: Ascend P8, Ascend P7, Honor6Plus, Honor6, Mate 10          mate 10Pro, Mate 9, Mate 9Pro, AscendG7, nove 2, nove 2Plus          P10, P9, P10 Plus, P9 Plus</p> <p><b>Support Android4.3 and above for selected Android devices with Bluetooth 4.0 and above</b></p>

## Step 2 Install the batteries



1/ – Open the cover of battery

2/ – Install the batteries in the right direction

3/ – Close the cover

## Step 3 Install App

Use your device to scan QR Code or search OKOK from App Store/Google Play to install App.

1/ Switch on Bluetooth

2/Phone or Tablet connects to Net (WiFi, 3G or 4G)




#### Step 4 Pair the devices, register and complete user information


Open APP, begin to weigh and pair



Register as a user



Complete personal information



#### Step 5 Weighing

Enter into "Dynamic Weight" interface



Stand on the scale and synchronize weight data onto APP



## Step 6 View the test result

Enter into “Weight Trend” interface



View test report



## Specification

Display: All data displayed when power on

Capacity: 3kg-180kgs

Division: 1D=100g/0.2lb

Unit: KG/LB/ST

Voltage: 3V

Battery: 2\*AAA battery

Low battery indication:  $2.4V \pm 0.2V$ , display “Lo” and power off

Current:  $\leq 25mA$

Standby current:  $\leq 10 \mu A$

Overload indication: overload more than 2.5kgs, total 182.5kgs display “Err”

Best working temperature:  $5 \sim 35^{\circ}C$

Storage temperature :  $-10 \sim 60^{\circ}C$

Moisture Environment:  $\leq 90\%RH$

## Instruction for OKOK App

- Instantly weighing

Open the OKOK APP, it will remind you to switch on Bluetooth. Turn Bluetooth on and the auto-pair process completed automatically. User begin weighing, the weight data will be synchronizing to APP. If the Bluetooth disconnected, the weight data will be saved into memory of the scale, then synced to APP while the Bluetooth connected.

- **Type weight data into APP**

The user can type weight data into APP to keep a track

- **Current weight**

Current weight interface displays weight information of the user and family members, and generates the report and give health advice based on the measurement. Including current weight interface, weight list interface, date report and notice.

- **Weight trend**

Including weekly report, monthly report and quarterly report, indicating the changes of weight through curves.



For fat scale, in addition to weighing data, it also includes body shape, obesity, weight, BMI, body fat rate, Fat weight, muscle rate. muscle weight. visceral fat, water, water content, basal metabolism. bone mass, Protein, fat free weight fifteen item data

- **Management of scale**

Management of scale including two parts: scale management and connection management

#### **Scale management**

Manage the process of pairing, bind and unbind

#### **Connection management**

Including Bluetooth status indication and Bluetooth connection

- **Account and user**

Including account management and user management

#### **Account Management**

Including registration, login in, change password, forget password, links to Facebook, Twitter...etc

#### **User management**

Manage user information. including nickname, birthday, height, gender. etc, also managing family member account and information (under the same registration account).

#### **Warnings**

1. Please keep the scale away from wet floor, don't use the scale if the glass is wet in case of slip.
2. Please stand on the middle of the scale in case of turnover.
3. The platform of the scale is made of glass, please avoid sharp or harsh things scratching the surface, excessive forces or sudden shocks, unexpected fallen-off will break the glass and injure the user.
4. The Bluetooth scale is not waterproof. Keep the surrounding environment dry. do not submerge it in water.
5. Use soft cloth with neutral cleaning liquids to clean the surface. keep away contacting with organic solvent nor chemicals.
6. When the batteries are in low power status, please remove the batteries in case of leaking liquids corrode the scale.

#### Notices

1. Please stand on the scale in the right direction.
2. Please stand on the scale no less than 10 seconds if you measure body fat.
3. Please confirm gender, height, weight data is correct.

#### Warning:

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

#### NOTE:


This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

#### FCC Statement:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

#### Documents / Resources

	<p><a href="#">TALEC WS2601 Bluetooth Smart Scale</a> [pdf] User Manual WS2601 Bluetooth Smart Scale, WS2601, Bluetooth Smart Scale, Smart Scale, Scale</p>
---	---

## References

- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.