

**TimerCap**

TC-38  
DISPLAYS  
TIME  
REMINDER  
CAP



## TimerCap TC-38 Displays Time Reminder Cap User Manual

[Home](#) » [TimerCap](#) » TimerCap TC-38 Displays Time Reminder Cap User Manual 

### Contents

- [1 TimerCap TC-38 Displays Time Reminder Cap](#)
- [2 INTRODUCTION](#)
- [3 SPECIFICATIONS](#)
- [4 WHAT'S IN THE BOX](#)
- [5 PRODUCT OVERVIEW](#)
- [6 FEATURES](#)
- [7 SETUP GUIDE](#)
- [8 CARE & MAINTENANCE](#)
- [9 TROUBLESHOOTING](#)
- [10 PROS & CONS](#)
- [11 WARRANTY](#)
- [12 FREQUENTLY ASKED QUESTIONS](#)
- [13 References](#)
- [14 Related Posts](#)

**TimerCap**<sup>®</sup>

**TimerCap TC-38 Displays Time Reminder Cap**



## INTRODUCTION

It's easy to keep track of your medication plan with the TimerCap TC-38 Displays Time Reminder Cap. This clever cap shows you a warning every time the bottle is opened, which helps you remember when to take your medicine. This cap was made by TimerCap LLC and is small, cheap, and simple to use. It's perfect for people who have trouble remembering to take their pills. The TimerCap TC-38, which costs only \$13.95, is an effective and affordable way to make sure that medications are taken on time. It was first released on December 11, 2013, and comes with an LR41 battery that can power it for 12 to 18 months. This simple but useful cap shows the time every time the bottle is opened, which makes it a great tool for keeping track of daily medications.

## SPECIFICATIONS

<b>Brand</b>	TimerCap
<b>Price</b>	\$13.95
<b>Product Dimensions</b>	2 x 2 x 4 inches
<b>Weight</b>	1.25 ounces
<b>Item Model Number</b>	TC-38
<b>Battery Life</b>	12-18 Months
<b>Batteries</b>	1 LR41 battery required (included)
<b>Date First Available</b>	December 11, 2013
<b>Manufacturer</b>	TimerCap LLC

## WHAT'S IN THE BOX

- Displays Time Reminder Cap
- battery
- User Manual

## PRODUCT OVERVIEW



## FEATURES

- **No Programming Needed:** The TimerCap works on its own and doesn't need to be set up or programmed. The stopwatch will begin to count down the seconds, minutes, and hours as soon as the cap is closed.
- **Real-Time Time Display :** The cap shows the time instantly since the last time the pill bottle was opened. This helps users keep track of when they last took their medicine.
- **Easy Medication Management:** This is a great tool for people with long-term conditions, people who are recovering from surgery, or anyone else who needs an easy way to keep track of their medications.
- **Color-Coded for Easy Identification :** TimerCaps come in different colors that make it easier for users to tell the difference between medications, times of use, or family members (like mom, dad, kid) who are taking the medicine.
- **Portability:** The TimerCap's small size makes it easy to carry in a pocket, purse, or travel bag, so you can easily handle your medications even when you're out and about.
- **Improves Medication Compliance:** The clear show of time spent can help people stick to their medication schedules by reducing confusion.
- **Built-in Stopwatch:** The stopwatch starts automatically when the bottle is closed, so you can get an exact reading of how long it's been since the bottle was last opened.
- **Helps Stop Abuse:** The timer keeps track of time since the last dose, which keeps other people from trying to play with the medicine.
- **Long Battery Life:** The TimerCap runs on one LR41 battery, which is included. The battery lasts for 12 to 18 months, so you don't have to change it as often.
- **CRC (Child-Resistant Cap):** TimerCap has a Child-Resistant Cap that keeps kids from quickly getting to the medicine inside.
- **No Loud Alarms :** The TimerCap is different from other reminder devices because it doesn't have any loud alarms. Instead, it uses time tracking to provide a quiet and discreet reminder system.
- **Available in Multiple Quantities:** TimerCap comes in packs of two to four amber-colored bottles, each with its own time-tracking cap. This makes it easy to take more than one medicine.
- **Design That Is Small and Durable:** The TimerCap is made to last. It has a small and durable design that fits most pill bottles with a 38mm hole.
- **Easy to Use and Understand:** The TimerCap is made to be easy for people of all ages to use. It has a simple, straightforward interface that doesn't require any technical understanding.



## SETUP GUIDE

- **Insert Battery:** To power the timer, start by putting the LR41 battery that came with it into the cap. The battery area is simple to get to.
- **Connect to Pill Bottle:** Connect the TimerCap safely to the pill bottle, making sure the cap fits snugly and is placed correctly.
- **Close the Cap:** To start the time, just close the TimerCap. Since the bottle was last opened, the timer will start keeping track of time on its own.
- **Pick the Right Size:** Make sure that your pill bottle can fit the TimerCap. This cap fits 38mm holes, which are typical on many standard pill bottles.
- **No Need to Configure:** The TimerCap doesn't need to be programmed or configured in any other way; it starts working as soon as it's closed.
- **Check Timer Display:** Once the cap is closed, the timer will start showing how much time has passed. It's been exactly this long since the bottle was last opened.
- **Track Time:** The time shown on the screen can help you figure out when to take the next amount. The picture is a great way to remember not to take your medicine too early or too late.
- **If you need to change your medicine,** all you have to do is take off the old TimerCap and put on a new one

for the new bottle.

- **Use Color Coding to Help You Find Things:** Use TimerCaps of different colors to help you tell the difference between medicines, times of use, or people in your home.
- **Attach to Travel Bottles:** The TimerCap fits easily on travel-sized pill bottles, which makes it easier to keep track of medicines while you're away.
- **Watch the Battery Life:** The battery usually lasts between 12 and 18 months. Check it often and change it when it's time to make sure it keeps working.
- **Manual Reset:** To restart the clock, just open the cap and close it again. This will reset the timer.
- **Check for Compatibility:** Make sure your bottle has the right 38mm hole for the TimerCap before you use it.
- **Check the organization of your medications:** The TimerCap can help you sort various medications by time of day or person taking them (like mom, dad, child).
- **How to Store It:** Keep your TimerCap in a cool, dry place when not in use to make sure it stays in good shape for a long time.

## CARE & MAINTENANCE

- **Clean it often:** Use a damp cloth to wipe the TimerCap clean of dust, dirt, and any medicine leftovers.
- **Do Not Submerge in Water:** Do not put the TimerCap in water. To keep the mechanism from breaking, clean it with a wet cloth instead.
- **Replace Battery When It's Dead:** Keep an eye on the battery life and replace the LR41 battery when the screen starts to fade or stops working.
- **How to Store It:** Keep the TimerCap in a cool, dry place when not in use to keep it working.
- **Handle Carefully:** The TimerCap is a small device that should be handled carefully so that you don't break the cap or hurt the internal timer.
- **Check for Damage:** Make sure you regularly check the TimerCap for any damage that you can see, like cracks or problems with the screen.
- **Use Only on Bottles That Fit:** Make sure the TimerCap is only used on bottles with a 38mm hole to keep the seal from going wrong and to get accurate tracking.
- **Do Not Use Too Much Force:** If the TimerCap does not fit properly, do not force it onto the bottle. It could be bad to force the cap on.
- **Reset Timer If Needed:** To restart the clock, just open the cap and close it again. This will reset the timer.
- **Replace the TimerCap If It Gets Worn Out:** If the cap gets worn out or hard to use, get a new one to keep it working.
- **When setting up the stopwatch,** always make sure that the TimerCap is tightly attached to the bottle so that it starts counting correctly.
- **Use the Mute Option (If You Want to):** If you'd rather not see the timer's visual warning, use the mute option to turn off any other sounds.
- **Take Only One drug at a Time:** The TimerCap is made to hold only one drug at a time. Do not mix medicines that come in the same bottle.
- **Check Expiration Date:** Look at the medicine's expiration date on the bottle and make sure you use the TimerCap within the medicine's useful time.
- **Keep Out of Reach of Children:** The cap is child-resistant, but small children should always make sure it is securely closed and kept somewhere out of their reach to avoid accidentally eating it.

## TROUBLESHOOTING

Issue	Solution
Cap not working	Ensure the LR41 battery is correctly installed and has power.
Reminder not displaying	Check the battery life and replace the battery if necessary.
Cap doesn't fit the bottle	Make sure the cap is compatible with the bottle size you're using.
Time reminder appears late	Verify that the cap is properly aligned and tightly fitted on the bottle.
Time display is unclear	Clean the display area to remove any dirt or smudges.
Cap is loose or falls off	Ensure the cap is securely twisted onto the bottle.
Battery drains quickly	Check for corroded terminals or ensure you're using a new LR41 battery.
No time reminder after opening	Confirm the sensor inside is not obstructed and is functioning properly.
Cap is difficult to open	Make sure the cap is not over-tightened or stuck due to any residue.
Timer reminder keeps resetting	Verify that the battery is correctly installed and not loose.
Cap making strange noises	Inspect for any damage to the cap or internal components.
Reminder is out of sync	Try resetting the cap and reinstalling the battery to fix the time sync.
Cap doesn't stay on the bottle	Check if the bottle opening is the right size for the cap.
Cap is damaged	If there are any cracks or breaks, contact the manufacturer for a replacement.
Reminder display is too small	Consider using a larger font or clearer screen display for better visibility.

## PROS & CONS

### PROS:

- Displays time reminders to help you stay on track with your medication schedule.
- Affordable price of \$13.95, making it an accessible solution for many.
- Long battery life of 12-18 months, reducing the need for frequent replacements.
- Compact and lightweight design (1.25 ounces) fits most pill bottles comfortably.
- Easy to use, requiring no complicated setup to get started.

### CONS:

- Only works with certain bottle sizes, so it may not fit all pill containers.
- Battery replacement might be needed after 12-18 months, which could be inconvenient.
- The time reminder display may not be large enough for some users to read easily.
- Limited features, as it primarily functions as a simple reminder cap.
- May not be suitable for users who require more complex medication management tools.

## WARRANTY

The **TimerCap TC-38 Displays Time Reminder Cap** comes with a **1-year limited warranty** from the

manufacturer, **TimerCap LLC**. This warranty covers defects in material and workmanship under normal use. If the product malfunctions or breaks due to a manufacturing defect, TimerCap LLC will provide a replacement or refund within the warranty period. Be sure to keep your original purchase receipt as proof of purchase to ensure you can take advantage of the warranty.

## FREQUENTLY ASKED QUESTIONS

What is the price of the TimerCap TC-38 Displays Time Reminder Cap?

The TimerCap TC-38 Displays Time Reminder Cap is priced at \$13.95, providing an affordable solution for medication time reminders.

What are the dimensions of the TimerCap TC-38 Displays Time Reminder Cap?

The TimerCap TC-38 Displays Time Reminder Cap has product dimensions of 2 x 2 x 4 inches, making it compact and easy to use with most medication bottles.

How much does the TimerCap TC-38 Displays Time Reminder Cap weigh?

The TimerCap TC-38 Displays Time Reminder Cap weighs 1.25 ounces, making it lightweight and easy to handle.

What is the battery life of the TimerCap TC-38 Displays Time Reminder Cap?

The TimerCap TC-38 Displays Time Reminder Cap has a battery life of 12-18 months, ensuring long-term use before needing a battery replacement.

What type of battery does the TimerCap TC-38 Displays Time Reminder Cap use?

The TimerCap TC-38 Displays Time Reminder Cap requires 1 LR41 battery, which is included with the product.

When was the TimerCap TC-38 Displays Time Reminder Cap first available?



The TimerCap TC-38 Displays Time Reminder Cap was first available on December 11, 2013.

Who manufactures the TimerCap TC-38 Displays Time Reminder Cap?

The TimerCap TC-38 Displays Time Reminder Cap is manufactured by TimerCap LLC.

What is the purpose of the TimerCap TC-38 Displays Time Reminder Cap?

The TimerCap TC-38 Displays Time Reminder Cap is designed to help users remember when to take their medication by displaying a time reminder on the cap.

How does the TimerCap TC-38 Displays Time Reminder Cap work?

The TimerCap TC-38 Displays Time Reminder Cap has a built-in timer that displays the time, helping users track when they last took their medication.

How long does the battery last in the TimerCap TC-38 Displays Time Reminder Cap?

The TimerCap TC-38 Displays Time Reminder Cap's battery lasts between 12 to 18 months, depending on usage.

What is the model number of the TimerCap Displays Time Reminder Cap?

The model number of the TimerCap TC-38 Displays Time Reminder Cap is TC-38.

Where is the TimerCap TC-38 Displays Time Reminder Cap manufactured?

The TimerCap TC-38 Displays Time Reminder Cap is manufactured by TimerCap LLC, though the specific country of origin is not provided in the data.

Why is my TimerCap TC-38 Displays Time Reminder Cap not displaying the time?

If the TimerCap TC-38 Displays Time Reminder Cap is not displaying the time, ensure that the cap is securely attached to the bottle and that the internal battery is working. If the cap's display is blank, try replacing the battery or check if it has been installed incorrectly. Also, ensure that the cap is not damaged or malfunctioning.

Why is the time on my TimerCap TC-38 Displays Time Reminder Cap incorrect?

If the TimerCap TC-38 Displays Time Reminder Cap is showing the wrong time, check if the time has been set correctly. To adjust the time, press the Set button.

## References

- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.