



## WYZE WSCALX Scale X Digital Scale User Guide

[Home](#) » [Wyze](#) » WYZE WSCALX Scale X Digital Scale User Guide 



### WSCALX Scale X Digital Scale User Guide

Read this manual carefully and save it for future reference.

**Disclaimer:** The products, accessories, user interfaces, etc. illustrated in this guide are intended for reference only. Due to product updates, the actual product and the illustrations in this guide may vary.

#### WELCOME TO WYZE

BIENVENUE CHEZ WYZE BIENVENIDO A WYZE


Hi! Thanks for getting your hands on a Wyze Scale X, and if this is your first Wyze product, welcome to the Wyze family! You just made a commitment to a better, more connected lifestyle. This guide will help you get things set up.

Contents




- 1 IN THE BOX
- 2 TIPS
- 3 SETUP
- 4 UNITS AND FACTORY RESET
- 5 NOTES
- 6 LED DISPLAY
- 7 EXPLANATION OF SYMBOLS ON UNIT
- 8 GENERAL SAFETY & PRECAUTIONS
- 9 CARE & MAINTENANCE
- 10 WARRANTY INFORMATION
- 11 SPECIFICATIONS
- 12 Documents / Resources
  - 12.1 References
- 13 Related Posts

IN THE BOX

- Wyze Scale X × 1
- 1.5V AAA Batteries × 4
- Quick Start Guide × 1




TIPS

		
Place on a hard, flat surface	Don't wear socks	Don't let your l

## SETUP

1. Before first use, remove the plastic strip from the battery compartment.




**Note:** When inserting the batteries, line up the positive and negative signs to make sure the batteries are installed correctly.

2. Download the Wyze app and sign in or create an account.




3. In the Wyze app, tap the + plus sign, then Add Device > Lifestyle > Wyze Scale X. Follow the instructions in the Wyze app to complete setup.

To connect, enable your phone's Bluetooth and tap on the scale's surface to activate it.




4. After the connection is established, place the scale on a hard, flat surface and step onto the scale with bare feet to begin the measurement.

Make sure your feet are touching the 4 electrodes for the most accurate measurements.



## UNITS AND FACTORY RESET

On the back of your scale, find the UNIT button under the battery cover. Press it to change the units from lb to kg. To factory reset the scale, press and hold the same button for 3s until the LED display shows “Clr”.



## NOTES

- You can share Wyze Scale X with other family members by following the sharing instructions inside the Wyze app.
- Your scale may need to be calibrated to show accurate weight measurements. Place the scale on a hard, flat surface and re-calibrate it (step on and then step off, allowing the scale to return to “0.0 lb”).
- If you work out more than 10 hours a week and have a resting heart rate below 60 bpm, you may need to use Athlete Mode.
- Make sure your age, sex, and height entered during the setup process are correct to get the most accurate results.
- Open the Wyze app and connect to your Wyze Scale X before taking measurements to ensure the data is recorded on the correct profile.
- If your Wyze Scale X is shared by users with similar weights, make sure the scale is connected to the Wyze app of the user taking the measurement or data may not be recorded on the correct profile.
- People that are under 18 or over 65 years of age, that are pregnant, that have osteoporosis or very low bone density, that have edema, or that have artificial dialysis should use Wyze Scale X only as a reference, as the body fat measurement could be inaccurate due to deviations from standard averages.

## LED DISPLAY

Suspected error. Please clean the scale with a damp cloth if needed, then step on the scale again with bare feet.



The batteries are low.  
Please replace the batteries.




An error has occurred. Please step off of the scale, then step on again to repeat the measurement.



Critical low battery.

Please replace the batteries.

LED display will remain off after “Lo” is displayed 3 times.



Weight not in supported range.

Please measure again.




## EXPLANATION OF SYMBOLS ON UNIT



### Warning.

Please recycle where facilities exist. Waste electrical products should not be disposed of with household waste. Check with your local authority or retailer for recycling advice.



## CAUTION

People fitted with an implanted medical device such as a pacemaker, an electronic life support system such as artificial heart/lungs, or that use a portable electronic medical device such as an electrocardiograph, and women who are pregnant, must only use this product while wearing shoes, or in the “Only Measure Weight” setting. Please refer to the manufacturer’s instructions for your medical device for more information. Do not use this

product unless you are wearing shoes, or activate the “Only Measure Weight” setting, if you have an implanted medical device, have an electronic life support system, use a portable electrical medical device, or are pregnant. Pregnant women should also consult their healthcare provider before use.

To activate the “Only Measure Weight” setting:

1. Open the Wyze App and tap Wyze Scale.
2. Tap on the gear to open Scale Settings.
3. Tap “Users,” then select your email address to open the “Account Info” screen.
4. Tap the slider next to “Only Measure Weight,” then “Yes” to confirm.”

Some physical conditions could affect hydration levels that may lead to inaccurate results. Please consult your healthcare provider for more information.

- Only use 1.5V AAA batteries in the scale. Using the incorrect type of batteries creates a risk of
- Do not expose batteries to extremely high temperatures or fire. Doing so creates a risk of explosion or leakage of flammable liquid or gas.
- Do not manually crush, cut or deform batteries. Doing so creates a risk of explosion or leakage of flammable liquid or gas.
- Do not subject batteries to extremely low air Doing so creates a risk of explosion or leakage of flammable liquid or gas.

## **GENERAL SAFETY & PRECAUTIONS**

1. Do not stand on the edge of scale, as you may fall over or influence the accuracy of the measurement.
2. Do not use the scale on a tiled or wet floor, as this may result in a fall.
3. When measuring your heart rate with your phone’s camera lens, do not touch your phone’s camera flash. Your phone’s camera flash could overheat and burn your fingers.
4. Make sure the surface of the scale is clean before you step onto it, to prevent any slippage.
5. Treat your scale with care. Do not drop it or jump on it.

The scale is designed to be stood on; misuse or abuse may break the electronic sensors, cause injury, or result in inaccurate measurements.

6. Avoid using this scale near strong magnetic fields, such as microwave ovens, etc.
7. This device complies with part 15 of the FCC Rules.

Its operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

8. Changes or modifications not expressly approved by the party responsible for compliance could void the user’s authority to operate the equipment.
9. This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada’s licence-exempt RSS(s). Operation is subject to the following two conditions: (1)This device may not cause interference. (2)This device must accept any interference, including interference that may cause undesired operation of the device.
10. This equipment complies with IC RSS-102 radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
11. Bluetooth® associated logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such

marks by Wyze Labs Inc. is permitted under license. Other trademarks and trade names are those of their respective owners.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## **CARE & MAINTENANCE**

1. Avoid high temperatures and irradiation.
2. If the scale is stored close to freezing point, allow it to acclimate to room temperature before use.
3. Remove the batteries if the scale is not going to be used for 3 months or more.
4. Clean the scale with a soft damp cloth.

Do not use abrasive or solvent-based cleaners, and never immerse it in water.

5. The scale will maintain its safety and performance features for at least 10,000 measurements.

## **WARRANTY INFORMATION**

Your Wyze product is covered by a one-year limited warranty. You can view the limited warranty terms at <https://wyze.com/return-warranty-policy> or request a copy by contacting Wyze at +1 (206) 339-9646.

## **FCC WARNING**

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- This device may not cause harmful interference, and
- This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### **FCC RADIATION EXPOSURE STATEMENT**

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

### **ISED RADIATION EXPOSURE STATEMENT**

This equipment complies with IC RSS-102 radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

### **ISED RSS WARNING**

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

### **SPECIFICATIONS**

Product Name: Wyze Scale X Black

Model: WSCALX

Product Size: 300 mm x 300 mm x 27.5 mm

Product Weight: 3.96 lbs / 1.80 kg

Power: 4 x 1.5V AAA batteries (included)

Display Units: lb/kg

Measurement Range: Weight: 11 – 400 lbs / 5 – 180 kg

Operating Temperature: 41°F – 104°F (5°C – 40°C)

Operating Humidity: 20 – 85% RH

Storage and Transport Humidity: ≤90% RH

Storage and Transport Temperature: -4°F – 131°F (-20°C – 55°C)

Apple logo and App Store are trademarks of Apple Inc., registered in the U.S. and other countries and regions.

Google Assistant, Google Home, Google Play, and the Google Play logo are trademarks of Google LLC.

Wyze and Wyze Scale X are trademarks of Wyze Labs, Inc.

### **NEED HELP?**

Contact our Support Team: <https://wyze.com/support>


Join the Community: <https://wyze.com/community>

**Scan the QR code to learn more**



<https://support.wyze.com/hc/en-us>



	<p><a href="#">WYZE WSCALX Scale X Digital Scale</a> [pdf] User Guide</p> <p>WSCALX, 2AUIUWSCALX, Scale X, Digital Scale, Scale X Digital Scale, WSCALX Scale X Digital Scale</p>
---	---

References

- [W Join Our Community – Wyze](#)
- [W Refund policy – Wyze](#)