

Xclusive Tongue Scraper Set User Manual

Home » Xclusive » Xclusive Tongue Scraper Set User Manual



Contents

- 1 Xclusive Tongue Scraper Set
- 2 How to
- **3 Package Content**
- 4 Documents / Resources
 - **4.1 References**
- **5 Related Posts**



Xclusive Tongue Scraper Set



How to

- 1. Place the curved tongue cleaner in the middle of your tongue and move it with light pressure to the tip of your tongue.
- 2. Rinse the tongue cleaner and repeat a few times until all deposits have disappeared.
- 3. Repeat with the straight tongue cleaner for both sides.
- 4. Repeat with the straight open tongue cleaner for the back of the tongue.
- 5. Clean and dry the tongue cleaners thoroughly after use and then store them in the case provided.
- 6. To maintain good hygiene, it is advisable to additionally clean the tongue cleaners in boiling water or with a little alcohol after a few uses.

Package Content

- 1 x Curved tongue cleaner (tongue surface)
- 1 x Straight tongue cleaner (tongue side)
- 1 x Straight open tongue cleaner (back)
- 1 x Deluxe storage case
- 1 x Dutch-language instructions for use

Thank you for your purchase and your trust in us. We select our products with the greatest care. If you have any questions about our product, please send us an email at bol@xclusive-lifestyle.nl
*Caution! Keep out of reach of small children.



Documents / Resources



Xclusive Tongue Scraper Set [pdf] User Manual Tongue Scraper Set, Tongue, Scraper Set

References

• LS_Lifestyle.nl - Let's make our house feel more like a home.

Manuals+,