

BSN Medical ULCERCARE 2-PART SYSTEM

BSN Medical ULCERCARE 2-Part System

Model: 114481

1. INTRODUCTION

The BSN Medical ULCERCARE 2-Part System is designed to provide effective compression therapy for the management of venous leg ulcers. This system consists of an inner liner and an outer compression bandage, working together to deliver sustained, therapeutic compression. Please read this manual carefully before use to ensure proper application and care.

2. SAFETY INFORMATION

- Always consult a healthcare professional before using this product, especially if you have pre-existing medical conditions.
- Do not apply the system over open wounds, broken skin, or infected areas without direct medical supervision.
- Ensure the system is applied correctly to avoid excessive pressure or discomfort. Improper application can lead to complications.
- Discontinue use and consult your healthcare professional if you experience increased pain, numbness, tingling, swelling, or skin discoloration.
- This product is for external use only.

3. PACKAGE CONTENTS

Upon opening the package, please verify that all components are present and undamaged:

- 1 x ULCERCARE Inner Liner
- 1 x ULCERCARE Outer Compression Bandage

4. SETUP AND APPLICATION

Follow these steps for proper application of the ULCERCARE 2-Part System:

4.1. Preparing the Limb

1. Ensure the leg is clean and dry.
2. If prescribed, apply any topical medications or dressings to the ulcer area before applying the compression system.

4.2. Applying the Inner Liner

1. Roll the inner liner down to the heel.
2. Place your foot into the liner, ensuring the heel pocket is correctly positioned.
3. Gently pull the liner up the leg, smoothing out any wrinkles. The liner should extend just below the knee.
4. Ensure the liner fits snugly but comfortably, without excessive tightness or bunching.



Image: The BSN Medical ULCERCARE 2-Part System, showing the inner liner and outer compression bandage. This image illustrates the product components before application.

4.3. Applying the Outer Compression Bandage

1. Starting from the base of the toes, wrap the outer compression bandage around the foot and ankle.
2. Continue wrapping upwards in a spiral motion, overlapping each layer by approximately 50%.
3. Maintain consistent, firm tension as you wrap, ensuring even compression from the ankle to just below the knee.
4. Secure the end of the bandage with the provided fasteners or tape, ensuring it is firm but not constricting.

5. WEARING INSTRUCTIONS

The ULCERCARE 2-Part System is typically worn for a prescribed period, often for several days, before being changed. Your healthcare professional will advise you on the specific duration and frequency of changes. During wear, monitor your comfort level and check for any signs of discomfort or adverse reactions as mentioned in the Safety Information section.

6. CARE AND MAINTENANCE

6.1. Cleaning Instructions

To maintain the integrity and hygiene of your ULCERCARE system:

- **Hand Wash Only:** Wash the inner liner and outer bandage by hand in lukewarm water with a mild detergent.
- Rinse thoroughly to remove all soap residue.
- Do not wring or twist. Gently squeeze out excess water.
- Lay flat to air dry away from direct heat or sunlight. Do not tumble dry, iron, or bleach.

6.2. Storage

Store the system in a clean, dry place at room temperature when not in use. Keep away from sharp objects that could damage the fabric.

7. TROUBLESHOOTING

If you encounter any issues while using the ULCERCARE 2-Part System, consider the following:

- **Discomfort or Excessive Tightness:** Re-evaluate the application. The system should be snug but not painful. If discomfort persists, remove the system and consult your healthcare professional.
- **System Slipping Down:** Ensure the inner liner is pulled up smoothly and the outer bandage is applied with consistent, firm tension, overlapping correctly.
- **Skin Irritation:** Check for proper skin preparation and ensure no wrinkles are present. If irritation continues, consult your healthcare professional.

8. SPECIFICATIONS

Attribute	Detail
Brand	BSN Medical
Model Number	114481
ASIN	B00VS44V98
Unit Count	1.00 Count
Product Care Instructions	Hand Wash Only

9. WARRANTY AND SUPPORT

For specific warranty information or product support, please contact BSN Medical directly through their official website or customer service channels. Keep your purchase receipt as proof of purchase.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question