[Manuals.plus](#) /

› [TUMAKOU](#) /

› TUMAKOU RF Radiofrequency Device User Manual (Model RF01)

TUMAKOU RF01

TUMAKOU RF Radiofrequency Device User Manual (Model RF01)

Model: RF01 | Brand: TUMAKOU

INTRODUCTION

The TUMAKOU RF Radiofrequency Device is designed for at-home skin care, utilizing radiofrequency technology to promote collagen regeneration, smooth wrinkles, and improve skin elasticity. This device offers a convenient solution for facial and body skin rejuvenation, helping to reduce signs of aging and enhance overall skin appearance.

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Image: The TUMAKOU RF Radiofrequency Device, featuring its main unit and two interchangeable probes for different treatment areas.

WHAT'S IN THE BOX

Please check the package contents upon receipt to ensure all items are present:

- Main Unit
- Facial Probe
- Body Probe
- Power Cord
- User Manual

Image: All components included in the TUMAKOU RF device package.

DEVICE OVERVIEW

The TUMAKOU RF device features a main unit with a control panel and two interchangeable probes. The control panel allows for adjustment of energy levels and treatment time.

Probes

- **Facial Probe:** Designed for delicate areas such as the face and around the eyes.
- **Body Probe:** Suitable for larger body areas including the neck, abdomen, arms, hands, thighs, and buttocks.

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Image: Illustration of the two distinct probes for comprehensive body and facial treatment.

Control Panel

The device features an intuitive control panel with:

- **Energy Level Selector:** Three adjustable energy levels (1, 2, 3) to customize treatment intensity.
- **LCD Time Display:** Shows remaining treatment time (10, 20, or 30 minutes).

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Image: Detailed view of the device's control panel, highlighting energy level and time settings.

SETUP

1. **Install the Probe:** Select either the facial or body probe and securely attach it to the main unit. Ensure it clicks into place.

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Image: Step 1 - Installing the desired probe onto the device.

2. **Connect Power:** Plug the white handle cord and the black power cord into their respective ports on

the machine. Then, plug the power cord into a suitable electrical outlet.

Comment Utiliser

Image: Step 2 - Connecting the power and handle cords.

OPERATING INSTRUCTIONS

Preparation

- Cleanse the skin thoroughly before treatment.
- Apply a conductive gel or serum evenly to the treatment area. This is crucial for effective and comfortable use.

Treatment Steps

1. **Power On:** Press the power button to turn on the device.
2. **Select Treatment Time:** Press the 'TIME' button to choose a treatment duration of 10, 20, or 30 minutes. The selected time will be displayed on the LCD screen.

Comment Utiliser

Image: Step 3 - Selecting the desired treatment time.

3. **Select Energy Level:** Press the 'LEVEL' button to choose an appropriate energy level (1, 2, or 3) for your skin. Start with Level 1 and gradually increase if comfortable.

Comment Utiliser

Image: Step 4 - Selecting the appropriate energy level.

4. **Begin Treatment:** Gently move the probe in circular motions over the treated area. Maintain constant contact with the skin.
5. **Completion:** The device will automatically stop after the selected time. Turn off the device and unplug it.

Recommended Usage

- Use the device once every 1-2 days.
- A full treatment cycle consists of 10 sessions.
- Visible improvements may be observed after 6-12 uses. Consistent use is key for optimal results.

Expected Results

Individual results may vary. Generally, visible skin improvement can be seen after 6 to 10 weeks of consistent use (10 to 30 minutes per treatment, once every 2-3 days). Continued use can lead to longer-lasting effects.

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
Image: Before and after comparison illustrating potential skin improvements.

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Image: Diagram explaining the bipolar RF technology and its effect on skin layers.

MAINTENANCE

- **Cleaning:** After each use, wipe the probe and main unit with a damp cloth to remove any gel residue. Do not use abrasive cleaners or immerse the device in water.

- **Storage:** Store the device in a cool, dry place, away from direct sunlight and moisture. Keep it out of reach of children.

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Image: A summary of useful tips for optimal device performance and user experience.

1. Always apply conductive gel during use for best results.
2. If you experience overheating or discomfort, choose a lower energy level.
3. After use, clean the device with a damp cloth.
4. It is normal for the heat sensation to decrease after a period of use, as your skin adapts. The effect is still present.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not heat up or feels cold.	Insufficient conductive gel; incorrect energy level; device malfunction.	Ensure sufficient conductive gel is applied. Increase energy level if comfortable. If problem persists, contact customer support.
No visible results after several uses.	Inconsistent use; insufficient treatment duration/intensity; individual skin response.	Ensure consistent use (once every 1-2 days). Use for recommended duration (10-30 mins) and highest comfortable energy level. Visible results may take 6-12 uses.
Skin feels too hot or uncomfortable.	Energy level is too high for your skin sensitivity.	Reduce the energy level to a comfortable setting. Ensure adequate conductive gel is applied.

SPECIFICATIONS

Feature	Detail
Model Number	RF01
Dimensions (L x W x H)	21.5 x 21.5 x 13 cm (8.46 x 8.46 x 5.12 inches)

Feature	Detail
Item Weight	2.15 Kilograms (4.74 lbs)
Material	Metal
Power Source	Corded Electric
Applicable Skin Type	All skin types
Included Components	Facial probe, Body probe
GTIN	794587324337

WARRANTY AND SUPPORT

The TUMAKOU RF Radiofrequency Device comes with a **1-year warranty** from the date of purchase. This warranty covers manufacturing defects and ensures a risk-free experience.

For any questions, issues, or support needs regarding your device, please contact TUMAKOU customer service. Refer to the contact information provided in your product packaging or on the official TUMAKOU website.

We are committed to providing high-quality products and excellent customer service.