

[Manuals.plus](#) /

› [CurrentBody](#) /

› CurrentBody Skin LED Hair Regrowth Device User Manual

CurrentBody LED Hair Regrowth Device

CurrentBody Skin LED Hair Regrowth Device User Manual

Model: LED Hair Regrowth Device

1. INTRODUCTION

The CurrentBody Skin LED Hair Regrowth Device is an advanced at-home solution designed to address hair thinning and promote new hair growth. Utilizing clinically proven red LED light therapy, this device targets the scalp to stimulate hair follicles and improve overall hair health. It is FDA-cleared and offers a non-invasive, natural, and pain-free approach to hair regrowth.



Image: A woman wearing the CurrentBody Skin LED Hair Regrowth Device, demonstrating its comfortable fit and the illuminated red lights.

2. SAFETY INFORMATION

Please read all safety information carefully before using the CurrentBody Skin LED Hair Regrowth Device. Failure to follow these instructions may result in injury or device malfunction.

- **Consult your doctor:** Prior to using this device, consult your doctor if you are taking prescription or over-the-counter drugs known to cause light sensitivity.
- **Skin conditions:** Do not use if you have any suspicious or cancerous lesions on your scalp.
- **Eye conditions:** Do not use if you have any photosensitive disorders or any genetic conditions of the eye.
- **Potential risks:** Immediately after treatment, skin may appear red. Skin should return to normal within 1 to 2 hours. Temporary headaches may immediately follow the light treatment. Headaches can be treated with an appropriate painkiller and rehydration.
- **Intended Use:** The CurrentBody Skin LED Hair Regrowth Device is indicated to treat androgenetic alopecia and promote hair growth in males with Norwood-Hamilton Classifications of IIa to V patterns of hair loss, and in females with Ludwig-Savin Scale I-1 to I-4, II-1, II-2, or frontal patterns of hair loss,

both with Fitzpatrick Skin Types I to IV.

3. PRODUCT OVERVIEW

The CurrentBody Skin LED Hair Regrowth Device utilizes specific red light wavelengths (620nm - 660nm) to penetrate the scalp and stimulate hair follicles. This process helps to improve scalp health, increase circulation, and reactivate dormant follicles for thicker, healthier hair.



Image: Front view of the device with red lights.



Image: Close-up of the device's red LED lights.

The device features 120 medical-grade LED bulbs, delivering 18 J/cm² of precise red light. It is designed for full head coverage and includes an integrated display to show treatment time remaining.

Why 620 – 660nm?




Image: Infographic detailing the hair concerns targeted by LED technology.

4. SETUP

Follow these steps to set up your CurrentBody Skin LED Hair Regrowth Device:

1. **Unboxing:** Carefully remove the device and its accessories from the packaging. Ensure all components are present: the LED Hair Regrowth Device, charging stand, and power adapter.
2. **Initial Charge:** Place the device on its charging stand. Connect the power adapter to the charging stand and plug it into a power outlet. Allow the device to fully charge before first use. The indicator light on the device will show charging status.
3. **Familiarization:** Take a moment to familiarize yourself with the device's design, including the internal LED array and the control panel.



Image: The device partially unboxed, showing its initial presentation.



Image: The device placed on its dedicated charging stand.


5. OPERATING INSTRUCTIONS

For optimal results, use the CurrentBody Skin LED Hair Regrowth Device daily for 10 minutes.

1. **Preparation:** Ensure your hair is clean and dry before use.
2. **Placement:** Gently place the device on your head, ensuring it sits comfortably and covers your entire scalp.
3. **Activation:** Press the power button to start the treatment. The red LED lights will illuminate, and the integrated display will show the remaining treatment time.
4. **Treatment:** Relax during the 10-minute session. The device is designed to be comfortable for daily wear.
5. **Automatic Shut-off:** The device will automatically turn off after 10 minutes, indicating the completion of your treatment.
6. **Post-Treatment:** Remove the device and place it back on its charging stand for storage and to ensure it's ready for the next use.



Image: A user wearing the device during a treatment session.



Most people treat hair concerns at the ends. But real results start at the root.

Red light LED therapy works beneath the surface to improve scalp health by hydrating, calming redness and rebalancing oil.

This creates the ideal environment for visibly thicker, stronger hair.

Image: Highlighting the device's focus on scalp health for hair regrowth.

Consistent daily use is key to achieving visible results, which typically appear within 16 weeks, with significant improvements in hair thickness and reduced shedding observed in clinical studies.

You'll see the difference:

Healthier scalp

Thicker, fuller hair in 12 weeks

Less grease & redness



Image: Expected benefits from consistent use, including healthier scalp and thicker hair.

6. MAINTENANCE

Proper care and maintenance will ensure the longevity and effectiveness of your device.

- **Cleaning:** After each use, gently wipe the interior and exterior surfaces of the device with a soft, dry, or slightly damp cloth. Do not use abrasive cleaners, alcohol, or harsh chemicals. Ensure the device is completely dry before storing or charging.
- **Storage:** Store the device on its charging stand in a cool, dry place, away from direct sunlight and extreme temperatures. Keep it out of reach of children.
- **Charging:** Always use the provided charging stand and power adapter. Do not attempt to charge the device with incompatible chargers.

7. TROUBLESHOOTING

If you encounter any issues with your CurrentBody Skin LED Hair Regrowth Device, please refer to the following troubleshooting guide:

Problem	Possible Cause	Solution
Device does not turn on	Low battery or not charged.	Place the device on the charging stand and ensure it is properly connected to power. Allow it to charge fully.
LED lights are dim or flickering	Low battery.	Recharge the device. If the problem persists after charging, contact customer support.
Device shuts off unexpectedly during treatment	Low battery or treatment cycle completed.	Check the battery level. If not low, the 10-minute treatment cycle may have finished. Recharge if necessary.

Problem	Possible Cause	Solution
Discomfort during use	Improper placement or sensitivity.	Adjust the device for a comfortable fit. If redness or irritation occurs, discontinue use and consult a healthcare professional.

If the issue is not resolved by the above steps, please contact CurrentBody customer support for further assistance.

8. SPECIFICATIONS

- **Model:** LED Hair Regrowth Device
- **Brand:** CurrentBody
- **LEDs:** 120 medical-grade LED bulbs
- **Wavelength:** Red light (620nm - 660nm)
- **Energy Output:** 18 J/cm²
- **Treatment Time:** 10 minutes daily
- **Regulatory:** FDA-cleared
- **UPC/GTIN:** 731093153136
- **ASIN:** B0CQRK9FH3

9. WARRANTY AND SUPPORT

The CurrentBody Skin LED Hair Regrowth Device comes with a manufacturer's warranty. Please refer to the warranty card included in your product packaging for specific terms and conditions, including coverage duration and claim procedures.

For product support, technical assistance, or warranty claims, please visit the official CurrentBody website or contact their customer service department. Contact details can typically be found on the CurrentBody website or within your product documentation.

Online Support: www.currentbody.com