

Aopvui FZ201-Curve

Aopvui IPL Permanent Hair Removal Device User Manual

Model: FZ201-Curve

INTRODUCTION

Thank you for choosing the Aopvui IPL Permanent Hair Removal Device. This manual provides essential information for the safe and effective use of your device. Please read it thoroughly before first use and keep it for future reference. This device utilizes advanced IPL (Intense Pulsed Light) technology to target hair follicles, helping to prevent hair regrowth and achieve long-lasting smooth skin.

IMPORTANT SAFETY INFORMATION

Your safety is paramount. Please adhere to the following guidelines:

- Always wear the provided protective glasses during operation to shield your eyes from the intense light flashes.
- Do not use the device on tattoos, dark spots, moles, or sunburned skin.
- Avoid using the device if you have a skin disease, are pregnant, or are undergoing medical treatment that makes your skin sensitive to light.
- This device is not suitable for all skin tones and hair colors. Refer to the compatibility chart below to ensure it is safe and effective for you.
- Keep the device out of reach of children.
- Ensure the device is clean and dry before each use.
- This device is not chargeable; it requires continuous power from an electrical outlet during use.

Check if this Device is Right for you?

SKIN TONES



Very Fair



Fair



White-Olive



Light Brown



Brown



Dark

HAIR COLORS



Black



Brown



Dark Blonde



Light Blonde



Red



Grey



Safe and Suitable



Not Suitable

Figure 1: Skin Tone and Hair Color Compatibility. This chart indicates which skin tones and hair colors are safe and suitable for IPL treatment. Green checkmarks denote suitability, while red crosses indicate it is not suitable. Generally, lighter skin tones and darker hair colors yield the best results.



Figure 2: Usage Steps and Contraindications. This image illustrates the six steps for using the device, including shaving, connecting power, wearing glasses, adjusting energy, placing on skin, and moisturizing. It also clearly lists conditions where the device should not be used, such as on tattoos, skin diseases, sunburn, mustache area (for men), and by pregnant women.

PRODUCT OVERVIEW

The Aopvui IPL device is designed for convenient and effective at-home hair removal. It features a durable design with a high flash count for long-term use.



Figure 3: Device and Included Accessories. The image displays the Aopvui IPL device, a pair of protective glasses, a razor, and the power adapter, all included in the package.

Key Features:

- **IPL Technology:** Targets hair follicles to inhibit regrowth.
- **9 Energy Levels:** Adjustable intensity from 1 to 9 (up to 15J) to suit different skin sensitivities and treatment areas.
- **2 Modes:** Automatic mode for large areas and Manual mode for precise treatment on smaller, sensitive areas.
- **High Flash Count:** Over 999,900 flashes, providing a long lifespan for the device.
- **Narrow Flash Window:** Designed for efficient hair removal with less missing area.
- **Integrated Safety Features:** Includes eye protection and a smart sensor for safe operation.

600-1200nm Red Light Tech

**Targets hair root precisely
without damaging your skin.**

A Wavelength of 600nm Above Is Absorbed
Selectively By The Melanin of Hair Follicles.




Figure 4: Red Light Technology. This diagram illustrates how the 600-1200nm red light technology precisely targets the hair root without damaging the skin, as the wavelength is selectively absorbed by the melanin in hair follicles.

Narrower Flash Window Frame

**More efficient hair removal
with less missing hair removal area.**

IPL devices with larger flash window
frame may miss some skin areas, resulting
in incomplete hair removal session.



Figure 5: Narrow Flash Window Frame. The image highlights the narrow flash window frame of the device, which contributes to more efficient hair removal by minimizing missed areas during treatment.

SETUP

1. **Prepare the Skin:** Ensure the treatment area is clean, dry, and free of any lotions or creams. Shave off any visible hair from the area you intend to treat.
2. **Connect Power:** Plug the power adapter into the device and then into a suitable electrical outlet. The device is not battery-operated and must be plugged in during use.
3. **Wear Protective Glasses:** Always put on the provided protective glasses before turning on the device to protect your eyes from the intense light.

ATTENTION! IMPORTANT NOTICE:



This Device is not chargeable because it has no battery.



This Device requires Electricity to power on and should be plugged into electric socket during use.

Figure 6: Power Connection. This image demonstrates the correct way to connect the power adapter to the device and plug it into an electrical socket, emphasizing that the device requires continuous power.

OPERATING INSTRUCTIONS

1. **Power On:** Press and hold the power button to turn on the device.
2. **Adjust Energy Level:** Start with the lowest energy level (Level 1) for your first treatment. Press the power button briefly to cycle through the 9 energy levels. Increase the level gradually as your skin adapts and tolerates the intensity.
3. **Select Mode:**
 - **Automatic Mode:** For large areas like legs, arms, chest, and back. Press and hold the flash button until the device enters automatic mode. The device will flash continuously when placed on the skin.
 - **Manual Mode:** For small or sensitive areas such as the upper lip, bikini line, armpits, and face. Press the flash button once for each flash.
4. **Apply to Skin:** Place the flash window firmly and completely against your skin. The device will only flash when full contact is made.
5. **Treat the Area:** Move the device slowly across the treatment area, ensuring each section is covered once. Avoid overlapping flashes excessively.
6. **Power Off:** After completing your session, press and hold the power button to turn off the device.

9 ENERGY LEVELS **13.1J**

Always start from a lower level



Figure 7: Energy Level Adjustment. This image illustrates the 9 adjustable energy levels of the device, ranging from Level 1 (7J) to Level 9 (13.1J), advising users to always start from a lower level.

AUTO MODE

Suitable for Large areas such as



Chest



Arms



Legs



Back



MANUAL MODE

Suitable for small areas such as



Upper Lip



Bikini Line



Armpit



Face




Figure 8: Auto and Manual Modes. This image demonstrates the application of the device in Auto Mode for large body areas like the chest, arms, legs, and back, and Manual Mode for smaller, sensitive areas such as the upper lip, bikini line, armpit, and face.

RECOMMENDED USAGE SCHEDULE

Consistent use is key to achieving optimal results. Follow this general schedule:

Phase	Frequency	Expected Outcome
Weeks 1-2	3 times per week	Initial reduction in hair growth.
Weeks 3-4	2 times per week	Noticeable slowing of hair regrowth.
Weeks 5-8	1-2 times per week	Significant hair reduction, up to 94% for some users.
After 8 Weeks (Maintenance)	Once every 2 months	Long-lasting smooth skin.

Individual results may vary. Clinical tests show that 80% of users experience 94% hair reduction after 7-8 treatments.

HAIR REMOVAL PROCESS

Amazing results within 8weeks




Figure 9: Hair Removal Process. This visual timeline illustrates the progressive reduction in hair growth over 8 weeks of consistent use, leading to complete hair removal and smooth skin.

VISIBLE RESULTS WITH LESS TIME

Before



VS

After



Figure 10: Visible Results. This image provides a visual comparison of an armpit before and after using the IPL device, demonstrating the visible reduction in hair and smoother skin achieved with consistent treatment.

MAINTENANCE

- **Cleaning:** After each use, ensure the device is unplugged and cool. Gently wipe the flash window and the device body with a soft, dry cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place, away from direct sunlight and moisture. Keep it in its original

packaging or a protective case to prevent damage.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Not properly plugged in; power outlet issue.	Ensure the power adapter is securely connected to the device and the electrical outlet. Try a different outlet.
Device does not flash.	Flash window not in full contact with skin; device in manual mode and flash button not pressed.	Ensure the flash window is pressed firmly and completely against the skin. If in manual mode, press the flash button.
Skin feels too hot or uncomfortable.	Energy level is too high for your skin sensitivity.	Reduce the energy level. Always start with the lowest setting and increase gradually. Ensure skin is clean and dry.
No visible results after several sessions.	Inconsistent use; incorrect energy level; unsuitable skin/hair type.	Ensure consistent use as per the recommended schedule. Increase energy level if comfortable. Verify your skin tone and hair color are suitable using the compatibility chart.

SPECIFICATIONS

- **Model:** FZ201-Curve
- **Technology:** IPL (Intense Pulsed Light)
- **Flashes:** 999,900+
- **Energy Levels:** 9 (up to 15J)
- **Modes:** 2 (Automatic, Manual)
- **Wavelength:** 600-1200nm
- **Power Source:** Corded Electric
- **Package Dimensions:** 11.02 x 7.32 x 2.28 inches
- **Weight:** 1.76 Pounds
- **Manufacturer:** Aopvui

WARRANTY AND SUPPORT

Aopvui is committed to providing high-quality products and customer satisfaction. Your purchase includes the IPL laser hair removal device, a razor, a pair of glasses, and a power adapter. If any accessories are missing or if you encounter any issues with your device, please contact Aopvui customer support for assistance.

For after-sales support, please refer to the contact information provided with your product packaging or visit the official Aopvui store on Amazon: [Aopvui Store](#).

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

My Aopvui FZ201-Curve device won't turn on and flashes the error codes E3 and E4. What does that mean and how can I fix it?

Based on the available manual for the Aopvui FZ201-Curve, specific alphanumeric error codes like E3 and E4 are not explicitly defined in the standard troubleshooting table. However, when an IPL device fails to turn...

[Read full answer](#)