

[Manuals.plus](#) /

› [NVBOTY](#) /

› NVBOTY LED Face Mask Light Therapy User Manual

NVBOTY 2344

NVBOTY LED Face Mask Light Therapy User Manual

Model: 2344

Brand: NVBOTY

INTRODUCTION

This manual provides instructions for the safe and effective use of the NVBOTY LED Face Mask Light Therapy device. Please read this manual thoroughly before first use and retain it for future reference.



Image: The NVBOTY LED Face Mask Light Therapy device, remote control, and eye shields. The mask is shown on a person's face, illustrating its application.

IMPORTANT SAFETY INFORMATION


- Ensure the USB plug of the LED Face Mask Light Therapy is fully inserted into the USB-A port of the remote control. Improper connection may prevent the device from starting or cause some lights to malfunction.
- Do not use the red light therapy for face while charging the remote control.
- Statements regarding light therapy devices have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. Consult a healthcare professional for any medical concerns.
- Perform a skin sensitivity test before the first use.
- Avoid direct eye exposure to the LED lights. Use the provided eye shields.
- Keep out of reach of children.
- Do not use if you have any pre-existing skin conditions that may be exacerbated by light therapy.

PACKAGE CONTENTS

- NVBOTY LED Face Mask Light Therapy Device
- Rechargeable Remote Control
- Dedicated Eye Shields
- USB Charging Cable (for remote)
- User Manual (this document)

SETUP

1. Step 1: Prepare Your Skin



TIP: Perform a skin sensitivity test before the first use.



Step 1

Clean and dry your face and hands ensuring there is no makeup or dirt.



Step 2

Fully connect the light therapy face mask plug to the remote control.



Step 3

Place the mask on your face and adjust it to the desired angle.



Step 4

Select the mode and desired time, then start the session.



Step 5

Remove the LED light therapy mask and apply skincare products to your skin.

The best recommended usage method for the spectral mask instrument is as follows: For the first month use it once daily. Starting from the second month, use it 3-4 times a week, with each session lasting 20 minutes.

Image: Step 1 of the usage guide, showing a woman cleaning her face. This image also contains steps 2-5.

Clean and dry your face and hands thoroughly, ensuring no makeup or dirt remains on your skin.

2. Step 2: Connect the Device

Fully connect the light therapy face mask plug to the USB-A port of the remote control. Ensure a secure connection.

3. Step 3: Position the Mask

Place the mask on your face and adjust it to a comfortable and desired angle. The soft, skin-friendly silicone material is designed to conform to your face. Ensure the dedicated eye shields are properly positioned to block light from entering your eyes.

OPERATING INSTRUCTIONS

The NVBOTY LED Face Mask Light Therapy features a rechargeable remote control for easy operation.

Powering On/Off

- Press and hold the power button on the remote control to turn the device on or off.

Selecting a Mode




Image: Visual representation of the four light therapy modes available on the device, detailing their primary functions and light colors. The device offers 4 distinct modes:

- **Infrared (850nm) & Red Light (630nm) Mode:** Intended to smooth skin, enhance elasticity, and firm the face.
- **Orange Light (600nm) Mode:** Intended to reduce pigmentation and soothe redness.
- **Blue Light (460nm) Mode:** Aims to balance oil production for clearer skin.
- **4-Color Cycle Therapy Mode:** Cycles through all available light colors for comprehensive treatment.

Press the "Mode Button" on the remote control to cycle through the available light therapy modes.

Setting the Timer

Timer Function for Safe Sleeping

10Min / 15Min / 20Min(recommended) / 25Min / 30Min

Set a Time



Remote Control



Image: A user wearing the LED face mask while resting, demonstrating the timer function for convenient use. The remote control is shown with timer options.

The system defaults to a 10-minute timer. Short presses of the "Time Button" on the remote allow you to adjust the session duration to 15, 20, 25, or 30 minutes. The device will automatically shut off when the timer ends.

For optimal results, use the device 3-4 times a week, with each session lasting 20 minutes.

Post-Session Care

After your session, remove the LED light therapy mask and apply your preferred skincare products to your skin.

MAINTENANCE

- **Cleaning:** Gently wipe the mask and remote control with a soft, damp cloth after each use. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place away from direct sunlight and extreme temperatures.
- **Charging:** The remote control has a 2000mAh rechargeable battery. Connect the remote to a suitable USB power

source using the provided charging cable when the battery is low.

TROUBLESHOOTING

Device not starting or some lights not turning on:

Ensure that the USB plug of the LED Face Mask Light Therapy is fully inserted into the USB-A port of the remote control. A poor connection can prevent proper operation. Verify the remote control battery is charged.

Remote control not responding:


Ensure the remote control battery is charged. If the issue persists, try restarting the remote control by powering it off and on.

Irregular light patterns:

Check for any physical damage to the mask or connections. If the issue persists, contact customer support.

SPECIFICATIONS

Feature	Detail
Model Number	2344
Brand	NVBOTY
LED Count	400 LEDs
Wavelengths	Infrared (850nm), Red Light (630nm), Orange Light (600nm), Blue Light (460nm)
Remote Battery	2000mAh Lithium Ion (included)
Timer Settings	10, 15, 20, 25, 30 minutes
Material	Soft, skin-friendly silicone
Package Dimensions	15.94 x 10.47 x 2.01 inches
Weight	2.14 Pounds



Better Penetration with Near-Infrared Light Technology

The near-red light wave penetrates the dermis, a layer not reachable by skincare products or other types of light. This process tightens the skin, stimulates cell vitality in the dermis and provides excellent anti-aging effects. It also enhances microcirculation, slows down skin aging and the formation of fine lines, and gets rid of excess subcutaneous fat.

Image: A cross-section diagram of skin layers (Epidermis, Dermis, Hypodermis) showing the varying penetration depths of blue, orange,

red, and infrared light wavelengths.

Four All-in-One Light Sources

More Even Distribution, Wider Coverage

Skin irradiation area increased by **1.47** times




Image: An illustration highlighting the device's four-in-one light sources, emphasizing even distribution and wider coverage for skin irradiation.

WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the product packaging or contact NVBOTY customer service directly. Contact details can typically be found on the manufacturer's website or through your purchase platform.

Manufacturer: NVBOTY

© 2024 NVBOTY. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question