

JOVS Blacken

JOVS Blacken DPL Facial Skincare Device User Manual

1. INTRODUCTION

The JOVS Blacken DPL Facial Skincare Device is designed for at-home skin rejuvenation, targeting concerns such as dark spots, blemishes, and uneven skin tone. Utilizing a unique 500-650nm narrow-spectrum light, this device aims to improve skin texture and firmness. It features a sapphire cooling surface for comfort and a 3-layer safety system for user protection. This manual provides comprehensive instructions for safe and effective use of your device.



Image 1.1: The JOVS Blacken DPL Facial Skincare Device, protective goggles, and packaging.

2. SAFETY INFORMATION

Before using the device, please read this user manual carefully, including all contraindications, warnings, and safety information. Failure to follow these instructions may result in injury or device malfunction.

- Always wear the provided protective goggles during operation to shield your eyes from the light.
- Do not use the device on broken, irritated, or infected skin.
- Avoid using the device on areas with tattoos, permanent makeup, or moles.
- Consult a healthcare professional before use if you have any pre-existing skin conditions, are pregnant, breastfeeding, or taking photosensitive medications.
- Keep the device out of reach of children.
- Ensure the device is clean and dry before each use.
- Do not immerse the device in water or any other liquid.

3. WHAT'S IN THE BOX

Your JOVS Blacken DPL Facial Skincare Device package includes the following items:

- 1 x JOVS Blacken Pro Device
- 1 x US Power Plug
- 1 x Pair of Protective Goggles
- 1 x Detailed User Guide (this manual)

In the box



JOVS Blacken



Us Plug



Goggles



User Guide

Image 3.1: All items included in the JOVS Blacken DPL Facial Skincare Device package.

4. DEVICE OVERVIEW

The JOVS Blacken device is ergonomically designed for comfortable handling. It features a sapphire cooling surface that provides immediate skin cooling during use, enhancing comfort and protecting the skin from heat. The device utilizes a precise narrow-spectrum light (500-650nm) to target specific skin concerns.



Sapphire cooling surface
pain-free



Narrow-spectrum light
of 500-650nm



Intelligent continuous flash
makes for easier handling



4 weeks treatment time to
brighten spot reducing



Image 4.1: Key features of the JOVS Blacken device, including sapphire cooling and narrow-spectrum light technology.

The device offers three distinct modes to address various skincare needs:

- **BRI./PORE Mode:** Designed for full-face rejuvenation, this mode helps refine pores, improve overall skin texture, and enhance radiance.
- **REDNESS Mode:** Focuses on reducing facial redness and irritation, promoting a more balanced skin tone.
- **POWER Mode:** Provides targeted care for stubborn spots and pigmentation, aiming to reduce visible blemishes.




Image 4.2: Illustration of the three operational modes and their targeted skin areas.

5. SETUP

1. **Cleanse Your Face:** Before each treatment, ensure your face is clean and free of makeup, lotions, or oils.
2. **Connect Power:** Plug the provided US power plug into the device and then into a power outlet.
3. **Power On:** Press and hold the start button for 2 seconds to turn on the device.

How To Use



Image 5.1: Visual guide for the initial setup and usage steps.

6. OPERATING INSTRUCTIONS

Follow these steps for effective use of your JOVS Blacken DPL Facial Skincare Device:

1. **Select Mode:** Short press the start button to cycle through and select your desired treatment mode (BRI./PORE, REDNESS, or POWER).
2. **Wear Goggles:** Put on the included protective goggles to protect your eyes from the intense light.
3. **Apply to Skin:** Place the device head gently against your skin. The device will automatically detect skin contact.
4. **Initiate Flash:** Short press the flash button to emit a single flash. For continuous flashing, hold the device against your skin and glide it slowly across the treatment area.
5. **Treatment Frequency:**

- **Beginners:** Start with 4-6 flashes per area, 3-4 times per week.
- **Adapted Skin:** Once your skin adapts, you may increase to 8-10 flashes per area, potentially daily if no irritation occurs.

6. **Post-Treatment Care:** After completing your session, apply a hydrating mask or moisturizer to soothe and nourish your skin.

Important Note: If you experience dryness or irritation, reduce both the frequency of use and the number of flashes per area.



Image 6.1: Proper application of the JOVS Blacken device on the face.

Official Product Video: JOVS Blacken DPL Facial Device

Your browser does not support the video tag.

Video 6.1: This video demonstrates the JOVS Blacken DPL Facial Device in use, showcasing its design and application for various skin treatments. It provides a visual guide to the device's operation and features.

7. MAINTENANCE

Proper maintenance ensures the longevity and effectiveness of your device:

- **Cleaning:** After each use, gently wipe the device head with a soft, dry cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place, away from direct sunlight and moisture. Keep it in its original packaging or a protective case when not in use.
- **Power Cord:** Inspect the power cord regularly for any signs of damage. If damaged, discontinue use and contact customer support.

8. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Not properly connected to power, or power outlet issue.	Ensure the power plug is securely connected to both the device and a working power outlet. Try a different outlet.
No flashes emitted.	Device not in proper contact with skin, or incorrect mode selected.	Ensure the device head is fully in contact with the skin. Verify the correct mode is selected.
Skin irritation or discomfort.	Excessive use, sensitive skin, or improper post-treatment care.	Reduce frequency of use and number of flashes. Ensure sapphire cooling is active. Apply hydrating mask after treatment. If irritation persists, discontinue use and consult a dermatologist.
No visible results after 4 weeks.	Inconsistent use, or individual skin response varies.	Ensure consistent use as recommended (2-3 times weekly). Results may vary; continue use for longer periods as advised in the operating instructions. Ensure proper mode selection for your concerns.

9. SPECIFICATIONS

- **Model:** Blacken
- **Brand:** JOVS
- **Light Spectrum:** 500-650nm (Narrow-spectrum DPL)
- **Cooling Technology:** Sapphire Cooling Surface
- **Power Source:** US Power Plug (included)
- **Dimensions:** 13.35 x 9.13 x 4.37 inches (Package)
- **Weight:** 9.03 ounces (Package)
- **Manufacturer:** JOVS
- **ASIN:** B0DJK4YRLR

10. WARRANTY AND SUPPORT

For warranty information, product support, or any inquiries not covered in this manual, please contact JOVS customer service directly. Refer to the contact information provided on the product packaging or the official JOVS website.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question