

[Manuals.plus](#) /

› [SUNLUXY](#) /

› SUNLUXY BCE1 3-in-1 Body Sculpting Machine User Manual

SUNLUXY BCE1

SUNLUXY BCE1 3-in-1 Body Sculpting Machine User Manual

Model: BCE1

IMPORTANT SAFETY INFORMATION

Please read all instructions carefully before using the SUNLUXY BCE1 Body Sculpting Machine. Retain this manual for future reference.

- Do not use on broken, irritated, or infected skin.
- Consult a physician before use if you have any medical conditions, including but not limited to heart conditions, skin disorders, or are pregnant.
- Keep out of reach of children.
- Do not immerse the device in water. Clean with a damp cloth only.
- Do not use for more than 3 minutes per area, with a maximum of 10 minutes total per daily session.
- Discontinue use if you experience excessive discomfort, pain, or adverse reactions.
- The device automatically powers off after 10 minutes for safety.

PRODUCT OVERVIEW

The SUNLUXY BCE1 is a portable 3-in-1 body sculpting machine designed for skin tightening, cupping therapy, and red light treatment. It features adjustable intensity levels and interchangeable suction cups for various body areas.

Why Choose Sunluxy?

SunLuxy



- ✓ **3-IN-1 Massage**
Red Light + Cupping + Lifting
- ✓ **2 Modes / 3 Levels**
Adjustable, Comfortable & Effective
- ✓ **Double Cup**
Suitable for multiple areas
- ✓ **Safe & Easy**
10-Min Auto Shut-Off

Other



- ✗ **Limited functionality**
Does not accommodate diverse usage scenarios
- ✗ **Reduced compatibility**
Cannot be applied to multiple body parts
- ✗ **No power failure protection**
Insufficient safety safeguards
- ✗ **Bulky and cable-dependent**
Not suitable for quick portability

Image: The SUNLUXY BCE1 3-in-1 Body Sculpting Machine, a white handheld device with a clear cupping attachment and red light elements.

SETUP

1. Charging the Device

Before first use, fully charge the device. Connect the provided Type-C cable to the charging port located at the bottom of the device and to a suitable USB power adapter (not included). The charging indicator light will show the charging status.

Simple & Quick




01 Apply essential oil or cream evenly to the target area.



02 Select the suction cup (large cup for abdomen, thighs, and buttocks; small cup for arms and



03 Choose Cupping or Lifting mode and adjust to a comfortable intensity.



04 Place the massager on the target area for 3 minutes.

Image: Diagram showing the device's components, including the Type-C charging port at the bottom.

2. Attaching Suction Cups

The device comes with two interchangeable suction cups: a large cup for areas like the abdomen, thighs, and hips, and a small cup for arms and more precise areas. To attach a cup, align it with the device head and twist clockwise until it clicks securely into place. To remove, twist anti-clockwise.



Image: The SUNLUXY BCE1 device shown with both large and small interchangeable suction cups.

3. Inserting Filters

Place a filter into the center of the device head before attaching a suction cup. This helps maintain hygiene and device performance.

OPERATING INSTRUCTIONS

1. Power On/Off

Press and hold the Power button (located in the center of the device) to turn the device on or off.

2. Mode Selection

The device has two main modes: Body and Face. Use the designated button to switch between these modes. In Face Mode, the cupping function is disabled for gentle care.

3. Function and Intensity Adjustment

The device offers Cupping and Lifting functions, each with three adjustable intensity levels. Press the respective function buttons to cycle through the intensity levels (Level 1, Level 2, Level 3) or to turn the function off.

4. Application Steps

1. **Prepare Skin:** Apply essential oil, cream, or body lotion evenly to the target area to ensure smooth movement of the device and enhance absorption.
2. **Select Cup:** Choose the appropriate suction cup (large for body, small for arms/face) and attach it.
3. **Select Mode:** Power on the device, select Body or Face mode, and choose either Cupping or Lifting function. Adjust to a comfortable intensity level.
4. **Apply Device:** Firmly place the device against your skin. Keep a slow and continuous movement. Treat each area for up to 3 minutes. The device will automatically shut off after 10 minutes.
5. **Release Suction:** If the suction is too strong or you wish to move the device, press the 'RELEASE' button in the center to quickly release the suction.

10 Minutes to a Better You



Image: Visual guide demonstrating the four simple steps for using the device: apply lotion, select cup, choose mode/intensity, and apply to skin.

5. Recommended Usage

Use for 10 minutes daily for a complete session. Treat each specific area for no more than 3 minutes. Consistent use is recommended for optimal results.

Video: A demonstration of the SUNLUXY BCE1 device being used on various body parts, highlighting its features and ease of use for a healthy lifestyle.

MAINTENANCE

1. Cleaning the Device

After each use, power off the device. Remove the suction cup and filter. Wipe the device head and suction cup with a clean, damp cloth. Do not use harsh detergents or immerse the device in water. Ensure all parts are dry before reassembling or storing.

2. Storage

Store the device in a cool, dry place away from direct sunlight and high humidity. Keep it in its original packaging or a protective case to prevent damage.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on	Low battery; Device not powered on correctly	Charge the device fully; Press and hold the Power button
Weak or no suction	Suction cup not properly attached; Filter clogged; Skin not adequately moisturized	Reattach suction cup securely; Clean or replace filter; Apply more lotion/oil to skin
Skin redness or temporary marks	Normal reaction to cupping/suction; Intensity too high	This is normal and should subside; Reduce intensity level; Do not treat the same area for too long
Device stops unexpectedly	Automatic 10-minute shut-off; Low battery	This is a safety feature; Recharge the device

SPECIFICATIONS

- **Product Name:** 3-IN-1 Sculpting Machine
- **Model Number:** BCE1
- **Red Light Wavelength:** 630+5nm
- **Cupping (with suction cups):** 20-45kPa
- **Battery:** 3.7V 1400mAh
- **Charging Time:** 2-4 hours
- **Full Charge Usage Time:** ≥60 minutes
- **Product Dimensions:** 6.02 x 2.36 x 1.18 inches
- **Item Weight:** 15.2 ounces (0.43 Kilograms)
- **Power Source:** Battery Powered
- **Material:** Plastic
- **Manufacturer:** GUANGDONG ACE-TEC CO.,LTD

WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the product packaging or contact SUNLUXY customer service directly. Keep your purchase receipt for any warranty claims.