

[Manuals.plus](#) /

› [medicube](#) /

› medicube Booster Pro X2 Facial Device User Manual

medicube Booster Pro X2

medicube Booster Pro X2 Facial Device User Manual

1. INTRODUCTION

The medicube Booster Pro X2 is an advanced 7-in-1 facial beauty device designed to enhance various aspects of skin care. It features dual-head technology, 7 distinct modes, and 6 adjustable intensity levels to address concerns such as skin glow, pore appearance, skin definition, and overall radiance. This manual provides detailed instructions for safe and effective use of your device.




Image: The medicube Booster Pro X2 device, showcasing its design and key features like dual-head technology, 7 modes, and 6 intensity levels.

7 Total Skin Solutions

Care for Real Result

DUAL CUSTOMIZING

● Dual Mode

MASK BOOSTING

● Mask Mode

2X GLOW BOOSTING

● Booster Mode

2X BOUNCE BOOSTING

● MC Mode

2X SKIN-DEFINING

● Derma Shot Mode

2X PORE REFINING

● Air Shot Mode

AI PERSONALIZING

● AI Mode*

Image: An illustration depicting the seven distinct skin solutions provided by the Booster Pro X2 device.

2. SAFETY INFORMATION

For external use only. Do not use on broken, irritated, or infected skin. Consult a physician before use if you have any pre-existing skin conditions, are pregnant, or have electronic medical implants. Keep out of reach of children.

3. SETUP

3.1 Charging the Device

Before first use, fully charge the device. Use a USB Type-C cable with an input voltage of DC 5V and a maximum of 2A. The device includes a charging cable.

Button Guide



1. MODE / INTENSITY

Turn on with the Button ①, switch modes with Button ②, and adjust intensity with the Button ①.

2. LANGUAGE

Hold the Buttons ① and ②, select a language with the Button ②, then press the Button ① to save.

Image: The medicube Booster Pro X2 device, its user manual, and the included USB-C charging cable.

3.2 App Connection

Connect the device to the medicube app for remote control, AI-assisted care recommendations, and customized Dual Mode combinations. Scan the QR code in the physical manual or search for the 'AGE-R' app in your device's app store.

Video: A quick guide demonstrating how to use the AGE-R App, including connection and navigation for personalized skincare routines.

4. OPERATING INSTRUCTIONS

4.1 Button Guide

The device features two main buttons for operation:

1. **Power/Intensity Button (1):** Turns the device on/off and adjusts intensity levels.
2. **Mode Button (2):** Switches between the different operating modes.

To change the language, hold both buttons (1 and 2) simultaneously, then use Button 2 to select the language and

NEW

DUAL MODE



Better Multi Care

Two Modes in Less Time

*Clinical Trial Tested by Global Institute of Dermatological Sciences on 20 participants to examine Booster Mode + Derma Shot Mode (2026.01.28–2026.02.11), and on 21 participants each to examine Booster Mode + MC Mode (2026.02.04–2026.02.19) and Derma Shot Mode + MC Mode (2026.02.05–2026.02.19). / Used with medicube PDRN Pink Peptide Serum twice daily for 2 weeks. / Results may vary.

Glow Skin Definition

⬆ Derma Shot X2


⬇ Booster X2



Double Bounce

⬆ MC X2

⬇ Derma Shot X2



Youthful Glow Boosting

⬆ MC X2

⬇ Booster X2



Image: A detailed view of the Booster Pro X2's control buttons and instructions for basic operation and language selection.

4.2 Modes and Usage

The Booster Pro X2 offers 7 distinct modes, each targeting specific skin concerns. It is recommended to use the device on clean skin with a suitable skincare product (serum, cream, or sheet mask) for optimal results, unless otherwise specified.

- **Booster Mode (Glow Boosting):** Enhances skin radiance and hydration. Apply a hydrating serum or cream before use. Glide the device gently over the face in upward motions.

BOOSTER MODE

2X GLOW BOOSTING

Radiance care for a luminous glow

24Hr Lasting Hydration **76.3%** ^{**} ↑

Skin Radiance **28.3%** ^{*} ↑

Skincare Hydration **69.2%** ^{*} ↑

*Clinical Trial Tested by Global Institute of Dermatological Sciences on 21 participants to examine Booster X2 Mode (2026.01.13-2026.01.27) / After 1 use with medicube PDRN Pink Peptide Serum and Booster X2 Mode Level 6 / Results may vary. / Images Provided by Global Ambassadors.
**Clinical Trial Tested by Global Institute of Dermatological Sciences on 21 participants to examine Booster X2 Mode (2026.01.13-2026.01.27) / After 24 hours with medicube PDRN Pink Peptide Serum and Booster X2 Mode Level 6 / Results may vary. / Images Provided by Global Ambassadors.

Before use

✓

After

Image: The device in Booster Mode, emitting an orange glow as it is used for glow boosting.

- MC Mode (Bounce Boosting):** Improves skin texture and elasticity, targeting fine lines and wrinkles. Apply a firming serum or cream. Use upward lifting motions.

Video: A quick guide to the Booster Pro X2, demonstrating its various modes and functions.

MC MODE

2X BOUNCE BOOSTING

Skin texture Care for Healthy-Looking Bounce

Forehead Line **28.3%** ^{*} ↓

Smile Line **69.2%** ^{*} ↓

Neck Line **76.3%** ^{*} ↓

Eye Area Line **16.8%** ^{*} ↓

*Clinical Trial Tested and Images Provided by Global Institute of Dermatological Sciences on 20 participants to examine MC X2 Mode (2026.01.14-2026.01.29) / Used with medicube PDRN Pink Peptide Serum and MC X2 Mode Level 6 twice daily for 2 weeks. / Results may vary. / Images Provided by Global Ambassadors.
**Clinical Trial Tested and Images Provided by Global Institute of Dermatological Sciences on 20 participants to examine MC X2 Mode (2026.01.14-2026.01.29) / After 1 use with medicube PDRN Pink Peptide Serum and MC X2 Mode Level 6 twice daily for 2 weeks. / Results may vary. / Images Provided by Global Ambassadors.

Before use

✓

After

Image: The device in MC Mode, emitting a green glow as it is used for bounce boosting.

- Derma Shot Mode (Skin Definition):** Focuses on enhancing facial contours and definition. Apply a sculpting or firming product. Use upward and outward strokes.

DERMA SHOT MODE

2X SKIN-DEFINING

Skin Definition Care for Enhanced Facial Balance

Under-Eye **21.7%** ^{*} ↑

Jawline **7.73%** ^{*} ↑

Cheek **5.72%** ^{*} ↑

*Clinical Trial Tested by Global Institute of Dermatological Sciences on 22 participants to examine Derma Shot X2 Mode (2026.01.20-2026.02.03) / Used with medicube PDRN Pink Peptide Serum and Derma Shot X2 Mode Level 6 twice daily for 2 weeks. / Results may vary. / Images Provided by Global Ambassadors.

Before use

✓

After

Image: The device in Derma Shot Mode, emitting an orange glow as it is used for skin definition.

- **Air Shot Mode (Pore Refining):** Targets pore appearance and skin texture. This mode is typically used on clean, dry skin without products. Glide over areas with visible pores.

AIR SHOT MODE

2X PORE REFINING
Refining care for pores & skin texture

Blackhead & Whitehead **10.5%*** ↓
Pores **31.4%**** ↓
Dead Skin **80.8%*** ↓
Excess Sebum **28.1%*** ↓

*Clinical Trial Tested and Images Provided by Global Institute of Dermatological Sciences (GIDS) to examine Air Shot X2 Mode (2026.01.21-2026.02.04) / Used with medicube PDRN Pink Peptide Serum and Air Shot X2 Mode Level 5 / Work daily for 4 weeks. / Results may vary. **Clinical Trial Tested and Images Provided by Global Institute of Dermatological Sciences (GIDS) to examine Air Shot X2 Mode (2026.01.21-2026.02.04) / After 1 use with medicube PDRN Pink Peptide Serum / Air Shot X2 Mode Level 6 / Results may vary.

Before use

After

Image: The device in Air Shot Mode, emitting a blue glow as it is used for pore refining.

- **Dual Mode:** Combines two modes simultaneously for a comprehensive and time-saving treatment. Select desired dual modes via the medicube app. Apply appropriate skincare products for the combined modes.

NEW DUAL MODE

NEXT-GEN DUAL CUSTOMIZATION

Combine two care modes at once for time-saving, multi-benefit customized care.

*Dual Mode Settings Available in App and Booster Mode and Derma Shot Mode are provided as default modes. / LED indications can only be selected in a single color mode.

Glow Skin Definition
Derma Shot X2
Booster X2

Double Bounce
MC X2
Derma Shot X2

Youthful Glow Boosting
MC X2
Booster X2

Image: The device operating in Dual Mode, combining two distinct care modes for enhanced treatment.

- **Mask Mode (Sheet Mask Boost):** Designed to optimize the absorption of sheet mask essences. Apply a sheet mask to your face, then gently glide the device over the mask.



Image: The device being used over a sheet mask to enhance the absorption of the mask's essence.

- **AI Mode:** Provides personalized skincare routines based on your skin condition and goals, accessible through the medicube app. Follow the app's recommendations for product application and device usage.




Image: The medicube app displaying the AI Mode interface, offering personalized skincare routine recommendations.

4.3 Layering Modes

Different modes can be incorporated into your skincare routine based on your needs:

- **Without Skincare:** Air Shot Mode (2-3 uses weekly).
- **With Skincare:** Booster, MC, Derma Shot, and Dual Modes (1-3 uses daily).
- **With Sheet Mask:** Mask Mode (1 use daily).



Image: An infographic illustrating the recommended layering of Booster Pro X2 modes based on skincare routine steps.

5. MAINTENANCE

To ensure the longevity and hygiene of your medicube Booster Pro X2, follow these maintenance guidelines:

- After each use, gently wipe the device head with a soft, damp cloth to remove any residual skincare products.
- Do not use abrasive cleaners, alcohol, or harsh chemicals, as these can damage the device.
- Ensure the device is completely dry before storing it.
- Store the device in a cool, dry place, away from direct sunlight and extreme temperatures.

6. TROUBLESHOOTING

If you encounter any issues with your medicube Booster Pro X2, refer to the following common troubleshooting steps:

- **Device not turning on:** Ensure the device is fully charged. If the battery is low, connect it to the charger and allow sufficient time for charging.
- **No sensation or weak sensation during use:** Ensure adequate contact with the skin. Apply a sufficient amount of conductive skincare product (for modes requiring it). Increase the intensity level if comfortable.
- **App connection issues:** Ensure Bluetooth is enabled on your smartphone and the device is within range. Restart both the device and the app.
- **Unusual noise or odor:** Discontinue use immediately and contact customer support.

If the problem persists after attempting these steps, please contact medicube customer support for further assistance.

7. SPECIFICATIONS

Feature	Detail
Model Name	Booster Pro X2
Dimensions	10.08 x 3.98 x 2.83 inches
Weight	1.21 Pounds
Battery	1 Lithium Ion battery (included)

Feature	Detail
Charging Input	DC 5V, max 2A (USB Type-C)
Modes	7 (Booster, MC, Derma Shot, Air Shot, Dual, Mask, AI)
Intensity Levels	6
Manufacturer	APR
Model Number	VAPMEJP26098R

8. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the official medicube website or contact their customer service directly. Keep your proof of purchase for any warranty claims.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

What is the minimum age requirement for the medicube Booster Pro X2, and can an older teenager safely use the EMS (Derma Shot) mode?

The official user manual for the medicube Booster Pro X2 does not specify a strict minimum age for use; however, it does include a safety warning to "Keep out of reach of children."Regarding the...

[Read full answer](#)

