AQUAPICK AQ-300

Aquapick AQ-300 Oral Care Dental Irrigator Water Flosser Instruction Manual

Model: AQ-300

Introduction

Welcome to the Aquapick AQ-300 user manual. This document provides essential information for the safe and effective operation, maintenance, and troubleshooting of your new oral care dental irrigator. The Aquapick AQ-300 is designed to provide thorough cleansing and support oral hygiene, including prevention of peri-implantitis, care for orthodontic patients, prevention of tooth decay, and halitosis treatment.

Important Safety Information

Please read all instructions carefully before using the Aquapick AQ-300. Failure to follow these guidelines may result in injury or damage to the product.

- Do not use with broken or damaged parts.
- Keep away from children.
- Do not use while sleeping or unconscious.
- Unplug from the power source when not in use.
- If the power cord is damaged, stop using the device and replace the cord immediately.
- Do not use if you have any medical conditions that make you prone to infection or if you are taking any medications that affect your ability to heal. Consult your dental professional if you have concerns.
- Always follow the instructions provided with the product.

Package Contents

Verify that all items are present in your Aquapick AQ-300 package:

- Aquapick AQ-300 Main Unit
- Water Reservoir
- Standard Jet Tips (multiple, typically 4)

- Power Cord
- Instruction Manual (this document)



Image: Four standard jet tips for the Aquapick AQ-300 water flosser. Each tip is clear with a distinct colored band (green, red, blue, orange) at its base for easy identification among multiple users. These tips are essential for directing the water jet during oral irrigation.

Product Overview

The Aquapick AQ-300 is a powerful oral irrigator designed for effective plaque removal and gum health improvement. It features a high-pressure pulsating water jet to clean between teeth and below the gumline, where traditional brushing and flossing may not reach. The unit includes a water reservoir, a handle with a tip ejection button, and a pressure control system.

Setup

Follow these steps for initial setup of your Aquapick AQ-300:

1. **Unpack:** Remove all components from the packaging.
2. **Place Unit:** Position the main unit on a flat, stable surface near an electrical outlet and a sink.
3. **Fill Reservoir:** Remove the water reservoir from the base and fill it with lukewarm water. You may also add a small amount of mouthwash if desired. Place the reservoir back onto the base, ensuring it is seated correctly.

4. **Attach Tip:** Select a standard jet tip. Insert the tip into the handle until it clicks into place. To remove, press the tip ejection button on the handle.
5. **Plug In:** Connect the power cord to the unit and then plug it into a suitable electrical outlet (AC 90~240V, 50~60Hz).

Operating Instructions

For optimal oral hygiene, use your Aquapick AQ-300 daily:

1. **First Use (Priming):** If using for the first time or after extended storage, direct the tip into the sink. Turn the unit ON and set the pressure to high until water flows consistently.
2. **Adjust Pressure:** The Aquapick AQ-300 offers adjustable pressure settings. Start with the lowest setting and gradually increase to a comfortable level. The maximum pressure is 7 kgf/cm².
3. **Position Tip:** Lean over the sink and place the tip in your mouth, aiming the water jet at your gumline at a 90-degree angle. Close your lips slightly to prevent splashing, but allow water to flow into the sink.
4. **Begin Flossing:** Turn the unit ON. Guide the tip along your gumline, pausing briefly between teeth. Ensure you clean both the front and back surfaces of your teeth. The unit provides 1,800 pulsations per minute for effective cleaning.
5. **Duration:** Continue until the reservoir is empty, which typically takes about 1.5 seconds for 600ml of water.
6. **Power Off:** Turn the unit OFF after use.
7. **Empty Reservoir:** Empty any remaining water from the reservoir to prevent bacterial growth.

Maintenance

Regular cleaning and maintenance will ensure the longevity and optimal performance of your Aquapick AQ-300.

Cleaning the Unit

- Always unplug the unit before cleaning.
- Wipe the exterior of the main unit with a soft, damp cloth. Do not immerse the unit in water.
- Avoid using abrasive cleaners or harsh chemicals.

Cleaning the Reservoir

- Remove the reservoir from the base.
- Wash the reservoir with warm, soapy water. Rinse thoroughly.
- For deeper cleaning, you can run a solution of white vinegar (1 part) and water (2 parts) through the unit periodically to remove mineral deposits. Follow with a full reservoir of clean water to rinse.

Cleaning the Tips

- Remove the tip from the handle.
- Rinse the tip under running water.
- For thorough cleaning, soak tips in a solution of hydrogen peroxide or denture cleaner for 5-10 minutes. Rinse well before use.
- Replace tips every 3-6 months, or sooner if they become worn or damaged.

Storage

- Ensure the unit is unplugged and empty of water before storing.
- Store the unit and tips in a clean, dry place.

Troubleshooting

If you encounter issues with your Aquapick AQ-300, refer to the following table for common problems and solutions.

Problem	Possible Cause	Solution
Unit does not turn on.	Not plugged in; Power outlet issue; Unit malfunction.	Ensure power cord is securely plugged into a working outlet. Test the outlet with another device. If still not working, contact customer support.
Water does not flow.	Reservoir empty; Tip clogged; Unit not primed; Air lock.	Fill reservoir. Remove and clean tip. Prime the unit by running it on high with the tip in the sink until water flows. Ensure reservoir is seated correctly.
Low water pressure.	Pressure setting too low; Mineral buildup; Clogged tip.	Increase pressure setting. Clean the unit and tips for mineral buildup (refer to Maintenance section). Replace tip if necessary.
Leaking from unit.	Reservoir not seated correctly; Damaged hose/seal.	Ensure reservoir is firmly in place. Inspect for visible damage. If leaking persists, contact customer support.

Specifications

Technical specifications for the Aquapick AQ-300 Oral Care Dental Irrigator:

Feature	Specification
Power	AC 90~240V, 50~60Hz
Capacity	24W
Pressure	Up to 7 kgf/cm²
Pulsation Rate	1,800 / minute
Working Time (600ml)	Approximately 1.5 seconds
Dimensions (H x L x W)	188 x 150 x 90 mm (7.4 x 5.91 x 3.54 inches)
Reservoir Volume	600ml
Weight	900g (1.98 Pounds)
Power Source	Corded Electric
Product Benefits	Cleansing

Warranty and Support

For warranty information, technical support, or service inquiries regarding your Aquapick AQ-300, please refer to the warranty card included with your product or contact Aquapick customer service directly. Keep your purchase receipt as proof of purchase.