Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

ACCU-ANSWER TC01

ACCU-ANSWER Cholesterol Test Strips TC01 User Manual

Model: TC01 | Brand: ACCU-ANSWER

1. Introduction

This manual provides essential instructions for the safe and effective use of ACCU-ANSWER Cholesterol Test Strips TC01. These test strips are designed for use with the ACCU-ANSWER 4-in-1 Test Meter to measure total cholesterol levels. Please read this manual thoroughly before using the test strips to ensure accurate results and proper operation.

2. Important Safety Information

- For *in vitro* diagnostic use only. Do not ingest.
- Keep out of reach of children.
- Store test strips in their original vial with the cap tightly closed.
- Do not refrigerate.
- Use within 3 months after first opening the strip vial.
- Do not use expired test strips. Check the expiration date on the vial.
- Dispose of used test strips and lancets according to local regulations for biohazardous waste.
- Consult a healthcare professional for interpretation of test results.

3. Package Contents

Verify that all items are present in your package:

- 25 x ACCU-ANSWER Cholesterol Test Strips (in a vial)
- 25 x Lancets



Figure 1: Package contents showing one vial of 25 test strips and 25 individual lancets.

4. Setup and Preparation

1. **Gather Supplies:** Ensure you have your ACCU-ANSWER 4-in-1 Test Meter, a Cholesterol Test Strip, a lancet, and an alcohol swab (not included).
2. **Wash Hands:** Wash your hands thoroughly with soap and warm water. Dry them completely.
3. **Prepare Lancet Device:** Insert a new lancet into your lancing device according to the device's instructions.
4. **Prepare Test Strip:** Remove one test strip from the vial. Immediately close the vial cap tightly to protect the remaining strips from moisture.
5. **Insert Test Strip:** With the meter turned off, insert the test strip into the test strip port of your ACCU-ANSWER 4-in-1 Test Meter. The meter should turn on automatically and display a test strip symbol.



Figure 2: ACCU-ANSWER Cholesterol Test Strips vial and packaging. Note the tight-fitting cap for strip preservation.

5. Operating Instructions: Performing a Cholesterol Test

1. **Obtain Blood Sample:** Use the lancing device to prick the side of your fingertip. Gently squeeze your finger to form a round drop of blood.
2. **Apply Blood Sample:** Touch the tip of the test strip to the blood drop. The blood will be drawn into the strip automatically. Do not apply blood to the top of the strip. Ensure the blood sample completely fills the reaction area.
3. **Wait for Results:** The meter will begin a countdown. Do not remove the test strip during this time.
4. **Read Results:** After the countdown, your cholesterol level will be displayed on the meter's screen. Record your result if necessary.
5. **Dispose of Used Materials:** Carefully remove the used test strip and lancet. Dispose of them in a puncture-resistant container (sharps container) according to local regulations.

For detailed instructions on using the ACCU-ANSWER 4-in-1 Test Meter, please refer to its specific user manual.

6. Maintenance and Storage

6.1 Storage of Test Strips

- Store the test strip vial in a cool, dry place at room temperature (typically 4°C - 30°C or 39°F - 86°F).
- Do not store in direct sunlight or high humidity.
- Always keep the vial cap tightly closed immediately after removing a test strip.
- Do not transfer test strips to another container.
- Use all test strips within 3 months of opening the vial, or before the expiration date printed on the vial, whichever comes first.

6.2 Cleaning

The test strips themselves do not require cleaning. Ensure your ACCU-ANSWER 4-in-1 Test Meter is kept clean according to its manual.

7. Troubleshooting

If you encounter issues, refer to the troubleshooting section of your ACCU-ANSWER 4-in-1 Test Meter manual. Common issues related to test strips include:

- **Error Message:** Ensure the test strip is inserted correctly and that it is not expired.
- **Insufficient Blood Sample:** Apply a larger, well-formed blood drop to the strip. Do not add more blood after the initial application.
- **Inaccurate Results:** Verify that the test strips are stored correctly and are not expired. Ensure proper technique for blood application. If issues persist, contact customer support.

8. Specifications

Feature	Detail
Product Name	Cholesterol Test Strips
Model	TC01
Brand	ACCU-ANSWER
Quantity	25 Test Strips, 25 Lancets
Storage Temperature	4°C - 30°C (39°F - 86°F)
Use After Opening	Within 3 months
UPC	705794433354

9. Warranty and Customer Support

9.1 Warranty Information

ACCU-ANSWER Cholesterol Test Strips are consumable items. Please check the packaging for any specific warranty information regarding manufacturing defects. Generally, test strips are guaranteed to be free from defects until their expiration date, provided they are stored and used according to the instructions.

9.2 Customer Support

For technical assistance, questions about your ACCU-ANSWER Cholesterol Test Strips, or to report any issues, please contact your retailer or the manufacturer directly. Refer to the product packaging or the ACCU-ANSWER website for contact details.

© 2024 ACCU-ANSWER. All rights reserved.
This manual is for informational purposes only. Always consult a healthcare professional for medical advice.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0CNYJ9JCB

More documents and user manuals: **https://manuals.plus/**