AURAGLOW B0CRGSK7Y

Auraglow Teeth Whitening Kit: Wired LED Accelerator Light User Manual

Introduction

This manual provides instructions for the safe and effective use of your Auraglow Teeth Whitening Kit. This kit is designed for at-home teeth whitening, utilizing a wired LED accelerator light and 35% carbamide peroxide gel.

Key Features:

- LED accelerator light for enhanced whitening.
- Plugs into various devices (USB, Lightning, Type-C, Micro USB).
- Professional-level whitening results.
- Fast-acting formula, whitens in 30 minutes.
- Includes 12 whitening treatments.
- Utilizes dental laser whitening technology.
- Complete kit with gel, light, and cover case.

What's in the Box

The Auraglow Teeth Whitening Kit includes the following components:

- 1 x Wired LED Teeth Whitening Light (with USB, Lightning, Type-C, and Micro USB adapters)
- 2 x Syringes of 35% Carbamide Peroxide Whitening Gel
- 1 x LED Light Cover Case



Image: The kit components, including the LED light, gel syringes, and storage case.

How to Use

Follow these steps for optimal teeth whitening results:

1. **Prepare Your Teeth:** Brush and floss your teeth thoroughly before each treatment. Ensure your teeth are clean and dry.
2. **Apply Whitening Gel:** Open the cap of one whitening gel syringe. Carefully apply a thin, even layer of the 35% carbamide peroxide gel along the top and bottom inner surfaces of the LED mouthpiece tray. Avoid overfilling to prevent gel from contacting gums.
3. **Connect and Insert LED Light:** Plug the wired LED light into your smartphone (Lightning, Type-C, Micro USB) or any standard USB power source. The LED light will illuminate. Carefully insert the mouthpiece into your mouth, ensuring your teeth are fully covered by the gel-coated tray.
4. **Begin Treatment:** Keep the mouthpiece in your mouth for 30 minutes. The LED light will remain active during this period.
5. **Complete Treatment:** After 30 minutes, remove the mouthpiece. Rinse your mouth thoroughly with water to

remove any remaining gel. Clean the LED mouthpiece with water and store it in its cover case.



Image: Visual guide illustrating the four steps of using the whitening kit: filling the tray, plugging in the light, wearing for 30 minutes, and enjoying results.

Care and Maintenance

To ensure the longevity and effectiveness of your Auraglow Teeth Whitening Kit:

- After each use, rinse the LED mouthpiece thoroughly with cool water. Do not use hot water as it may damage the tray.
- Store the mouthpiece in its protective cover case in a cool, dry place.
- Keep the whitening gel syringes sealed and stored in a cool, dark place, away from direct sunlight.
- Avoid submerging the entire wired LED light unit in water. Only clean the mouthpiece portion.

Troubleshooting

LED Light Not Turning On:

- Ensure the LED light is securely plugged into a compatible power source (smartphone or USB port).
- If using with a smartphone, check your phone's settings to ensure external accessories are permitted (e.g., for iPhone, check "Settings > Face ID & Passcode > USB Accessories" or "Settings > General > Transfer or Reset iPhone > Reset > Reset Location & Privacy").
- Try plugging the device into a different USB port or smartphone to rule out a faulty power source.

Uneven Whitening Results:

- Ensure the whitening gel is applied evenly across all tooth surfaces that come into contact with the mouthpiece.
- Make sure the mouthpiece is properly seated in your mouth, covering all visible teeth.
- Consistent daily use for the recommended duration is crucial for best results.

Tooth Sensitivity:

- If you experience sensitivity, reduce the treatment time or frequency.
- Ensure gel is not applied directly to gums. If gel contacts gums, wipe it off immediately.
- Consider using a sensitive toothpaste during your whitening regimen.

Specifications

Feature	Detail
Brand	AURAGLOW
Active Ingredients	35% Carbamide Peroxide
Product Benefits	Teeth Whitening
Item Form	Gel
Number of Items	1 Kit (includes 2 gel syringes)
Power Source	Wired (USB, Lightning, Type-C, Micro USB compatible)
Treatment Duration	30 minutes per session
Origin	Made in USA (gel)

Important Safety Information

- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.
- Consult your dentist before use if you have sensitive teeth, gum disease, or any existing dental work (crowns, veneers, fillings).
- Do not use if pregnant or breastfeeding.
- Keep out of reach of children.
- Avoid contact of the whitening gel with eyes or skin. Rinse immediately with water if contact occurs.
- Do not swallow the whitening gel.

- Discontinue use if irritation or excessive sensitivity occurs.

Warranty and Support

Auraglow offers a 60-day satisfaction guarantee for this product. For warranty claims, product support, or further assistance, please contact Auraglow customer service through their official website or the retailer where the product was purchased. Please retain your proof of purchase.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question