HAUBAT Laser Cap

HAUBAT Laser Cap: FDA Cleared Laser Hair Growth Device Instruction Manual

Model: Laser Cap

1. Product Overview

The HAUBAT Laser Cap is a hand-held, cap-shaped low-level laser therapy (LLLT) device designed to promote hair growth in both women and men. It emits laser light directly to the scalp, ensuring comprehensive coverage for effective treatment.

This device utilizes advanced laser and red light therapy technology at an optimal wavelength to stimulate hair follicles across the hairline, temples, sides, and lower crown. The medical-grade padding ensures a comfortable fit during use.



Image 1.1: The HAUBAT Laser Cap, illustrating the red light emission during operation.

2. Safety Information

Please read all safety instructions carefully before using the HAUBAT Laser Cap. Failure to follow these instructions may result in injury or device malfunction.

- **Eye Safety:** Do not look directly into the laser light. The device is designed for scalp treatment only.
- **Medical Conditions:** Consult a healthcare professional before use if you have any pre-existing medical conditions, are pregnant, or are undergoing other medical treatments.
- **Children:** Keep out of reach of children. This device is not intended for use by individuals under 18 years of age.
- **Damage:** Do not use the device if it appears damaged. Contact customer support for assistance.
- **Cleaning:** Refer to the Maintenance section for proper cleaning instructions. Do not immerse the device in water.
- **Power Source:** Use only the provided power adapter or a compatible power bank. Ensure the power source meets the device's specifications.

3. Indications for Use

The HAUBAT Laser Cap is indicated to promote hair growth in individuals with androgenic alopecia (pattern hair loss) who meet specific classifications:

- **Males:** Norwood-Hamilton classifications of IIa-V.
- **Females:** Ludwig-Savin Classifications of I-II.
- **Skin Phototypes:** Applicable for individuals with Skin Phototypes I-IV.



Image 3.1: Visual guide for applicable users based on hair loss classification and skin phototype.

4. Package Contents

Upon opening the package, please verify that all components are present and undamaged:

- HAUBAT Laser Cap

- Wall Power Adapter
- User Manual (this document)

Note: A portable power bank is not included but can be used as an alternative power source.

5. Setup

1. **Unpack:** Carefully remove the HAUBAT Laser Cap and all accessories from the packaging.
2. **Inspect:** Check the device for any visible damage. Do not use if damaged.
3. **Connect Power:**
    - For wall power: Connect the provided wall adapter to the cap's power port, then plug the adapter into a standard electrical outlet.
    - For portable power: Connect a compatible portable power bank (not included) to the cap's power port using a suitable USB cable. Direct plug-in is recommended for optimal performance.
4. **Position:** Place the cap comfortably on your head, ensuring it covers the areas of the scalp you wish to treat. The medical-grade padding should gently hug your head.




Image 5.1: Proper placement of the HAUBAT Laser Cap during use, demonstrating comfort and ease of wear.

6. Operating Instructions

1. **Power On:** Once connected to a power source, press the power button located on the cap to activate the laser lights. The lights will illuminate, indicating the device is operational.
2. **Treatment Session:** Wear the cap for the recommended treatment duration of 30 minutes per session. The device is designed to provide consistent laser light during this period.
3. **Power Off:** The device may automatically turn off after the session. If not, press the power button again to turn off the laser lights.
4. **Disconnect:** Disconnect the cap from its power source after each use.

7. Treatment Schedule

For optimal results, adhere to the following treatment schedule:

- **Frequency:** Use the HAUBAT Laser Cap 3 times per week, spaced out (e.g., every other day).
- **Duration:** Each treatment session should last 30 minutes.
- **Recommended Course:** A typical course of treatment is 16 weeks.
- **Expected Results:** Visible improvement in hair thickness and growth can typically be observed within 3 to 6 months of consistent use. Individual results may vary.




Image 7.1: Visual representation of the recommended treatment schedule and laser wavelength (650nm).




Image 7.2: Examples of hair growth results observed after consistent use of the laser cap over several months.

8. Maintenance

Proper care and maintenance will extend the life of your HAUBAT Laser Cap.

- **Cleaning:** Ensure the device is powered off and disconnected from the power source before cleaning. Gently wipe the interior and exterior surfaces of the cap with a soft, dry cloth. Do not use abrasive cleaners, solvents, or immerse the device in water.
- **Storage:** Store the cap in a cool, dry place away from direct sunlight and extreme temperatures. Keep it in its original packaging or a protective case when not in use to prevent dust accumulation and damage.

9. Troubleshooting

If you encounter any issues with your HAUBAT Laser Cap, please refer to the following common troubleshooting steps:

Problem	Possible Cause	Solution
Device does not turn on.	Not properly connected to power source. Power source issue. Device malfunction.	Ensure the power adapter or power bank is securely connected to the cap and the power outlet/power bank is functional. Try a different outlet or power bank. If the issue persists, contact customer support.
Lights are dim or flickering.	Insufficient power. Loose connection.	Verify the power connection is secure. If using a power bank, ensure it is adequately charged. Use the wall adapter for the most stable performance.
Cap feels uncomfortable.	Improper positioning.	Adjust the cap's position on your head to ensure a comfortable fit. The internal padding is designed to conform gently.

If you experience any issues not listed here or if the problem persists after attempting these solutions, please contact HAUBAT customer support.

10. Specifications

- **Model:** Laser Cap
- **Manufacturer:** HAUBAT
- **ASIN:** B0F22FZR5V
- **UPC/GTIN:** 792522470477
- **Wavelength:** 650nm (Red Light)
- **Treatment Duration:** 30 minutes per session
- **Power Source:** Wall adapter or compatible portable power bank
- **Package Dimensions:** 8.58 x 8.43 x 3.82 inches
- **Weight:** 1.08 Pounds



Image 10.1: Key features of the HAUBAT Laser Cap, including its lightweight design, universal fit, and FDA clearance.

11. Warranty and Support

HAUBAT is committed to providing high-quality products. For information regarding warranty coverage, product registration, or technical support, please refer to the warranty card included with your purchase or visit the official HAUBAT website.

If you require assistance, please contact our customer support team with your product model and purchase details ready. Contact information can typically be found on the manufacturer's website or on the product packaging.

© 2025 HAUBAT. All rights reserved.
This manual is for informational purposes only. Specifications are subject to change without notice.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question