



## aidapt VM924 One Way Drinking Straw Instructions

[Home](#) » [Aidapt](#) » aidapt VM924 One Way Drinking Straw Instructions 

### Contents

- [1 aidapt VM924 One Way Drinking Straw](#)
- [2 INSTRUCTIONS FOR USE](#)
- [3 CARE & MAINTENANCE](#)
- [4 Documents / Resources](#)
  - [4.1 References](#)



**aidapt VM924 One Way Drinking Straw**



## **INTRODUCTION**

Thank you for deciding to purchase the One-way Drinking Straw from Aidapt. The One-way Drinking Straw is manufactured from the finest materials available. When used correctly it is designed to give many years of reliable, trouble-free service.

## **BEFORE USE**

Carefully remove all packaging. Avoid using any knives or other sharp instruments as this may damage the surface of the product.

Check the product for any visible damage. If you see any damage or suspect a fault, do not use your product, but contact your supplier for support.

## **INTENDED USE**

The One-Way Drinking Straw makes drinking easier for people with weak oral musculature or those who have difficulty swallowing. The straw is a dining aid featuring a one-way valve which stays filled with liquid even after removing the straw from the lips, this reduces the air ingested and the effort required by the user, leading to less coughing and choking while drinking thin liquids.

NB. The cup (pictured above) is not included.

## **INSTRUCTIONS FOR USE**

The One-Way Drinking Straw is simple to use:

- Slip the drinking straw into the cup or glass, with the adjustable clip on the outside edge of the cup.
- Use the drinking straw as normal.

NB. Not recommended for use with hot drinks.

## CLEANING

Clean your One-way Drinking Straw using a non-abrasive cleaner or mild detergent with a soft cloth. Abrasive cleaners e.g. AJAX and/ or abrasive cleaning pads may damage the product beyond repair and should not be used. Always ensure that you wipe dry the equipment after cleaning.

## REISSUE

If you reissue or are about to reissue this product, please thoroughly check all components for their safety. If in any doubt, please do not issue or use, but immediately contact your supplier for service support.

## CARE & MAINTENANCE

Please make a safety check of the product at regular intervals or if you have any concerns.

## IMPORTANT INFORMATION

The information given in this instruction booklet must not be taken as forming part of or establishing any contractual or other commitment by Aidapt Bathrooms Limited, Aidapt (Wales) Ltd or its agents or its subsidiaries and no warranty or representation concerning the information is given.


Please exercise common sense and do not take any unnecessary risks when using this product; as the user you must accept liability for safety when using the product.

Please do not hesitate to contact the person who issued this product to you or the manufacturer (detailed below) if you have any queries concerning the assembly/use of your product.

Aidapt Bathrooms Ltd, Lancots Lane, Sutton Oak, St Helens, WA9 3EX

- Telephone: +44 (0) 1744 745 020
- Fax: +44 (0) 1744 745 001
- Web: [www.aidapt.co.uk](http://www.aidapt.co.uk)
- Email: [accounts@aidapt.co.uk](mailto:accounts@aidapt.co.uk)
- [adaptations@aidapt.co.uk](mailto:adaptations@aidapt.co.uk)
- [sales@aidapt.co.uk](mailto:sales@aidapt.co.uk)

## Documents / Resources

 <p>aidapt One-Way Drinking Straw VM924 Usage and Maintenance Instructions</p>	<p><a href="#">aidapt VM924 One Way Drinking Straw</a> [pdf] Instructions VM924, VM924 One Way Drinking Straw, One Way Drinking Straw, Drinking Straw, Straw</p>
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------

## References

-  [Aidapt Mobility Aids & Independent Living Solutions - Buy Mobility Aids Online](#)

- [M<sub>H</sub> Search - Manual-Hub.com](#)

Manuals+.