



## Apps CuddleCot Cooling Cot User Guide

[Home](#) » [APPs](#) » Apps CuddleCot Cooling Cot User Guide 

### Contents

- [1 Apps CuddleCot Cooling Cot](#)
- [2 Product Usage Instructions](#)
- [3 Documents / Resources](#)
  - [3.1 References](#)
- [4 Related Posts](#)



**Apps CuddleCot Cooling Cot**

**Product Information:**

CuddleCot

The CuddleCot is a device designed to keep a deceased baby cool, allowing families to spend more time with their baby before saying goodbye. It consists of a cooling unit, a cooling mat, and various accessories.

**Key Features:**


- Insulating foil to be placed on the surface of a cot or Moses basket
- Cooling pad for maintaining a cool temperature
- Double connectors for attaching the cooling mat to the cooling unit
- Single connector for connecting the cooling unit to the hose
- Filler cap for adding biocide and water
- Blankets for covering and insulating the baby

**Important Note:**

Please read the complete user manual contained in the CuddleCot box or visit our website at [www.CuddleCot.com/manuals](http://www.CuddleCot.com/manuals) for further instructions and troubleshooting.

**Product Usage Instructions****CuddleCot Quickstart Guide**

1. Place the insulating foil on the surface of the cot or Moses basket, ensuring that the silver side is facing up.
2. Put the cooling pad on top of the foil and cover it with a thin sheet or insert it into a pillowcase.



3. Connect the double connectors of the hose to the cooling mat and the single connector to the cooling unit.




4. Remove the filler cap from the cooling unit, add two drops of biocide, and fill it three-quarters full with purified water.



5. Plug in the cooling unit and switch it on. Check the water level and top up if necessary.
6. Cover the baby with blankets and tuck them in to provide insulation.

For further instructions and troubleshooting see instructions contained in the CuddleCot box or visit our website [www.CuddleCot.com/manuals](http://www.CuddleCot.com/manuals)

## Documents / Resources

 <p><b>CuddleCot Quickstart Guide</b></p> <ol style="list-style-type: none"> <li>1. Place cooling mat on surface of cot/bassinet before baby goes in.</li> <li>2. Place cooling mat on top of cot/bassinet, place in room and plug in.</li> <li>3. Connect the double connectors of the hose to the cooling mat and the single connector to the cooling unit.</li> <li>4. Remove the cap, add two drops of biocide and fill with purified water.</li> <li>5. Plug in the cooling unit and switch it on.</li> <li>6. Check the water level and top up if necessary.</li> <li>7. Cover the baby with blankets and tuck them in to provide insulation.</li> </ol> <p>For further instructions and troubleshooting see instructions contained in the CuddleCot box or visit our website <a href="http://www.CuddleCot.com/manuals">www.CuddleCot.com/manuals</a></p>	<p><b><a href="#">Apps CuddleCot Cooling Cot [pdf] User Guide</a></b>  CuddleCot, CuddleCot Cooling Cot, Cooling Cot</p>
---	--

## References

- [Set-Up Manuals - CuddleCot](#)