



ARROW AK-17702-CDC Two-Lumen CVC Instruction Manual

[Home](#) » [Arrow](#) » ARROW AK-17702-CDC Two-Lumen CVC Instruction Manual 



AK-17702-CDC Two-Lumen CVC Instruction Manual



Contents

- [1 Contents:](#)
- [2 Documents / Resources](#)
- [2.1 References](#)
- [3 Related Posts](#)

Contents:

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
- 1: Spring-Wire Guide: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm

- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint
- Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® II Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- Safety Scalpel: #11
- Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle
- ¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65



WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).


Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (14 Ga) Proximal (18 Ga.)	0.60 0.43	5337 1453

* Priming volumes are approximate and are done without accessories.

+Flow rate values are approximate and are determined using deionized water at 100 cm head height.



Arrow International LLC
A subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA



www.teleflex.com

LBLO53650 ROO (2021-11)

7 Fr. **2** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

Documents / Resources

	ARROW AK-17702-CDC Two-Lumen CVC [pdf] Instruction Manual AK-17702-CDC Two-Lumen CVC, AK-17702-CDC, Two-Lumen CVC
	ARROW AK-17702-CDC Two-Lumen CVC [pdf] Instruction Manual AK-17702-CDC, Two-Lumen CVC, Two-Lumen, Lumen

References

-  p65warnings.ca - [p65warnings Resources and Information.](http://p65warnings.ca)

Manuals+.