

[Skip to content](#)

Manuals+

User Manuals Simplified.



ARROW AK-22123-CDC Lumen CVC for High Volume Infusions User Manual

July 22, 2022 July 23, 2022 [Leave a comment on ARROW AK-22123-CDC Lumen CVC for High Volume Infusions User Manual](#)

[Home](#) » [Arrow](#) » **ARROW AK-22123-CDC Lumen CVC for High Volume Infusions User Manual** 

Contents [hide](#)

[1 ARROW AK-22123-CDC Lumen CVC for High Volume Infusions](#)

[2 Contents](#)

[3 Rx only](#)

[4 Documents / Resources](#)

[4.1 References](#)

[4.2 Related Manuals / Resources](#)

ARROW

ARROW AK-22123-CDC Lumen CVC for High Volume Infusions



Contents

- Arrow+ard Blue® Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 16 cm
- Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- Injection Needle: 25 Ga. x 1" (2.54 cm)
- Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- Pressure Transduction Probe
- Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- Syringe: 3 mL Luer-Lock
- Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®1 One-Step Solution with Hi-Lite Orange™ Tint
- SharpsAway® II Locking Disposal Cup
- SharpsAway® Disposal Cup
- Drape: 41" x 55" with 4" fenestration
- Needle Holder
- Safety Scalpel: #11
- Gauze Pad: 2" x 2" (5 cm x 5 cm)
- Gauze Pad: 4" x 4" (10 cm x 10 cm)
- Suture: 3-0 Silk with Curved Needle

1A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex. Store between 20 – 25°C (68 – 77°F). Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



Lumen

Priming
Volume*
(mL)

Gravity
Flow
Rate†
(mL/hr)

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.35	2581
Medial (12 Ga.)	1.2	9719
Proximal (12 Ga.)	1.1	12666

Priming volumes are approximate and are done without accessories. † Flow rate values are approximate and are determined using deionized water at 100 cm head height. Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

Documents / Resources



[ARROW AK-22123-CDC Lumen CVC for High Volume Infusions](#) [pdf] User Manual
AK-22123-CDC Lumen CVC for High Volume Infusions, Lumen CVC for High Volume Infusions,
High Volume Infusions, Volume Infusions



References

- [Harm.P65Warnings.ca](http://www.Harm.P65Warnings.ca)

Related Manuals / Resources

[ARROW CK 04711 One Lumen CVC User Guide](#)



CK 04711 One Lumen CVC User Guide Contents: 1: One-Lumen Catheter:14 Ga.(2.2mmOD)x16cm,Side Holes 1: Spring-Wire Guide:.032" (0.81mm) dia.x17-3/4"(45cm) (Straight...

[ARROW AK-14402 Small Vein Two-Lumen CVC User Manual](#)



AK-14402 Small Vein Two-Lumen CVC ARROW AK-14402 Small Vein Two-Lumen CVC Contents 1: Two-Lumen Catheter: 4 Fr....

[AK-25123-F Arrow-Howes Large-Bore Multi-Lumen CVC Kit User Guide](#)

Arrowg+ard Blue® Large-Bore Three-Lumen CVC for High Volume Infusions Contents: 1: 1: 1: 1: 1: 1: 1: 1: 1:...



[BRIGHTK XHP70-S High Lumen Flashlights User Guide](#)

XHP70-S Product Features Product Name High Lumen XHP70-S Flashlights LED Light Source CREE XHP70.2 Body Material Aluminum Alloy...

Leave a comment

Your email address will not be published.

Comment

Name

Email

Website

Save my name, email, and website in this browser for the next time I comment.

[Manuals+](#),

- [home](#)
- [privacy](#)