



AK-25123-F Arrow-Howes Large-Bore Multi-Lumen CVC Kit User Guide

[Home](#) » [Arrow](#) » AK-25123-F Arrow-Howes Large-Bore Multi-Lumen CVC Kit User Guide 



Contents

- [1 Arrowg+ard Blue® Large-Bore Three-Lumen CVC for High Volume Infusions](#)
- [2 Documents / Resources](#)
 - [2.1 References](#)
- [3 Related Posts](#)

Arrowg+ard Blue® Large-Bore Three-Lumen CVC for High Volume Infusions

Contents:

1:	Arrowg+ard Blue® Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 20 cm
1:	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J"
1:	Tip on Other) with Arrow Advancer
1:	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1:	Injection Needle: 25 Ga. x 1" (2.54 cm)
1:	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1:	Pressure Transduction Probe
1:	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire
1:	Introduction Syringe
1:	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1:	Syringe: 3 mL Luer-Lock
1:	Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
1:	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1:	3 mL Applicator 2% CHG and 70% IPA Chloraprepel One-Step Solution with Hi-Lite Orange™ Tint Sh
1:	arpsAway® Disposal Cup
1:	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
2:	Scalpel: #11
5:	Gauze Pad: 2" x 2" (5 cm x 5 cm) Gauze Pad: 4" x 4" (10 cm x 10 cm)
1:	Suture: 3-0 Braided Silk with Straight Needle
1A	registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. California Prop. 65



WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F). Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic. Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs.

 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 G a.)	0.40	2457
Medial (12 G a.)	1.3	10151
Proximal (12 G a.)	1.2	11039

* Priming volumes are approximate and are done without accessories. † Flow rate values are approximate and are determined using deionized water at 100 cm head height.





Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL052097 R01 (2022-02)



AK-25123-F

12 Fr. 3Lumen 20cm catheter length .035inch dia. spring-wire guide

ARROW®

Documents / Resources

	ARROW AK-25123-F Arrow-Howes Large-Bore Multi-Lumen CVC Kit [pdf] User Guide Arrow-Howes, Large-Bore, Multi-Lumen, CVC, AK-25123-F
--	---

References

- p65warnings.ca - [p65warnings Resources and Information.](#)

[Manuals+.](#)