



ARROW CDA-34052-HPK1A Pressure Injectable Two-Lumen PICC User Guide

[Home](#) » [Arrow](#) » ARROW CDA-34052-HPK1A Pressure Injectable Two-Lumen PICC User Guide 

Contents

- [1 ARROW CDA-34052-HPK1A Pressure Injectable Two-Lumen PICC](#)
- [2 SPECIFICATION](#)
- [3 Pressure Injectable](#)
- [4 Documents / Resources](#)
- [5 Related Posts](#)

ARROW[®]

ARROW CDA-34052-HPK1A Pressure Injectable Two-Lumen PICC

1: **Surgical Apparel:** Bouffant Cap
 1: **Surgical Apparel:** Mask
 1: **Tape:** Steri-Strip®2
 1: HemoHopper® Fluid Receptacle
 1A trademark of Becton, Dickinson and Company.
 2A registered trademark of 3M Company.
 All components are CE 2797 unless otherwise noted.

Not made with natural rubber latex.

Contraindications: The Pressure Injectable PICC is contraindicated wherever there is presence of device-related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.

5Fr.
 2Lumen
 40cm
 catheter
 length
 .018 inch dia.
 spring-wire
 guide

Pressure Injectable



Lumen

Priming
Volume*
(mL)

Gravity
Flow
Rate†
(mL/hr)

Pump
Flow Rate††
(mL/hr)

MAX Pressure
Injection
Flow Rate**
(mL/sec)

Distal (18 Ga.)
 Proximal (18 Ga.)

0.44
 0.44

568
 568

4720
 4680

5
 5

Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.

Documents / Resources

	ARROW CDA-34052-HPK1A Pressure Injectable Two-Lumen PICC [pdf] User Guide CDA-34052, HPK1A, Pressure Injectable Two-Lumen, Injectable Two-Lumen PICC, Injectable
--	--