

ARROW CS-42703-E Three Lumen CVC Instructions

PB-10005-000 Dimensions: Width: 6-3/4" (17,2 cm) Length: 11-5/16" (28,7 cm)

T04200100A6 Dimensions: 6 5/8" x 10 1/2"

Product of XXX
Packaged in XXX

ARROW®



REF **CS-42703-E** (IPNXXXXXX)

7 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Contents

- [1 Arrowg+ard Blue Plus® Three-Lumen CVC](#)
- [2 Documents / Resources](#)
 - [2.1 References](#)
- [3 Related Posts](#)

Arrowg+ard Blue Plus® Three-Lumen CVC

Contents:

- 1: Arrowg+ard Blue Plus®1 Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹ Licensed under US Patent No. 7,329,412.

All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Not made with natural rubber latex.


Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



CAS # 7440-48-4

 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.38	2490
Medial (18 Ga.)	0.37	1162
Proximal (18 Ga.)	0.39	1328

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



LotNumber



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber



EU Authorized Representative and Importer:



Teleflex Medical




IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL055933 R01 (2024-03)

Documents / Resources

	<p>ARROW CS-42703-E Three Lumen CVC [pdf] Instructions CS-42703-E Three Lumen CVC, CS-42703-E, Three Lumen CVC, Lumen CVC, CVC</p>
---	--

References

- [User Manual](#)

Manuals+ Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.