




ARROW MP-15703-JA Three-Lumen CVC Instruction Manual


[Home](#) » [Arrow](#) » ARROW MP-15703-JA Three-Lumen CVC Instruction Manual 

Contents

- [1 ARROW MP-15703-JA Three-Lumen CVC](#)
- [2 Documents / Resources](#)
- [3 Related Posts](#)

ARROW®

ARROW MP-15703-JA Three-Lumen CVC



- 1: Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm, Stylet
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 23-5/8" (60 cm) (Angled Soft Tip on One End – Straight Stiff Tip on Other) with Arrow Advancer and hub
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over
- 22 Ga. TW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 2-3/4" (7 cm) XTW

- 1: Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 1-3/8" (3.5 cm) XTW
- 1: Syringe: 10 mL Luer-Slip
- 1: Syringe: 5 mL Luer-Slip
- 1: Syringe: 2 mL Luer-Slip
- 1: Tissue Dilator: 8 Fr. x 9.7 cm
- 3: Dust Cap: Non-Vented
- 1: Second Site™ Adjustable Hub: Fastener
- 1: Second Site™ Adjustable Hub: Catheter Clamp
- 1: Medicine Cup: 60 mL
- 1: Sharps Away® Disposal Cup
- 1: Drape: 35" x 50" (89 cm x 127 cm) with 4" (10 cm) fenestration, with adhesive
- 1: Needle Holder
- 1: Scalpel: #11
- 2: Stopcock: 4-way
- 1: Prep Sponge Swab
- 6: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: OpSite® 4-3/4" x 4" (12 cm x 10 cm)
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Tray: Prep
- 1: Y-hub
- 2: Guidewire Introducer

1A registered trademark of T.J. Smith & Nephew, Limited.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex. Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.



Lumen

**Priming
Volume***
(mL)

**Gravity
Flow
Rate†**
(mL/hr)

Distal (16 Ga.)	0.47	2656
Medial (18 Ga.)	0.4	1245
Proximal (18 Ga.)	0.42	1328

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

Documents / Resources



[ARROW MP-15703-JA Three-Lumen CVC](#) [pdf] Instruction Manual
MP-15703-JA Three-Lumen CVC, MP-15703-JA, Three-Lumen CVC