



## PB-40002-100B Pressure Injectable Arrowg plus and Owner's Manual

[Home](#) » [Arrow](#) » PB-40002-100B Pressure Injectable Arrowg plus and Owner's Manual 

# ARROW

PB-40002-100B



**PR-45563-HPHNL**

(IPNXXXXXX)

**6** Fr.

**3** Lumen

**55** cm  
catheter  
length

**.018** inch dia.  
spring-wire  
guide

**Pressure Injectable Arrowgtard Blue Advance® Three-Lumen PICC with 130  
cm Hydrophilic-Coated Nitinol Guidewire**

### Contents

[1 Contents:](#)

[2 Rx only](#)

[3 Documents /  
Resources](#)

[3.1 References](#)

### Contents:

1: Three Lumen TaperFree® Cathetert}+ with Arrowgtard Blue Advance® Antimicrobial/Antithrombogenic Protection": 6 Fr. (2.11 mm OD) x 55 cm Pressure Injectable and Contamination Guard

1: GlideThru™ Peel-Away Sheath: 6 Fr. x 4" (10 cm) Radiopaque over 6 Fr. Dilator

1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – Straight Stiff Tip on Other)

1: Guidewire, Marked: .018" (0.46 mm) dia. x 130 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil  
 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW  
 1: Syringe: 10 mL Luer-Lock  
 3: Dust Cap: Non-Vented  
 1: Catheter Trimmer  
 1: Safety Scalpel: #11  
 1: Patient ID Card  
 1: Patient Information Booklet  
 1: Paper Tape Measure  
 1: Dressing: STATLOCK®2 PICC Plus Catheter Stabilization Device


\*Licensed under US Patent No. 7,329,412 and 11,219,706.

2A registered trademark of C. R. Bard, Inc.

All components are CE 2797 unless otherwise noted.

## Rx only

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU) California Prop. 65

 **WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

+++MR Safe


Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 30°C (86°F).

Fluid path components are non-pyrogenic.

Contraindications: The Pressure Injectable Arrowgard Blue Advance antimicrobial/antithrombogenic catheter is contraindicated:

- for patients with known hypersensitivity to chlorhexidine
- in the presence of device related infection in the intended insertion vessel or catheter pathway
- in the presence of thrombosis in the intended insertion vessel or catheter pathway


 Lumen	Priming Volume* (mL)	Gravity Flow Rate † (mL/hr)	Pump Flow Rate†† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Pink (18 Ga.)	0.70	990	5440	6
Blue (19 Ga.)	0.52	240	1930	not rated
White (19 Ga.)	0.49	270	2070	not rated


\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head eight.

†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.


\*\* Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



 EU Authorized  
Representative and Importer:

**EC REP** Teleflex Medical  
IDA Business and  
 Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**All components are CE 2797 unless otherwise noted.**


 Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA  
Product of XXX  
Packaged in XXX

**Teleflex®**


[www.teleflex.com](http://www.teleflex.com)

LBL070055 R01 (2024-05)

## Documents / Resources

	<p><b><a href="#">ARROW PB-40002-100B Pressure Injectable Arrowg plus ard</a> [pdf] Owner's Manual</b> PB-40002-100B, PR-45563-HPHNL, LBL070055 R01, PB-40002-100B Pressure Injectable Arrowg plus ard, PB-40002-100B, Pressure Injectable Arrowg plus ard, Injectable Arrowg plus ard, Arrowg plus ard</p>
---	---

## References

-  [p65warnings.ca](http://p65warnings.ca) - [p65warnings Resources and Information.](#)
- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.