



ARROW PSI Kit for use with 7-7.5 Fr. Catheters User Guide


[Home](#) » [Arrow](#) » ARROW PSI Kit for use with 7-7.5 Fr. Catheters User Guide 

Contents

- [1 ARROW PSI Kit for use with 7-7.5 Fr. Catheters](#)
- [2 Features](#)
- [3 Instructions](#)
 - [3.1 Rx only](#)
- [4 Contact us](#)
- [5 Documents / Resources](#)
- [6 Related Posts](#)



ARROW PSI Kit for use with 7-7.5 Fr. Catheters



Features

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer

Needle

- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow®

Raulerson Spring-Wire Introduction Syringe

- 1: Syringe: 3 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Drape: 24" x 36" with 4" fenestration
- 1: Scalpel: #11
- 1: Stopcock: 3-way High-Flow
- 2: Prep Sponge Swab

- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Nylon with Straight Needle
- 1: Tray: Prep

Instructions

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Caution: This Product Contains Dry Natural Rubber.

Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.

Contact us

Arrow International LLC
Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA
Teleflex
www.teleflex.com
LBLO58336 R01 (2022-02)

Documents / Resources

	ARROW PSI Kit for use with 7-7.5 Fr. Catheters [pdf] User Guide Catheters, PSI Kit, 7-7.5 Fr. Catheters, PSI
---	---