

[Manuals.plus](#) /

> [DMI](#) /

> DMI Inflatable Hair Washing Basin (Model 540-8085-0000) - Instruction Manual

DMI 540-8085-0000

DMI Inflatable Hair Washing Basin

Model: 540-8085-0000

Brand: DMI

INTRODUCTION

This instruction manual provides essential information for the safe and effective use of your DMI Inflatable Hair Washing Basin. Designed for individuals who require hair washing while lying down, such as bedridden patients, the disabled, injured, elderly, or children, this portable basin offers a convenient solution. Please read all instructions carefully before use.

PACKAGE CONTENTS

- 1 Inflatable Bed Shampoo Basin
- 1 Inflatable Pillow
- 1 Drain Tube
- 1 PVC Patch Repair Kit
- 1 Instruction Sheet (this document)

NOTE: Please inspect all parts carefully for cracks prior to use.

SETUP INSTRUCTIONS

1. **Locate Inflation Valves:** Find the inflation valves on both the upper and lower portions of the basin.
2. **Inflate Basin:** Inflate the basin by mouth or with an air pump (air pump not included). For extra comfort, inflate the pillow section as well.
3. **Position Basin:** With the user lying flat, gently place the basin under the user's head with their neck resting gently in the outside opening.
4. **Connect Drain Tube:** Connect the larger end of the drain hose to the drain outlet on the underside of the basin. Ensure the drain hose is fully inserted until it cannot go in any further.
5. **Position Drain:** Place the other end of the drain hose into a wastewater container (e.g., a large basin or sink).

Ensure the wastewater container is positioned lower than the basin for proper drainage.



Image: The DMI Inflatable Hair Washing Basin positioned under a user's head for comfortable hair washing.

OPERATING INSTRUCTIONS

6. **Wash Hair:** Wash the user's hair as usual, ensuring water collects in the basin.
7. **Drain Water:** When finished shampooing, open the drain plug to allow all wastewater to completely drain from the basin into the prepared container.
8. **Dry Hair:** Gently dry the user's hair and remove the basin from underneath their head.
9. **Rinse Basin:** Rinse the basin and wipe out any excess water.

EASY TO CLEAN

Connect drain tube to remove water



Image: Hands gently washing hair within the basin, demonstrating the deep design.



Image: A hand opening the drain plug to release water from the basin.

MAINTENANCE

- **Cleaning:** The basin is made of heavy-duty white vinyl and is easy to clean. After each use, rinse thoroughly with water and a mild soap if necessary.
- **Drying:** Leave the basin inflated for several hours until it is completely dry before deflating and folding for storage. This prevents mold and mildew growth.
- **Storage:** Store the deflated and dry basin in a cool, dry place away from direct sunlight and sharp objects.
- **Repair Kit:** In case of small punctures, use the included PVC Patch Repair Kit according to its instructions.

PRODUCT FEATURES

- **Portable Design:** Quickly inflates, deflates, and stores compactly for reuse, making it ideal for home use.
- **Comfortable Support:** Features a built-in pillow for head and neck support, enhancing user comfort during washing.

- **Deep Basin:** Extra deep design prevents water spillage and allows for thorough washing, suitable for various hair types including locs.
- **Durable Construction:** Made from heavy-duty white vinyl for durability and easy cleaning.
- **Convenient Drainage:** Includes a 40-inch drain tube for easy and mess-free water removal.



INCLUDES PILLOW
For Added Head & Neck Support

Image: The inflatable pillow provides added head and neck support for user comfort.

DOUBLE LAYER TUBING

Prevents
spillage
during use.



Image: Double layer tubing design helps prevent water spillage during use.

SPECIFICATIONS

Feature	Detail
Model Number	540-8085-0000
Product Dimensions	24 x 20 x 8 inches (approximate, as per product description)
Material	Heavy-duty White Vinyl
Drain Tube Length	40 inches
Weight	11.2 ounces



Image: Visual representation of the basin's dimensions.

TROUBLESHOOTING

- **Water Leakage:** Ensure all inflation valves are securely closed and the drain plug is firmly in place. Check for any punctures and use the repair kit if needed.
- **Slow Drainage:** Verify that the drain tube is not kinked or obstructed. Ensure the wastewater container is positioned lower than the basin.
- **Difficulty Inflating:** Ensure the inflation valves are fully open during inflation. If using a pump, check its connection.

OFFICIAL PRODUCT VIDEOS

VERY Comfortable inflatable shampoo bed. Easy assemble!

Your browser does not support the video tag.

Video Description: This video demonstrates the ease of assembly and comfort of the DMI Inflatable Hair Washing Basin. It highlights the inflation process and how the basin supports the user's head and neck during hair washing.

Apple cider vinegar rinse with Duro-med portable inflatable shampoo bowl

Your browser does not support the video tag.

Video Description: This video showcases the DMI Inflatable Hair Washing Basin being used for an apple cider vinegar rinse, demonstrating its practical application for hair treatments while lying down.

WARRANTY AND SUPPORT

For warranty information or product support, please refer to the manufacturer's official website or contact DMI customer service directly. Keep your purchase receipt for any warranty claims.

Manufacturer: DMI

ASIN: B0009STN7O