

[Manuals.plus](#) /

› [Drive Medical](#) /

› Drive Medical Aluminum Rollator Walker R726BK Instruction Manual

Drive Medical R726BK

Drive Medical Aluminum Rollator Walker (Model R726BK) Instruction Manual

Comprehensive guide for setup, operation, and maintenance of your Drive Medical Rollator Walker.

PRODUCT OVERVIEW

The Drive Medical Aluminum Rollator Walker is designed to provide mobility support for individuals requiring assistance with walking. It features a lightweight yet sturdy aluminum frame, a comfortable padded seat with a backrest, and convenient storage. This rollator is suitable for both indoor and outdoor use, offering stability and ease of movement.



Figure 1: Drive Medical Aluminum Rollator Walker, Black model.

SAFETY INFORMATION

- **Brakes:** Always engage the loop locks on the handles by pressing them down to lock the wheels and prevent unintended movement, especially when sitting or standing. The serrated edges on the brakes provide a firm hold.
- **Weight Capacity:** This rollator supports users up to 300 lb. Do not exceed this weight limit.
- **Stability:** Ensure all four wheels are on a stable, level surface before applying weight or sitting.
- **Storage Pouch:** Do not overload the zippered storage pouch under the seat. Distribute weight evenly.
- **Assembly:** Ensure all components are securely fastened before use. Refer to the setup section for proper assembly.

FEATURES

- **Padded Seat and Backrest:** Provides comfortable resting options during use. The backrest is removable and hinged, allowing it to be folded up or down.
- **Lightweight Aluminum Frame:** Weighs approximately 15.6 lb, making it easy to maneuver and transport.
- **6-Inch Wheels:** Equipped with soft-grip tires for smooth movement across various indoor and outdoor surfaces.
- **Adjustable Ergonomic Handles:** Easily customize handle height with a simple lever for a comfortable and personalized fit.
- **Side-to-Side Folding Design:** Allows for one-handed folding, enabling the rollator to stand upright when folded for convenient storage and transport.
- **Under-Seat Storage Pouch:** A zippered pouch provides secure and private storage for personal items.

SETUP

The Drive Medical Aluminum Rollator Walker comes largely pre-assembled. Follow these steps to prepare it for use:

1. **Unfolding:** Carefully remove the rollator from its packaging. Gently pull the sides of the frame apart until it fully extends and locks into place. Ensure the seat is properly positioned.
2. **Handle Installation/Adjustment:** Locate the adjustment levers or knobs on the handle posts. Press or turn to release, slide the handles to the desired height (typically wrist height when standing upright), and secure them firmly.
3. **Backrest Attachment:** If the backrest is not already attached, secure it to the designated points on the frame. It is designed to be hinged for convenience.
4. **Brake Check:** Before first use, test both hand brakes to ensure they engage and disengage smoothly and effectively lock the wheels.

Customized Fit



Easily adjust the handle height between 32"-37" inches for a comfortable reach



Figure 2: Adjusting the handle height for personalized comfort.

OPERATING INSTRUCTIONS

1. **Walking:** Stand upright behind the rollator, holding both handles. Push the rollator forward a comfortable distance, then step forward. Maintain a steady pace and keep your body close to the rollator for optimal balance.
2. **Using Brakes:** To slow down or stop, gently squeeze the hand brake levers. To lock the wheels for sitting or standing, push the brake levers down until they click into the locked position. Pull them up to release.
3. **Sitting:** Ensure the rollator is on a level surface. Engage both hand brakes by pushing them down to lock the wheels. Turn around and slowly lower yourself onto the padded seat. Use the backrest for support.
4. **Standing Up:** From a seated position, lean forward slightly and push off the handles to stand up. Once standing, release the brakes.
5. **Folding for Storage/Transport:** To fold the rollator, lift the padded seat by pulling the strap located in the center of the seat. The frame will collapse inward. The side-to-side folding design allows it to stand upright for easy storage or lifting into a vehicle.

Your browser does not support the video tag.

Video 1: Demonstration of the Drive Medical Foldable Rollator Walker with Seat, showing its features and folding mechanism.

Enjoy Independent Mobility




Gives individuals with balance, strength, and other difficulties support and stability to navigate life more safely

drive

DeVilbiss
HEALTHCARE

Figure 3: User demonstrating outdoor mobility with the rollator walker.

All Day Comfort




1



Padded seat and backrest let you rest on the go

2



Contoured handgrips reduce stress on hands

Figure 4: The padded seat and backrest provide a convenient resting spot.

Easy to Store & Transport



At just 13 lb., the aluminum frame folds compact to fit in trunks and small spaces




Figure 5: The compact folding design allows for easy transport in a car trunk.

MAINTENANCE

- **Cleaning:** Wipe down the frame and seat regularly with a damp cloth and mild soap. Avoid abrasive cleaners.
- **Wheel Inspection:** Periodically check the 6-inch wheels for wear and tear. Ensure they spin freely and are free of debris.
- **Brake Adjustment:** If the brakes feel loose or do not engage properly, they may need adjustment. Consult a qualified technician or the manufacturer's support for guidance.
- **Fastener Check:** Regularly inspect all nuts, bolts, and fasteners to ensure they are tight and secure.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Rollator does not roll smoothly.	Debris in wheels, loose wheel fasteners, worn wheels.	Clean wheels, tighten fasteners, replace worn wheels if necessary.

Problem	Possible Cause	Solution
Brakes are not engaging or locking effectively.	Loose brake cables, worn brake pads, improper adjustment.	Check brake cable tension. If unsure, consult a professional for adjustment or replacement.
Rollator feels unstable.	Loose frame components, uneven handle height, improper unfolding.	Ensure frame is fully locked open. Check and tighten all visible fasteners. Verify handle heights are even and secure.

SPECIFICATIONS

Model Number: R726BK

Brand: Drive Medical

Color: Black

Material: Alloy Steel, Aluminum

Item Weight: 12 Pounds (approx. 5.44 kg)

Weight Capacity: Up to 300 Pounds (approx. 136 kg)

Product Dimensions (Overall): 24.5"D x 27.5"W x 32"H (approx. 62.2 cm D x 69.8 cm W x 81.3 cm H)

Wheel Size: 6 Inches

Dimensions





Supports up to 300 lbs.

Figure 6: Key dimensions of the rollator walker.

WARRANTY AND SUPPORT

Drive Medical is committed to providing quality medical equipment. For specific warranty details, product support, or to inquire about replacement parts, please contact Drive Medical customer service directly. Information regarding their products and support can often be found on their official website or through the retailer where the product was purchased.

Drive Medical has been manufacturing trusted medical equipment and homecare products since 2000, focusing on enhancing lives with products that support safety, independence, and comfort.

