

[Manuals.plus](#) /

> [Welch Allyn](#) /

> Welch Allyn LumiView Battery Pack Model 75200 Instruction Manual

Welch Allyn 75200

Welch Allyn LumiView Battery Pack Model 75200 Instruction Manual

Portable Power Source

1. INTRODUCTION

This instruction manual provides essential information for the safe and effective use of the Welch Allyn LumiView Battery Pack, Model 75200. This device is designed as a portable power source for compatible Welch Allyn LumiView products. Please read this manual thoroughly before operating the device and retain it for future reference.

2. SAFETY INFORMATION

WARNING: Failure to follow these safety instructions could result in fire, electric shock, injury, or damage to the product or other property.

- Use only with Welch Allyn LumiView compatible devices.
- Do not expose the battery pack to extreme temperatures, direct sunlight, or moisture.
- Do not disassemble, puncture, or modify the battery pack.
- Keep out of reach of children.
- Dispose of the battery pack according to local regulations. Do not dispose of in household waste.
- If the battery pack shows signs of damage, leakage, or overheating, discontinue use immediately and contact Welch Allyn support.

3. SETUP

3.1 Unpacking

Carefully remove the Welch Allyn LumiView Battery Pack from its packaging. Inspect the device for any visible damage. If damage is present, do not use the device and contact your supplier.

3.2 Initial Charge

Before first use, fully charge the battery pack. Connect the appropriate Welch Allyn charger (sold separately) to the charging port on the battery pack. Plug the charger into a suitable power outlet. The charging indicator (if present) will illuminate or change color to indicate charging status. Refer to your charger's manual for specific charging instructions.

3.3 Connecting to a Device

Align the connector of the LumiView device with the output port on the battery pack. Gently push until a secure connection is established. Ensure the connection is firm to prevent accidental disconnection during use.



Figure 1: Top view of the battery pack, illustrating the power switch and the output port for connecting to a LumiView device.

4. OPERATING INSTRUCTIONS

4.1 Powering On/Off

To power on the battery pack, locate the power switch (refer to Figure 1) and slide it to the 'ON' position. To power off, slide the switch to the 'OFF' position. Ensure the device is powered off when not in use to conserve battery life.

4.2 Charging the Battery Pack

When the battery level is low, connect the Welch Allyn charger to the charging port. A full charge typically takes several hours. Avoid overcharging by disconnecting the charger once the battery is fully charged, if the charger does not have an automatic cut-off feature.

4.3 Battery Status Indicators

The battery pack may feature an indicator light to show its charge status. Consult the specific LumiView device or charger manual for details on interpreting these indicators (e.g., solid light for charged, blinking for charging, red for low battery).

5. MAINTENANCE AND CARE

5.1 Cleaning

Clean the exterior of the battery pack with a soft, damp cloth. Do not use abrasive cleaners, solvents, or immerse the device in water. Ensure all ports are dry before connecting any cables.

5.2 Storage

Store the battery pack in a cool, dry place away from direct sunlight and extreme temperatures. For long-term storage, it is recommended to charge the battery pack to approximately 50% every few months to maintain battery health.

5.3 Battery Life

Rechargeable batteries have a limited number of charge cycles. If the battery pack's operating time significantly decreases, it may be time to replace it. Contact Welch Allyn support for replacement options.

6. TROUBLESHOOTING

Problem	Possible Cause	Solution
Battery pack does not power on.	Battery is depleted. Power switch is off.	Charge the battery pack. Ensure the power switch is in the 'ON' position.
Battery pack does not charge.	Charger not connected properly. Faulty charger or cable. Battery pack damage.	Ensure charger is securely connected. Try a different charger/cable if available. Contact support if issue persists.
LumiView device not receiving power.	Battery pack is off or depleted. Connection is loose.	Turn on the battery pack and ensure it is charged. Check the connection between the battery pack and the LumiView device.

If troubleshooting steps do not resolve the issue, contact Welch Allyn customer support for assistance.

7. SPECIFICATIONS

Model: 75200
Brand: Welch Allyn
Type: Portable Power Source / Battery Pack
Input: 5.3V
Output: 3.5V — 1.8Ah
Manufacturer: Welch Allyn
ASIN: B002MBE9KW
Unit Count: 1 Count
Date First Available: March 11, 2011



Figure 2: Rear view of the battery pack, displaying the model number 75200 and key electrical specifications.

8. WARRANTY AND SUPPORT

Welch Allyn products are manufactured to high-quality standards. For information regarding warranty coverage, please

refer to the warranty card included with your purchase or visit the official Welch Allyn website. For technical support, service, or replacement parts, please contact Welch Allyn customer service.

Welch Allyn Customer Support:

Website: www.welchallyn.com (Note: This is a placeholder URL. Please refer to official documentation for the correct support contact.)

Phone: Refer to your product documentation for regional contact numbers.

© 2024 Welch Allyn. All rights reserved. Information subject to change without notice.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question