

DonJoy 11-1268

DonJoy McGuire Knee Wrap-On Pad Instruction Manual

Model: 11-1268

[Introduction](#)[Safety Information](#)[Package Contents](#)[Setup](#)[Operating Instructions](#)[Maintenance](#)[Troubleshooting](#)[Specifications](#)[Warranty & Support](#)

1. INTRODUCTION

The DonJoy McGuire Knee Wrap-On Pad is designed for targeted cold therapy application to the knee. This pad features an open patella design and an integrated elastic wrap, allowing for direct application to the knee without the need for additional securing wraps. It is intended for use with compatible DonJoy IceMan Cold Therapy Units, such as the CLEAR3 or CLASSIC3 models (sold separately).

This manual provides essential information for the proper setup, operation, and maintenance of your McGuire Knee Wrap-On Pad to ensure effective and safe use.

2. IMPORTANT SAFETY INFORMATION

- Always consult with a healthcare professional before beginning any cold therapy treatment.
- Do not apply the cold pad directly to skin for extended periods without a protective barrier, as this may cause skin irritation or frostbite. Follow your healthcare provider's recommendations.
- Regularly inspect the pad for any signs of damage or leaks. Discontinue use if damage is observed.
- This product is not sterile and should not be used in applications requiring sterile conditions.
- Keep out of reach of children and pets.
- Do not modify the pad or its components. Use only with compatible DonJoy IceMan Cold Therapy Units.

3. PACKAGE CONTENTS

- 1 x DonJoy McGuire Knee Wrap-On Pad (Model 11-1268)
- Integrated elastic wrap for securing the pad

Note: DonJoy IceMan Cold Therapy Unit (e.g., CLEAR3 or CLASSIC3) is sold separately and is required for operation.

4. SETUP INSTRUCTIONS

Follow these steps to prepare your McGuire Knee Wrap-On Pad for use with a compatible DonJoy IceMan Cold Therapy Unit.

1. **Prepare the Cold Therapy Unit:** Fill your DonJoy IceMan Cold Therapy Unit with ice and water according to its specific instructions. Ensure the unit is properly assembled and ready for connection.
2. **Connect the Pad:** Locate the two connectors on the McGuire Knee Wrap-On Pad. Connect these securely to the corresponding ports on your DonJoy IceMan Cold Therapy Unit. Ensure a snug fit to prevent leaks.
3. **Purge Air from System (Important):** Before applying the pad, it is recommended to purge any air from the system for optimal performance.
 - Place the cold therapy unit higher than the pad (e.g., on a chair with the pad on the floor).
 - Turn on the cold therapy unit. Observe the hoses and pad for air bubbles. Allow the unit to run until all visible air bubbles are cleared from the system.
 - Once air is purged, the unit can be placed on the floor or at a level convenient for use.



5. OPERATING INSTRUCTIONS

Once the pad is connected and the system is purged of air, you can begin your cold therapy session.

1. **Position the Pad:** Place the McGuire Knee Wrap-On Pad directly onto the knee area requiring cold therapy. The pad is specifically designed with an open patella to conform to the knee's contours.
2. **Secure the Pad:** Use the integrated elastic wrap to secure the pad firmly but comfortably around your knee. Ensure the pad remains in place throughout the therapy session. Avoid overtightening, which could restrict circulation.
3. **Begin Therapy:** Turn on your DonJoy IceMan Cold Therapy Unit. The cold water will begin to circulate through the pad, providing therapeutic cooling to your knee.
4. **Monitor and Adjust:** Periodically check your skin under the pad for any signs of discomfort, excessive cold, or skin changes. Adjust the pad's position or the duration of therapy as advised by your healthcare professional.
5. **Duration of Use:** Follow your healthcare provider's recommendations for the duration and frequency of cold therapy sessions. Do not exceed recommended times to prevent adverse effects.

6. CARE AND MAINTENANCE

Proper care and maintenance will extend the life of your DonJoy McGuire Knee Wrap-On Pad.

- **Cleaning:** After each use, disconnect the pad from the cold therapy unit. Drain any remaining water. Clean the exterior of the pad with a mild soap and water solution. Rinse thoroughly and allow to air dry completely before storage. Do not submerge the entire pad in water for cleaning.
- **Drying:** Ensure the pad is completely dry, both internally and externally, before storing to prevent mold or mildew growth.
- **Storage:** Store the pad in a clean, dry place away from direct sunlight and extreme temperatures. Avoid folding or creasing the pad sharply, which could damage the internal tubing.
- **Inspection:** Regularly inspect the pad and hoses for any signs of wear, cracks, or leaks. If any damage is found, discontinue use and replace the pad.
- **Non-Sterile:** This pad is manufactured as a non-sterile product. It cannot be sterilized for surgical or sterile applications.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
Pad is not getting cold enough.	Insufficient ice in the unit; air in the system; unit not running properly.	Ensure the cold therapy unit is filled with adequate ice and water. Re-purge air from the system (see Setup). Check the cold therapy unit's manual for its specific troubleshooting.
Water is leaking from connections.	Loose connections; damaged O-rings or connectors.	Ensure all connections between the pad and the cold therapy unit are securely fastened. Inspect connectors for damage. If leaks persist, the pad or unit may require replacement.
Water is leaking from the pad itself.	Pad material or seams are damaged.	Discontinue use immediately. The pad is damaged and cannot be repaired. Replace the pad.

Problem	Possible Cause	Solution
Pad feels uncomfortable or causes skin irritation.	Direct contact with skin for too long; pad applied too tightly; allergic reaction.	Always use a protective barrier between the pad and skin. Loosen the elastic wrap if it feels too tight. If irritation persists, discontinue use and consult a healthcare professional.

8. PRODUCT SPECIFICATIONS

Feature	Detail
Model Number	11-1268
Product Dimensions	11 x 10 x 2 inches (27.94 x 25.4 x 5.08 cm)
Item Weight	11.2 ounces (0.32 kg)
Target Body Part	Knee
Special Features	Integrated Elastic Wrap, Open Patella Design
Compatibility	DonJoy IceMan Cold Therapy Units (CLEAR3, CLASSIC3)
Sterility	Non-sterile
Color	Gray

9. WARRANTY AND SUPPORT

Specific warranty information for the DonJoy McGuire Knee Wrap-On Pad is typically provided with the product packaging or can be obtained directly from DonJoy or your authorized distributor. Please retain your proof of purchase for warranty claims.

For technical support, product inquiries, or to report issues, please contact DonJoy customer service. Contact details can usually be found on the official DonJoy website or on the packaging of your cold therapy unit.

Manufacturer: DonJoy Braces

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question