

[Manuals.plus](#) /

› [Rayovac](#) /

› Rayovac Proline Advanced Mercury-Free Hearing Aid Batteries Size 312 User Manual

Rayovac IM312A-48ZM

Rayovac Proline Advanced Mercury-Free Hearing Aid Batteries Size 312 User Manual

INTRODUCTION

Welcome to the user manual for your Rayovac Proline Advanced Mercury-Free Hearing Aid Batteries, Size 312. These batteries are designed to provide reliable and long-lasting power for your hearing aids. Please read this manual carefully to ensure proper use, optimal performance, and safety.

IMPORTANT SAFETY INFORMATION

- Keep out of reach of children.
- If swallowed, promptly seek doctor; have doctor phone 800-498-8666.
- Do not dispose of in fire or charge – may explode or leak causing injury.
- Immediately dispose of used batteries.
- Store batteries in a cool, dry place, free of moisture and extreme temperatures.
- Do not remove the protective tab until ready to use.

WARNING: Batteries can be dangerous if misused. Always follow these safety guidelines.

PRODUCT OVERVIEW

The Rayovac Proline Advanced Mercury-Free Hearing Aid Batteries are zinc-air batteries, specifically designed for hearing aids. They are size 312 and come in child-resistant packaging.

- **Battery Type:** Zinc-Air
- **Size:** 312
- **Voltage:** 1.45 Volts
- **Reusability:** Single Use
- **Key Feature:** Mercury-Free formulation.



Image shows the front of the Rayovac Proline Advanced Hearing Aid Batteries packaging, highlighting "Size 312" and "Child Resistant Packaging".

SETUP AND INSTALLATION

1. **Opening Child-Resistant Packaging:** The batteries are supplied in child-resistant packaging. To access a battery, you will need scissors. Carefully cut around the individual battery compartment on the plastic wheel.

Your browser does not support the video tag.

This video demonstrates how to open the child-resistant packaging for Rayovac hearing aid batteries, showing the use of scissors to cut out an individual battery from the dispenser wheel.

2. **Battery Activation:** Once removed from the packaging, peel off the protective tab from the battery. For best results, allow the battery to sit un-tabbed for at least 60 seconds before inserting it into your hearing aid. This allows the zinc-air battery to fully activate.

NO HEARING AID BATTERY
LASTS LONGER*



RAYOVAC[®]

PRO LINE[®] ADVANCED
HEARING AID BATTERIES

Size
312

CHILD RESISTANT PACKAGING
CUT TO OPEN



6
PACK
Zinc Air
Batteries

Image shows a single Rayovac Size 312 battery with its orange protective tab, ready for activation.

3. **Inserting into Hearing Aid:** Open the battery compartment of your hearing aid. Insert the battery with the positive (+) side facing the correct direction as indicated in your hearing aid's manual. Close the battery compartment securely.

OPERATING INSTRUCTIONS

After installation, turn on your hearing aid. The battery should provide consistent power. Battery life varies depending on hearing aid usage, settings, and environmental factors.

Your hearing aid may provide an audible alert when the battery level is low, indicating it's time for a replacement.

MAINTENANCE AND STORAGE

- **Storage:** Keep unused batteries in their original, sealed packaging in a cool, dry place at room temperature. Avoid extreme heat, cold, or humidity.
- **Handling:** Always handle batteries with clean, dry hands.
- **Shelf Life:** These batteries have a long-lasting 4-year shelf life when stored correctly.
- **Used Batteries:** Place used batteries back into the dispenser wheel or a designated container to avoid confusion with new batteries.

Size 312 6 PACK
Zinc Air
Batteries

312CR-6PROUS
1.45V



HEARING AID BATTERIES



For best results, let a new battery sit un-tabbed for one minute before inserting into device.

Warning: (1) **Keep out of reach of children.** Keep in original package until use. If swallowed, promptly see doctor; have doctor phone 800-498-8666. (2) do not dispose of in fire or charge – may explode or leak causing injury. Immediately dispose of used batteries.



REPLACES ALL SIZE 312

©2024 Energizer. Trademarks and certain designs are owned by Energizer. Made in the United Kingdom for Energizer Brands, LLC, St. Louis, MO 63141



Cut to Open



RAYOVAC
PRO LINE[®] ADVANCED

*Based on ANSI/IEC
Performance Tests
for Size 312

Mercury Free
UPN-175875

BEST BEFORE



warnings.

TROUBLESHOOTING COMMON ISSUES

- **Hearing Aid Not Working After Battery Change:**

- Ensure the battery is inserted correctly with the positive (+) side facing the right direction.
- Check if the protective tab was fully removed.
- Allow at least 60 seconds for the battery to activate after removing the tab.
- Try a fresh battery from a new pack.

- **Short Battery Life:**

- Verify the battery size (312) is correct for your hearing aid.
- Ensure proper activation time (60 seconds) after removing the tab.
- Check hearing aid settings; high volume or streaming can drain batteries faster.
- Store batteries correctly to prevent premature discharge.

- **Difficulty Opening Packaging:**

- Use sharp scissors to carefully cut around the individual battery on the plastic wheel. Refer to the "Setup and Installation" section for a video demonstration.

If problems persist, consult your hearing aid professional or contact Rayovac customer support.

PRODUCT SPECIFICATIONS

Feature	Detail
Brand	Rayovac
Model Number	IM312A-48ZM
Battery Size	312
Battery Cell Type	Zinc
Voltage	1.45 Volts
Reusability	Single Use
Item Dimensions (Packaging)	5 x 6 x 0.5 inches
Battery Weight	7 Grams (per battery)
UPC	012800512164, 012800499953

DISPOSAL AND RECYCLING

Please dispose of used batteries responsibly according to local regulations. Do not throw them in household trash. Many communities offer battery recycling programs. Contact your local waste management facility for information on proper disposal and recycling options for zinc-air batteries.

WARRANTY AND SUPPORT

For product support or warranty inquiries, please contact Rayovac customer service. Refer to the product packaging or

Rayovac's official website for the most current contact information.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question