

**Medline EXTSRB3036AZ**

# Medline Ultrasorbs Premium Underpads Instruction Manual

## INTRODUCTION

---

Medline Ultrasorbs Premium Underpads are designed to provide superior skin dryness and comfort. They feature an air-permeable, moisture-proof backsheet and an absorbent core that effectively wicks moisture away from the skin. These underpads are suitable for various applications, including incontinence management, use with ICU/CCU patients, emergency room fluid absorption, labor and delivery, and low air loss therapy, ensuring enhanced patient dignity and skin care.



This image displays several folded Medline Ultrasorbs Premium Underpads, highlighting their compact form and the visible 'ULTRASORBS AP' branding on the top layer. The pads appear thick and designed for absorbency.

## FEATURES

- **Superabsorbent Polymer Core:** Wicks away and absorbs fluid to help protect patient or resident skin while controlling odors.
- **Soft Materials:** Provide protection against friction and shear, contributing to patient or resident comfort.
- **Air-Permeable:** Materials are suitable for use on low air loss mattresses and specialty support surfaces.
- **Non-Liquid Permeable Backsheet:** Protects underlying surfaces and helps keep them dry.
- **Lies Flat:** Designed to remain flat when wet without bunching, swelling, or disintegrating.

## SETUP AND PREPARATION

Before use, ensure the surface where the underpad will be placed is clean and dry. Carefully unfold the underpad completely, ensuring the absorbent side (typically the softer, non-plastic side) is facing upwards. Position the underpad to cover the desired area adequately.

## OPERATING AND USAGE

Position the Medline Ultrasorbs Premium Underpad directly beneath the individual, ensuring adequate coverage of the area prone to fluid loss. The underpad's design allows for effective absorption and containment of fluids, protecting bedding and other surfaces. For optimal performance, ensure the underpad lies flat and is not creased during use. Replace the underpad promptly once it becomes soiled or saturated to maintain hygiene and skin integrity.



This image shows an elderly woman seated in a wheelchair, smiling, with a younger woman standing behind her, also smiling. They are in a comfortable home setting, suggesting the application of Medline products for home care.

## MAINTENANCE AND CARE

Medline Ultrasorbs Premium Underpads are designed for single use. Once soiled, carefully fold the underpad inward to contain any fluids and dispose of it according to local waste disposal guidelines. Do not attempt to wash or reuse the underpads, as this will compromise their absorbency and protective properties, potentially leading to leakage and reduced hygiene.

## TROUBLESHOOTING

- Issue: Underpad is bunching or not staying flat.**  
*Solution:* Ensure the underpad is fully unfolded and smoothed out before placement. While these underpads are designed to lie flat, if the user is highly mobile, consider repositioning or ensuring the underpad is placed on a stable, flat surface.
- Issue: Fluid leakage.**  
*Solution:* Verify that the underpad is correctly positioned to cover the entire area of fluid loss. Ensure the appropriate size of underpad is being used for the expected volume of fluid. Replace the underpad promptly after significant fluid absorption, as over-saturation can lead to leakage.

## SPECIFICATIONS




Attribute	Detail
Product Dimensions	36 x 30 x 0.01 inches
Item Weight	12.8 ounces

Item Model Number	EXTSRB3036AZ
Date First Available	February 21, 2013
Manufacturer	Medline Industries Healthcare
ASIN	B00JIRAS94
Brand	Medline
Material	Bamboo
Age Range (Description)	Adult
Number of Layers	2
Unit Count	70.0 Count

## WARRANTY AND SUPPORT

For specific warranty information or product support regarding Medline Ultrasorbs Premium Underpads, please refer to the official Medline website or contact Medline customer service directly. General product inquiries and assistance can often be resolved through the manufacturer's dedicated support channels.

## Related Documents

	<p><a href="#">Medline Bariatric Bed MDR107004: Owner's Assembly and Operating Manual</a></p> <p>Comprehensive owner's assembly and operating manual for the Medline Bariatric Bed (Model MDR107004), including safety warnings, assembly instructions, operation details, maintenance procedures, troubleshooting guide, parts list, and limited warranty information.</p>
	<p><a href="#">Medline Physician Office and Patient Care Products Catalog</a></p> <p>A comprehensive catalog of medical supplies and equipment for physician offices and patient care, including wound care, personal care, mobility aids, and more, from Medline Industries.</p>
	<p><a href="#">Medline 2024 Sustainability Report: Driving Resilient and Responsible Healthcare</a></p> <p>Explore Medline's 2024 Sustainability Report detailing efforts to foster a resilient and sustainable healthcare supply chain, focusing on responsible products, environmental action, and community engagement.</p>



### [Medline Diagnostic and Laboratory Products Catalog](#)

A comprehensive catalog showcasing Medline's extensive range of diagnostic instruments, patient monitors, blood pressure equipment, thermometers, stethoscopes, defibrillators, electrodes, point-of-care testing devices, and laboratory supplies, designed for healthcare professionals.